

Neurotransmitters

Key Clinical Messages

What is the Neurotransmitter Test?

Vibrant's Neurotransmitter Test is a urine-based test that measures levels of the most common neurotransmitters and metabolites in the brain and nervous system. It includes inhibitory, excitatory, and other neurotransmitters along with diurnal rhythm of epinephrine and norepinephrine, providing insight into how the nervous system is communicating.

What are Neurotransmitters?

Neurotransmitters are endogenous chemical messengers that are responsible for communication in the nervous system. They are secreted by neurons in response to an electrical signal and cross a synapse to communicate with other cells, such as neurons, glands, or muscle cells.

Neurotransmitters are involved in many physiological processes in the central nervous system and in the periphery, including regulating mood, heart rate, sleep regulation, digestion, appetite, muscle movements, breathing, learning, memory, motivation, pain, and many more functions.

Neurotransmitters are classified as either excitatory, inhibitory, or neuromodulatory. Many neurotransmitters function both as excitatory and inhibitory depending on the location and circumstances.

Inhibitory Neurotransmitters: Acts as the "brakes" or "off switch" by decreasing/inhibiting action potentials. The following neurotransmitters are classified as inhibitory:

- GABA, Serotonin, Glycine, Taurine

Excitatory Neurotransmitters: Act as the "accelerator" or "on switch" by firing off action potentials. The following neurotransmitters are classified as excitatory:

- Glutamate, Epinephrine, Norepinephrine, Dopamine, Acetylcholine

Why Order the Neurotransmitters Test?


The neurotransmitter test is a non-invasive test that provides a snapshot of an individual's specific neurotransmitter levels. This test provides insight into how the nervous system is communicating.

Many treatment recommendations used with neurological or mental conditions are currently guided by a general diagnosis or a "trial and error" approach. The neurotransmitter test can provide clinically useful data to help guide you when monitoring and treating neurological or mental health conditions. It can also be used to assess medication efficacy and responsiveness or help wean medications or supplements.

Neurotransmitter imbalances can be caused by a host of factors. The neurotransmitter test allows you to assess all aspects of neurotransmitter synthesis and metabolism, highlighting where issues occur on the neurotransmitter pathways so interventions can be tailor-made.

Regulatory Statement:

This test has been laboratory developed and their performance characteristics determined by Vibrant America LLC, a CLIA-certified laboratory performing the test CLIA#:05D2078809. The test has not been cleared or approved by the U.S. Food and Drug Administration (FDA). Although FDA does not currently clear or approve laboratory-developed tests in the U.S., certification of the laboratory is required under CLIA to ensure the quality and validity of the tests.





Which Patients Benefit from This Test?

- ADHD/ADD
- Addictive behaviors
- Alzheimer's disease
- Altered pain response
- Anger
- Anxiety
- Appetite (poor/excess)
- Autism spectrum disorder
- Autoimmune diseases
- Autonomic nervous system disorders
- Bipolar disorder
- Brain fog
- Cancer
- Cardiovascular disease
- Chronic fatigue
- Constipation/ diarrhea
- Compulsive behaviors
- Dementia
- Depression
- Developmental disorders
- Difficulty concentrating
- Diabetes
- Dysmotility
- Eating disorders
- Epilepsy/seizure disorders
- Fatigue
- Fibromyalgia
- Gastroesophageal reflux disease
- Headaches/migraines
- Hormonal Imbalances
- Hyperactivity
- Huntington's disease
- Irritable bowel syndrome (IBS)
- Insomnia
- Irritable bowel disease (IBD)
- Irritability
- Low libido
- Low motivation
- Medication adjustments
- Memory impairments
- Mood disorders
- Movement disorders/ Motor dysfunction
- Multiple sclerosis
- Muscle twitching/spasms
- Obsessive compulsive disorder
- Panic attacks
- Parkinson's disease
- Schizophrenia
- Sensory processing disorders
- Tremors
- Weight imbalances
- Vomiting

Which Tests Pair Well with the Neurotransmitter Test?

- **Micronutrient test:** many enzymatic reactions for neurotransmitter synthesis require nutrient cofactors, therefore, it's extremely helpful to have information on nutrient levels when interpreting the neurotransmitter test. The most important micronutrients involved in neurotransmitter pathways include:
 - Vitamin B1, B2, B3, B5, B6, folate, choline, iron, magnesium, copper, vitamin C, vitamin D, amino acids
- **Gut Zoomer:** to assess for impaired digestion and absorption, which could lead to decreased amino acids and nutrients required for neurotransmitter synthesis.
- **Neural Zoomer Plus:** to assess for antibodies to different neurotransmitter receptors.
- **Methylation Panel:** to assess for impairments in the methylation cycle affecting neurotransmitter synthesis. This is important since many neurotransmitters rely on methylation support for either synthesis or degradation.
- **Total Tox Burden:** to assess for toxins that may interfere with neurotransmitter synthesis, receptor function, and neurotransmitter metabolism.
- **Hormone tests:** to assess for hormone imbalances since hormones can influence many neurotransmitter levels.
- **Wheat Zoomer:** to assess for gluten sensitivity that may affect the glutamic acid decarboxylase (GAD) enzyme, affecting glutamate and GABA levels.

Why Vibrant?

Vibrant is a **CLIA-certified lab** that utilizes reliable, high-quality methodologies to measure individual neurotransmitter levels.

Methodology

The Neurotransmitter test uses liquid chromatography with tandem mass spectrometry (LC-MS/MS). The analyte results are expressed by normalizing to the quantity of creatinine measured to account for urine dilution variations.

Test Prep

Collection: Four (4) urine specimen tubes.

- Requires four (4) urine collection samples to be collected throughout the day: first-morning urine, 2 hours after waking, evening (before dinner), night (before bed)
- For pediatric patients, a pediatric urine collection bag may be used
- Does not need to be refrigerated or frozen

Hydration: Do not drink more than 8 oz water 1 hour prior to each urine collection. Samples may be rejected if the urine is too dilute.

Fasting: Not required.

Diet: It's important to avoid certain foods because they can be a source of neurotransmitters that can artificially skew the results. The following foods must be avoided for 48 hours before collection:

- Bananas, pineapple, nuts, alcohol, protein powder, protein shakes

Medication Restrictions: None.

Dietary Supplement Restrictions: None.

Contraindications: Renal disease or other kidney abnormalities.



Reference Ranges

Reference ranges were established using a cohort of 1000 apparently healthy individuals. Reference ranges have been validated for patients aged <10 years, 10-14 years, and >14 years. The classification of red indicates a result that is outside the reference range and the classification of green denotes a result that is within the reference range.

Vibrant has a comprehensive interpretive guide for the Neurotransmitters Test to aid providers in the interpretation of their patient's results.



Neurotransmitters Summary Report

Test Name	Current	Previous	Result	Reference
Acetylcholine (mcg/g)	8.64			1.7-5.9
<p>Acetylcholine is the neurotransmitter used at the neuromuscular junction where it is released from the motor neurons of the nervous system in order to activate muscles. Acetylcholine functions in both the central nervous system (CNS) and the peripheral nervous system (PNS). In the CNS, cholinergic projections from the basal forebrain to the cerebral cortex and hippocampus support the cognitive functions of those target areas. In the PNS, acetylcholine activates muscles and is a major neurotransmitter in the autonomic nervous system. Acetylcholine has been implicated in learning and memory in several ways. In animals, disruption of the supply of acetylcholine impairs the learning of simple discrimination tasks. Acetylcholine is rich in food sources such as eggplant, bitter orange, common bean, foxglove, mistletoe, mung bean, nettle species, pea, radish, spinach, squash, wild strawberry.</p>				
5-HTP (mcg/g)	266.39			11.4-185.6
<p>5-Hydroxytryptophan (5-HTP), also known as oxitriptan, is a naturally occurring amino acid and chemical precursor as well as a metabolic intermediate in the biosynthesis of the neurotransmitter serotonin. 5-HTP is produced from the amino acid tryptophan through the action of the enzyme tryptophan hydroxylase. 5-HTP is normally rapidly converted to 5-HT by amino acid decarboxylase. 5-HTP is sold over the counter as a dietary supplement for use as an antidepressant, appetite suppressant, and sleep aid. Oral 5-HTP results in an increase in urinary 5-HIAA, a serotonin metabolite, indicating that 5-HTP is peripherally metabolized to serotonin, which is then metabolized.</p>				
HVA/DOPAC Ratio	8.47			2.6-8.3
<p>An elevated HVA/DOPAC ratio is associated with excessive supplementation of S-adenosyl methionine, methyltetrahydrofolate, methylcobalamin. A lowered HVA/DOPAC ratio is associated with a genetic deficiency of catechol-O-methyltransferase and/or a nutritional deficiency of S-adenosyl methionine.</p>				
Quinolinic acid/5-HIAA Ratio	0.12			0.32-1.1
<p>A low Quinolinic Acid/5-HIAA ratio is non-significant because quinolinic acid acts as a central nervous system toxin. Thus, if Quinolinic Acid and 5-HIAA are within normal limits, a low ratio is non-significant.</p>				




What Markers Are Included on the Neurotransmitter Test?


Excitatory Neurotransmitters	Inhibitory Neurotransmitters	Other Neurotransmitters
<ul style="list-style-type: none"> Glutamate Histamine PEA Dopamine DOPAC HVA Norepinephrine Epinephrine VMA Acetylcholine Aspartate Oxytocin 	<ul style="list-style-type: none"> Serotonin 5-HIAA GABA Glycine Taurine 	<ul style="list-style-type: none"> 3-Methoxytyramine Metanephrine Tryptamine Tyrosine Tyramine Serine 5-HTP L-DOPA Xanthurenic acid Quinolinic acid Kynurenic acid
Ratios	Diurnal Epinephrine	Diurnal Norepinephrine
<ul style="list-style-type: none"> Norepinephrine/Epinephrine HVA/VMA HVA/DOPAC Quinolinic acid/5-HIAA 	<ul style="list-style-type: none"> Epinephrine (1st Morning) Epinephrine (2nd Morning) Epinephrine (Evening) Epinephrine (Night) 	<ul style="list-style-type: none"> Norepinephrine (1st Morning) Norepinephrine (2nd Morning) Norepinephrine (Evening) Norepinephrine (Night)

Neurotransmitter Cascade


Glutamate/GABA, Glycine & Histamine



Serotonin & Metabolites



Catecholamines & Metabolites



Neurotransmitters & Metabolites

- HVA: Homovanillic acid
- NMN: Normetanephrine
- PEA: Phenethylamine
- VMA: Vanillylmandelic acid
- 5-HIAA: 5-hydroxyindole 3-acetic acid

Cofactors

- BH4: Tetrahydrobiopterine
- Cu: Copper
- Fe: Iron
- Mg: Magnesium
- Mn: Manganese
- MTHF: Methyltetrahydrofolate
- SAMe: S-adenosyl methionine

Enzymes

- AADC: Aromatic L-amino acid decarboxylase
- AANMT: Arylalkylamine N-methyltransferase
- AD: Aldehyde dehydrogenase
- AR: Aldehyde reductase
- COMT: Catechol-O-methyltransferase
- DBH: Dopamine beta hydroxylase
- GA: Glutaminase
- GAD: Glutamate decarboxylase
- GS: Glutamine synthetase
- HDC: Histidine decarboxylase
- HIOMT: Hydroxyindole-O-methyltransferase
- HNMT: Histamine N-methyltransferase
- MAO: Monoamine oxidase
- M6H: Melatonin 6 hydroxylase
- M6ST: Melatonin 6 sulfotransferase
- PHEH: Phenylalanine hydroxylase
- PNMT: Phenylethanolamine N-methyltransferase
- SERHMT: Serine hydroxymethyltransferase
- THRA: Threonine aldolase
- TRPH: Tryptophan hydroxylase
- TYRH: Tyrosine hydroxylase