



VibrantWellness

website: www.vibrant-wellness.com



Urine Hormones Interpretive Guide

1-866-364-0963 | info@vibrant-wellness.com
3521 Leonard Ct, Santa Clara, CA 95054

Table of Contents

Regulatory Statement	04
Introduction to Vibrant America's Urine Hormone Test	04
What are Hormone Metabolites?	04
Lab Methodology of the Urine Hormones Test	05
What Analytes are Included in the Urine Hormone Test?	05
The Steroid Hormone Cascade	06
The Estrogens	07
Estradiol (E2)	07
Treatment Strategies for Increased Estrogens	08
Treatment Strategies for Decreased Estrogens	10
Estrone (E1)	12
Estriol (E3)	13
E3/(E1+E2) Ratio	14
The Estrogen Metabolites	14
Phase I Estrogen Metabolism	15
The 2-OH Metabolites: 2-OH Estradiol and 2-OH Estrone	15
The 4-OH Metabolites: 4-OH Estradiol and 4-OH Estrone	17
16 α -OH Estrone – 16 α -OH-E1	19
2-OH (E1 + E2)/16 α -OH E1 Ratio	20
Phase 2 Estrogen Metabolism	21
The 2 Methoxyestrogens - 2-MeO Estradiol, 2-MeO Estrone	21
The 4 Methoxyestrogens - 4-MeO Estradiol and 4-MeO Estrone	22
Ratio of 4-MeO E1/4-OH E1 and 4-MeO E2/4-OH E2	23
Additional Estrogen Related Markers	24
Bisphenol A	24
Total Estrogen	24
Androgens	25
DHEA/ DHEA-S	26
Androstenedione (A4)	27
The Metabolites of Androstenedione: Androsterone and Etiocholanolone	28
Testosterone (T)	30
Treatment Considerations for High Testosterone and Low Testosterone	32
Epi-Testosterone (Epi-T)	35
T/Epi-T Ratio (Testosterone/Epi-Testosterone Ratio)	36
5 α -DHT (5 α -Dihydrotestosterone)	37
DHT Metabolites: 5 α ,3 α -Androstanediol and 5 β -Androstanediol	37
5 α ,3 α -Androstanediol.	38
5 β -Androstanediol	38

Progestogens-----	38
The Progesterone Metabolites -----	39
Beta(β)-Pregnanediol -----	39
Alpha-Pregnanediol (5 α -Pregnanediol) -----	40
Allopregnanolone -----	41
3 α -Dihydroprogesterone and 20 α -Dihydroprogesterone -----	42
β -Pregnanediol/E2 (β -pregnanediol/estradiol ratio) -----	42
Glucocorticoids and Glucocorticoid Metabolites -----	43
Total Cortisol -----	45
Which is More Accurate: Cortisol or Cortisone? -----	46
Cortisol/Cortisone Ratio -----	46
Cortisol and Cortisone Metabolites -----	47
Cortisol Metabolites: β -Tetrahydrocortisol (β -THF or THF) and α -Tetrahydrocortisol (α -THF or Allo-THF) -----	47
Cortisone Metabolite: β -Tetrahydrocortisone (β -THE) -----	48
Metabolized Cortisol (THF+THE) -----	48
The Mineralocorticoids: Deoxycorticosterone (11-Deoxycorticosterone) and Corticosterone-----	48
Deoxycorticosterone (11-Deoxycorticosterone) -----	48
Corticosterone -----	48
Oxidative Stress -----	49
8-OHDG -----	49
Diurnal Values of Free Cortisol and Cortisone -----	50
Urinary Free Cortisol -----	50
Urinary Free Cortisone -----	51
Interpretation of Cortisol Circadian Rhythm – HPA Axis Dysregulation -----	52
HPA Axis Function -----	52
Lifestyle Treatment Considerations for Glucocorticoids -----	54
Hypercortisolism – Treatment Objectives -----	56
Hypocortisolism – Treatment Objectives -----	56
Urinary Free Melatonin - 6-Hydroxymelatonin Sulfate (MT6s) -----	57
Urinary Creatinine -----	58
Disclaimer -----	59
References -----	60

Regulatory Statement

The content within this interpretive guide is not intended to be a stand-alone medical reference guide, nor is it intended to be a substitute for medical advice from a healthcare provider. The clients who receive Vibrant Wellness Urine Hormone test results are advised to consult their physician, and/or health care provider team for diagnosis and further follow up care.

This Urinary Hormones Interpretive Guide is intended to be used in tandem with Vibrant Wellness's Urine Hormone Test and this guide is provided to users pursuant to the Terms of Use Agreement (the "Terms") on its website www.vibrant-wellness.com.

Vibrant testing does not demonstrate absolute positive and negative predictive values for any disease state or condition. Its clinical utility has not been fully established. Vibrant validates the accuracy and precision of the testing but not of its clinical or diagnostic value. The Urine Hormone Test is advised to be used for wellness and informational purposes only.

Introduction to Vibrant America's Urine Hormone Test

Vibrant America's Urine Hormone test offers a comprehensive evaluation of urinary levels of steroid hormones along with the downstream compounds resulting from metabolism, known as urinary hormone metabolites. An analysis of hormone metabolites along with parent hormones can provide valuable insight into situational factors that influence hormone pathways such as nutrition and other lifestyle factors, genetics, environmental toxicities, and inflammation.

Vibrant's Urine Hormone test also offers adrenal hormone testing at 4 points throughout the day to provide diurnal values of free cortisol as well as free cortisone. Cortisol metabolites are also measured which provide a view of systemic cortisol production, metabolism, and related enzyme functions.

Furthermore, the Urine Hormone test provides diurnal 4-point melatonin values through measurement of the MT6s metabolite. Waking to bedtime analysis of melatonin output can be helpful for assessment of circadian rhythm balance as well as melatonin replacement.

Finally, the Urine Hormone test provides two additional tests, 8-OHDG, a marker of oxidative DNA damage and Bisphenol A (BPA) an estrogenic environmental toxicant. Both factors can help further inform detoxification needs and cancer prevention strategies.

What are Hormone Metabolites?

The process of metabolism of steroid hormones generates intermediate and conjugated compounds with active, systemic effects. Conjugated hormones are forms of hormones linked to glucuronide, sulfate, or glutathione and these are important for bioavailability, storage, cellular transport, and excretion.¹ Serum testing typically measures unconjugated forms of hormones, while urine testing captures conjugated forms which are deconjugated in the lab and free portions are added to free hormone totals.

Metabolites are products of various stages of hormone metabolism. Some metabolites are linked to increased harms such as carcinogenicity or excess hormonal effects. Other metabolites are linked to potentially protective benefits. Evaluation of both the enzyme pathways and metabolites provides opportunities for greater assessment of health risks as well as therapeutic lifestyle interventions.


Lab Methodology of the Urine Hormones Test

Vibrant America utilizes liquid chromatography-tandem mass spectrometry (LC-MS/MS) to measure parent steroid hormones (Estradiol, Estrone, Estriol, DHEAS, Testosterone, Cortisol, Cortisone, etc.) All metabolites are measured using gas chromatography/mass-spectrometry (GC-MS). These technologies allow for simultaneous measurements of multiple analytes and have a low limit of quantification which means precision and accuracy despite low levels of hormone in the sample.

What Analytes are Included in the Urine Hormone Test?

Estrogens	Estradiol Estrone Estriol 2-OH Estradiol 2-OH Estrone 4-OH Estradiol 4-OH Estrone	4-OH Estrone 16 α -OH Estrone 2-MeO Estradiol 2-MeO Estrone 4-MeO Estradiol 4-MeO Estrone Bisphenol A
Androgens	DHEA Testosterone Epi-Testosterone Androstenedione Androsterone	Etiocholanolone 5 α -DHT 5 α ,3 α -Androstanediol 5 α -Androstanediol DHEA-S
Progestogens	Progesterone Allopregnanolone Allopregnanediol 3 α -Dihydroprogesterone	20 α -Dihydroprogesterone β -Pregnanediol α -Pregnanediol
Glucocorticoids	Cortisol Cortisone β -Tetrahydrocortisol (β -THF)	α -Tetrahydrocortisol (α -THF) β -Tetrahydrocortisone (β -THE) Corticosterone
Other	Melatonin 8-OHdG	

The Steroid Hormone Cascade




- Androgens
- Estrogens
- Glucocorticoids
- Mineralocorticoids
- Progestogens

Enzyme Abbreviations

5α-R	5α-Reductase	11β-HSD	11β-Hydroxysteroid dehydrogenase
5β-R	5β-Reductase	17α-HSD	17α-Hydroxysteroid dehydrogenase
11β-OH	11β-Hydroxylase	17β-HSD	17β-Hydroxysteroid dehydrogenase
17α-OH	17α-Hydroxylase	20α-HSD	20α-Hydroxysteroid dehydrogenase
17,20-Lyase	Same enzyme as 17α-OH	AR	Aromatase
21-OH	21-Hydroxylase	AS	Aldosterone Synthase
3α-HSD	3α-Hydroxysteroid dehydrogenase	CYP	Cytochrome p450 (scc, 1A1, 1B1 & 3A4)
3β-HSD	3β-Hydroxysteroid dehydrogenase	COMT	Catechol-O-Methyl-Transferase

The Estrogens

Estrogens		
	Estradiol	4-OH Estrone
	Estrone	16 α -OH Estrone
	Estriol	2-MeO Estradiol
	2-OH Estradiol	2-MeO Estrone
	2-OH Estrone	4-MeO Estradiol
	4-OH Estradiol	4-MeO Estrone
	4-OH Estrone	Bisphenol A




The Parent Estrogens: Estradiol (E2) and Estrone (E1)

Estradiol (E2)

Estradiol, commonly known as E2 or 17 β -Estradiol, is the predominant and most biologically active estrogen in circulation in males and females. Estradiol is also the strongest estrogen and is 80 times more potent than estriol, the weakest estrogen.^{2,3} Estradiol plays a key role in the development of the female reproductive system, including breast development, and has non-reproductive roles in cognition and neuroprotection, lipid and glucose homeostasis, adipose distribution, cardiovascular health, pancreatic cell function, bone maintenance and wound healing.⁴ In males, E2 is predominantly produced by the testes.⁴ Males and postmenopausal women have significantly lower levels of estradiol than premenopausal women.

E2 originates from cholesterol, which is converted through a cascade of progesterone and androgen intermediates, finally resulting in estradiol through aromatization from testosterone.⁵ In pre-menopausal women this occurs primarily in the ovary. In postmenopausal women, estradiol is sourced directly from the peripheral tissues, predominantly from estrone, sourced from adrenal precursors, and converted to estradiol via 17 β -HSD.^{3,4}



Elevated E2

Known causes of increased estradiol as well as increased aromatase activity, which leads to increased estradiol, are listed in Table 1.0 below.

↑ **Table 1. Causes of Increased E2**

Known Causes of ↑ Estradiol	Known Causes of ↑ Aromatase Activity (which lead to ↑ E2)
Excess Aromatase Activity	Inflammation ⁶
Chronic Liver Disease ⁷	Age
Hyperthyroidism	Obesity
Dysregulated Gut Microbiome/Elevated Beta-Glucuronidase ⁶	Insulin
Cancer Tissue	Gonadotropin Use
Exogenous Estrogen Use	Alcohol
	Rare Genetic Causes /CYP19A1 Rearrangement ⁸

Conditions associated with prolonged high estradiol (E2)^{4,5}

- PCOS
- Uterine Cancer
- Gastric Cancer
- Multiple Sclerosis
- Obesity
- Endometriosis
- Prostate Cancer
- Thyroid Cancer/Goiter
- Gynecomastia
- Breast Cancer
- Pituitary Cancer
- Systemic Lupus Erythematosus (SLE)
- Oligospermia
- Male Hypogonadism
- Breast Cysts
- Gallbladder Disease
- Fibroids
- Ovarian Cancer

Treatment Strategies for Increased Estrogens

Lifestyle interventions such as diet, supplements and detoxification strategies can improve overall estrogen metabolism and clearance from the body through several mechanisms.

- Reduce active estradiol, estrone and xenoestrogens through reduced exposure, increased detoxification and increased sex hormone binding globulin activity
- Target microbiome support for enhanced estrogen clearance and reduction of beta-glucuronidase
- Decrease aromatization from androgens to estrogens, especially with excess adiposity

Table 2. Lifestyle Interventions for Estrogen Reduction

Lifestyle Considerations	Reduction of Estrogens and Xenoestrogens	Support Microbiome for Estrogen Clearance	Decrease Aromatase Activity
Specific Supplements	<p>DIM - up to 300mg per day⁹</p> <p>Sulforaphane^{10*}</p> <p>Vitamin D Repletion in deficient women (to blood levels above 32ng/mL) ↑ SHBG and ↓ bioavailable E2 and testosterone¹¹</p>	<p>Lactobacillus spp^{12*}</p> <p>Calcium D-glucarate 1500mg or higher/day*</p>	<p>Numerous flavonoids show preclinical evidence of decreased aromatase activity: *</p> <p>Resveratrol¹³</p> <p>Grape Seed extract¹⁴</p> <p>Citrus Peels/citrus flavonoids¹⁵</p> <p>Apigenin¹³</p> <p>Chrysin¹⁶</p>
Diet	<p>Brassica family vegetables – 1.5 servings/day¹¹</p> <p>Wheat sprout juice -100ml/day¹⁷</p> <p>Broccoli sprouts^{10*}</p>	<p>Fiber from fruit/veggies - 30 g/day¹⁸</p> <p>Flaxseed meal – 10 -30 g/day¹⁹, positive benefits with 50mg lignans per day</p> <p>Kefir^{12*}</p>	<p>White button mushrooms^{13*}</p> <p>Foods with apigenin and/or resveratrol: grapes, berries, celery, parsley, onions, oranges, chamomile, thyme, basil, oregano, red wine²⁰</p> <p>Kale, collards</p>
Exercise/Sauna	<p>Dry Sauna- 15 minutes 3x per week*</p> <p>Exercise –100-300minutes/week treadmill/aerobic reduces estrogen²¹</p>	<p>Moderate exercise benefits microbiome diversity²²</p>	<p>150 minutes moderate activity per week if weight loss is needed</p>
Misc. Lifestyle Factors	<p>Reduce or avoid exogenous HRT</p> <p>Avoid endocrine disruptors such as phthalates, BPA, growth hormones in meat, pesticides²³</p>		<p>Weight loss if obesity or PCOS is present²⁴</p> <p>Lifestyle support to reduce inflammation secondary to obesity²⁵</p>

*Dosage information for this indication is not well established, refer to empiric guidelines of safe and effective use

** Assess risk/benefit analysis for the specific individual prior to giving hormone replacement therapy

Further testing for further treatment avenues for consideration for estrogen excess:

1. **NutriPro or Methylation Panel** by Vibrant Wellness testing for genetic SNPs which may lead to estrogen metabolism issues (MTHFR, COMT)
2. **Gut Zoomer** test for beta glucuronidase, estrobolome characteristics
3. **Saliva Hormone** test for bioavailable hormone levels
4. **Environmental Toxins, Mycotoxins and/or Heavy Metals test-** for plasticizers, mold, metals and other environmental toxicant burden
5. **Vibrant America Cardiac Health and Diabetes Panel** for assessment of cardiometabolic function
6. **Neural Zoomer Plus or Neurotransmitter** test – progesterone is a neuroprotectant and neuromodulator, progesterone deficiency may affect these pathways

Decreased E2

Decreased levels of estradiol are most commonly caused by aging or ovariectomy, however there can be many causes from other mechanisms, such as decreased aromatase activity and increased prolactin levels, as outlined in Table 3.

↓ **Table 3. Causes of Decreased E2**

Systemic Causes of ↓ Estradiol	Causes of ↓ Aromatase Activity (which lead to ↓ E2)	Causes of ↓ Prolactin (which lead to ↓ E2)
Menopause/Aging	Aromatase Inhibitors	Antidepressants/Antipsychotics
Ovariectomy	Herbicides	Anticonvulsants
Compounds Which Cause Aromatase Inhibition	Agricultural Antifungals	Opiates
Compounds Which Increase Prolactin	Immunosuppressive Drugs	Hypothyroidism ²⁶
Hypogonadism	Antimalarials	Estrogens
Genetic Disorders/17α-hydroxylase/17,20-lyase deficiency ⁴	Cigarette Smoke ⁵	Anti-Androgens
	Anti-Mullerian Hormone*	Anti-Hypertensives
		H2 Receptor Antagonists

*A hormone produced by reproductive tissues, not currently measured at Vibrant Wellness

Conditions associated with prolonged low estradiol (E2)

- Alzheimer’s Disease
- Osteoporosis
- Osteoarthritis
- Parkinson’s Disease
- Diabetes Mellitus
- Eclampsia
- Colorectal Cancer

Treatment Strategies for Decreased Estrogens

Treatment Objectives - Consider which of the following single or combination of therapeutic approaches is the best fit for the patient with estrogen deficiency symptoms.

- Modulate and reduce vasomotor symptoms, insomnia, and mood disturbances
- Enhance quality of life during natural transitions in aging
- Prevent sequelae of severe estrogen deficiency conditions i.e., cognitive, bone, and cardiometabolic

Table 4. Lifestyle Interventions for Estrogen Deficiency

Lifestyle Considerations	Vasomotor Symptom and Mood Support
Specific Supplements	<p>Black Cohosh for mood and hot flashes 40-127mg/day²⁷</p> <p>Soy Isoflavones for vasomotor symptoms, bone loss, hypertension – trial doses ranged from 40-120mg isoflavones per day²⁸ Use organic, non-GMO isoflavones</p> <p>Maca 3.0 g/d for libido, depression in menopausal women; improvement in testosterone, but not estradiol, in women²⁹</p> <p>Kudzu (Pueraria spp.)- may aid vasomotor symptoms, bone loss, vaginal symptoms – 20 mg to 2.5g/day in most human clinical studies for menopausal symptoms. Topical doses 0.5–1.0 g/day in gel³⁰⁻³¹</p> <p>Red Clover for hot flashes $\geq 80\text{mg/d}$³²</p> <p>Red Ginseng 0.9-3mg/day for fatigue, memory, hot flashes³³</p> <p>Ashwagandha 300mg 2x/day for vasomotor symptoms, increased estradiol in perimenopausal women in one study³⁴</p> <p>Vitamin E 800 IU/day for hot flashes³⁵</p> <p>Fenugreek seed 600mg/day de-husked seed extract for vasomotor symptoms³⁶</p>
Diet	<p>Mediterranean diet to reduce menopausal symptoms³⁷ & cardiometabolic dysfunction³⁸</p> <p>Diet with increased fruits/vegetables, fiber & less sugar and fat improves vasomotor symptoms³⁷</p> <p>High fiber and soy isoflavones decrease vasomotor symptoms³⁷</p>
Exercise/Sauna	<p>150 minutes moderate aerobic + strength training³⁹ - support for bone health and general estrogen deficiency symptoms</p> <p>Yoga variable practices aid vasomotor and psychological health⁴⁰</p>
Stress Management	<p>Cognitive Behavioral Therapy for vasomotor symptoms and insomnia⁴¹</p> <p>Hypnosis for reduction in hot flash severity⁴²</p> <p>Aromatherapy/massage⁴²</p> <p>Mindfulness Meditation for distress from vasomotor symptoms⁴³</p> <p>Acupuncture for hot flashes⁴⁴</p>

Miscellaneous Support	<p>Nutrient and Bone Support via diet, supplements and exercise as indicated for osteoporosis prevention</p> <p>Cardiometabolic support via diet, exercise and supplements as indicated</p> <p>ERT or HRT - consider for severe symptoms or prevention of severe estrogen deficiency sequelae**</p>
-----------------------	--

* Use of some supplements could result in specific hormone increases with use. Always research a specific product prior to use with patients, especially those at higher risk for hormone related cancer.


** Fully research the individual patient's risk/benefit analysis regarding hormone related cancer and other side effects or ERT or HRT prior to recommending.

Further testing considerations for estrogen deficiency:

1. **Cardiac Health Panel** and **CardiaX tests**
2. **Salivary Hormones** test for bioavailable hormone levels
3. **APOE** test for Alzheimer's Disease
4. **Colorectal cancer screening (Not available at Vibrant)**
5. **DEXA** testing for osteoporosis

Estrone (E1)

Estrone (E1) is a weaker estrogen than estradiol (E2); by some estimates estrone has less than 10% the strength of estradiol.⁴⁵ While levels of estrone do not differ significantly in pre- vs. post-menopausal women,⁸ estrone is the predominant estrogen in post-menopausal women⁴⁶ by a factor of 100-fold compared to estradiol.⁴⁵ In post-menopausal women, estrone concentrations have been positively correlated with bone mineral density and breast cancer risk, and inversely correlated with colon cancer risk. Serum estrone levels are an important indicator of serum estradiol levels in post-menopausal women.⁴⁵




↑ Causes of Increased Estrone (E1):

- Hormone replacement therapy (HRT)
 - HRT can increase estrone, even if only estradiol is given, as estrone and estradiol can interconvert.⁴⁷
 - Advanced age and obesity will further elevate estrone values in users of HRT.⁴⁸
- Liver disease⁴⁹
- Hyperthyroidism⁵⁰
- Hormone producing tumors (ovary, adrenal gland, etc.) can also raise estrone^{51,52}
- Refer to Table 2. in E2 section for lifestyle considerations to reduce excess estrogen if clinically indicated.

↓ Causes of Decreased Estrone (E1):

- E1 levels appear to be similar in pre and postmenopausal women,⁸ however there is wide variation in hormone levels in the normal range.
- Factors such as DHEA and testosterone levels, SHBG levels, genetic predisposition, adiposity, and aromatase activity affect estrone levels.
- Low estrone levels, below reference ranges, often result from decreased aromatase activity.
 - Refer to the Table 3. in E2 section of this guide for substances that can decrease aromatase activity
 - Refer to Table 4. in the E2 section for lifestyle considerations to manage low estrogen symptoms if clinically indicated.

Estriol (E3)



What is Estriol (E3)?

Estriol or E3, originates from cholesterol and androgen intermediates, notably aromatization from testosterone and androstenedione. Estradiol is reversibly oxidized to estrone and most E3 is formed from estrone via CYP3A4 metabolism and 17β-HSD as seen in the image above.⁵³ Both E2 and E1 can be irreversibly converted to estriol in the liver. E3 is the weakest of the three estrogens: it dissociates rapidly from estrogen receptors.⁵⁴ Estriol is the predominant estrogen of pregnancy where it regulates uterine/placental blood flow and placental vascularization. Lab testing for E3 has been most used as a maternal screening for fetal anomalies.⁵⁴ E3 is used off-label for menopausal symptoms such as hot flashes, vaginal atrophy, and bone density. It has also been studied for use of immunomodulation and neuroprotection in multiple sclerosis and protection from atherosclerosis.⁵³

↑ Causes of Increased Estriol (E3):

- Pregnancy
- Use of bioidentical hormone replacement prescriptions⁵⁵
- Conditions that create elevated estradiol may also increase estriol, refer to E2 section for further reference.

↓ Causes of Decreased Estriol (E3):

- During pregnancy, low estriol is a marker of fetal compromise, genetic disorders, fetal adrenal deficiency, fetal growth restriction, congenital adrenal hyperplasia, aromatase deficiency and other fetal conditions.⁵⁶
- Low estradiol conditions can contribute to a low estriol status in non-pregnant women, for further reference refer to E2 section for causes and lifestyle considerations for low estrogen.

E3/(E1+E2) Ratio

What is the Ratio of E3/E1-E2?

The ratio of estriol (E3) to the sum of estrone (E1) + estradiol (E2) is known as the EQ or estrogen quotient. This equation was popularized by Henry Lemon, MD, a breast cancer researcher. The estrogen quotient quantifies the concept that a higher ratio of estriol relative to the two stronger estrogens has value for both prevention and outcomes in breast cancer. From his research with 24-hour urine samples in female participants, he suggested that a ratio value of E3/E1+E2 of 1 or 1.5 is more optimal.⁵⁷

Vibrant America's urine hormone test reference range for the ratio of E3/E1+E2 is based on average values found in samples from healthy women of varying ages. The high and low results from these reference ranges on the urine hormone test are not based on Dr. Lemon's findings in urinary estrogens. Thus, the ratio of E3/E1+E2 is provided for comparative use of estrogens only. Clinical application of the results should be based on the provider's discretion.

Lifestyle Considerations for Low Estriol Ratio - E3/E2+E1


Treatment objectives for a low E3/E2+E1 ratio

- Minimize exposure to exogenous estrogens
- Assess and correct for any pre-existing iodine deficiency
- Refer to Table 3. Lifestyle Interventions for Estrogen Reduction and Metabolism

Integrative medicine physicians have used iodine therapy empirically to balance effects of a low E3/E2+E1 ratio and potentially reduce risk of breast disease. This has been based on combining ideas from human studies regarding estriol risk and fibrocystic breast disease and breast cancer,^{57,58} as well as pre-clinical evidence that iodine/iodide can modulate estrogen metabolic pathways to favor inactivation of E2.^{59,60} Further clinical trials are needed.

Iodine interventions, beyond diet suggestions for deficiency, require physician oversight including measures such as initial thyroid and iodine testing (available at Vibrant), ongoing monitoring, and careful consideration of dose, form, and length of iodine usage. Iodine excess may precipitate hyperthyroidism, hypothyroidism, goiter, and/or thyroid autoimmunity.⁶¹


The Estrogen Metabolites



Phase I Estrogen Metabolism

While estrogen is crucial for female and some aspects of male development; the metabolites of estrogen can have negative as well as positive effects. Phase I metabolism of estrogen is an oxidative process that involves hydroxylation reactions to render the metabolites hydrophilic rather than lipophilic. The specific hydroxylation reactions that occur are either induced or inhibited by the cytochrome P-450 (CYP) enzyme system and affected by compounds that affect that system. These reactions produce estrogen metabolites. Modulation of the estrogen hydroxylation reactions is necessary, as some metabolites can promote disease through increased estrogenic effects and/or carcinogenesis.⁶² Carcinogenicity of estrogens can result from the trifecta of epigenetic effects, direct estrogen receptor effects and metabolite effects. Estrogen metabolite effects cause harm through metabolites binding to estrogen receptors, which cause effects of increased cell proliferation, and production of reactive quinones which cause DNA damage. There are 13 known estrogen metabolites that can be produced from hydroxylation at the 2, 4 or 16 carbon of the steroid ring.⁶³ These hydroxylation pathways of estrogen metabolism are termed the 2-hydroxylation pathway (2-OH), the 4-hydroxylation pathway (4-OH), and the 16-hydroxylation pathway (16-OH).⁶² The 2-hydroxyestrogen metabolites (2-OH) along with the 4-hydroxyestrogen metabolites (4-OH) are known as the catechol estrogens.

The 2-OH Metabolites: 2-OH Estradiol and 2-OH Estrone



Estrogen metabolism produces a larger amount of 2-OH and 16-OH metabolites and a lesser number of 4-OH metabolites. 2-OH metabolites are regulated by CYP1A1 enzymes.⁶² 2-OH is generally considered a safer metabolic pathway regarding the potential for harmful intermediate production. There are a few reasons for this theory. First, 2-OH metabolites have a significantly reduced binding affinity for the estrogen receptor. Furthermore, some early evidence indicates that 2-OH metabolites may act as tumor suppressors to inhibit breast cancer cell growth and proliferation.⁶⁴ Finally, increased 2-hydroxylation of parent estrogens, rather than 4- or 16-hydroxylation, is associated with reduced risk of breast cancer.⁶⁵ It is, however, known that the 2-OH pathway metabolites can form increased reactive oxygen species and can also result in the formation of a low level of depurinating adducts.^{66,67} Thus, though considered safer, 2-OH metabolites should be kept in balance. Prompt methylation of these intermediates in phase II metabolism via the COMT enzyme pathway is the optimal detoxification step.

Irregular levels of 2-OH intermediates suggest irregular levels of parent estrogens (E1, E2), and/or increased or decreased activity of enzyme pathways (CYP1A1 or COMT) which precede or follow 2-OH production.


↑ Potential Causes of Increased 2-OH Estradiol and 2-OH Estrone

- Increased endogenous hormone production or exogenous hormone exposure
- Increased CYP1A1 enzyme activity relative to CYP3A4 or CYP1B1 enzyme pathways
- Decreased COMT activity and decreased production of 2-MeO metabolites
 - Consider COMT support if there are increased 2-OH metabolites and sub-optimal levels of 2-MeO metabolites
 - Support for COMT function can be gained through optimization of nutrient cofactors and methyl donors, such as increasing intake of foods high in methionine, methyl vitamin B12 (methylcobalamin), pyridoxal-5-phosphate (P5P), betaine, choline, folate (5-methyltetrahydrofolate), and magnesium.⁶⁸
- Conditions that create elevated estrogens may also increase 2-OH metabolites
 - Refer to Table 2.0 in E2 section for lifestyle considerations to reduce excess estrogen if clinically indicated.

↓ Potential Causes of Decreased 2-OH Estradiol and 2-OH Estrone

- Decreased endogenous levels of estradiol and/or estrone
- Decreased induction of CYP1A1 enzyme activity
- Increased induction of alternate pathways of CYP3A4 or CYP1B1 enzymes
- Increased COMT activity and increased production of 2-MeO metabolites
- Are 2-OH metabolites decreased along with decreased E1 and/or E2 and estrogen deficiency symptoms?
 - Conditions that create decreased overall estrogen may also decrease 2-OH metabolites.
 - Refer to Table 4.0 in the E2 section for lifestyle considerations to manage low estrogen symptoms if clinically indicated.
- Are 2-OH metabolites decreased in combination with elevated E1 or E2, and/or increased 4-OH or 16-OH values?
 - Consider the following strategies to promote optimal estrogen metabolism and increase 2-OH pathway metabolism:
 - Refer to Table 2. in E2 section for lifestyle considerations to reduce overall estrogen burden
 - Interventions that have been studied to increase 2-OH metabolism include:
 - o Cruciferous vegetable compounds such as Indole-3-carbinol, DIM (diindolylmethane), and sulforaphane⁶⁹
 - o Flax⁷⁰
 - o Isoflavones (Soy, Kudzu, Clover)
 - o Resveratrol⁷¹
 - o Omega 3 Fatty Acids
 - o Coffee⁷²
 - o Weight Loss⁷³
 - o Moderate Exercise
 - o Rosemary⁷⁴

The 4-OH Metabolites: 4-OH Estradiol and 4-OH Estrone



4-OH metabolites are the product of estrogen metabolism governed by the CYP1B1 enzyme pathway. 4-OH metabolites are considered an unfavorable estrogen intermediate. 4-OH metabolites are easily oxidized to quinone and semi-quinone conjugates, which can react with DNA and lead to cancer initiating mutations.⁷⁵ Elevated levels of depurinating estrogen-DNA adducts have been found in women with breast, ovarian and thyroid cancer. Men with prostate cancer or non-Hodgkin lymphoma also have high levels of estrogen-DNA adducts.⁷⁶

In a recent comprehensive review, urine 4-OH metabolites were found to be significantly associated with breast cancer risk in premenopausal women. Furthermore, pre-clinical evidence links increased expression of CYP1B1, which produces 4-OH metabolites, with increased breast cancer risk.⁷⁷

↑ Potential Causes of Increased 4-OH Estradiol and 4-OH Estrone

- Elevated 4-OH-E1 can occur more commonly in pre/peri-menopausal women transitioning to menopause.
- Increased endogenous hormone production or exogenous hormone exposure
 - Conditions that create elevated estrogens may also increase 4-OH metabolites
 - Refer to Table 2.0 in E2 section for lifestyle considerations to reduce excess estrogen if clinically indicated.
- Downregulation of CYP 3A4 or CYP1A1 pathways of enzyme metabolism
- Decreased COMT activity and/or decreased methylation to 4-MeO metabolites
 - Consider lifestyle support to increase COMT activity
 - Support for COMT function can be gained through optimization of nutrient cofactors and methyl donors, such as increasing intake of foods high in methionine, methyl vitamin B12 (methylcobalamin), pyridoxal-5-phosphate (P5P), betaine, choline, folate (5-methyltetrahydrofolate), and magnesium.⁷⁸
- As 4-OH metabolites can react with DNA and initiate carcinogenic mutations, consider lifestyle support for harm mitigation from quinone conjugates/DNA adducts.
 - The substances below boost detoxification of quinone conjugates via increased activity of NQO1 (NAD(P)H Quinone Dehydrogenase 1 gene) and/or glutathione conjugation
 - o Resveratrol⁷⁶
 - o N-Acetyl Cysteine (NAC)⁷⁶
 - o Glutathione⁷⁶
 - o Sulforaphane⁷⁹

- Upregulation of CYP1B1 activity
 - Refer to CYP1B1 section below
 - o Weight loss, reduced pollutant exposure, and plant compounds such as stilbenes, flavonoids and more, have been shown to reduce expression of CYP1B1.



Potential Causes of Decreased 4-OH Estradiol and 4-OH Estrone

- Decreased estrone or estradiol can result in decreased 4-OH metabolites
 - Conditions that create decreased overall estrogen may also decrease 4-OH metabolites
 - Refer to Table 4.0 in the E2 section for lifestyle considerations to manage low estrogen symptoms if clinically indicated.
- In general, lower levels are preferred based on current understanding of estrogen metabolism pathways.



A Note about CYP1B1

CYP1B1 metabolizes environmental pollutants and plant flavonoids but also endogenous substances such as melatonin, retinoic acid and other fat-soluble vitamins, fatty acids, and steroid hormones.⁸⁰ It plays a role in homeostasis of glucose, weight, redox reactions, and blood pressure and supports eye development. Overexpression of CYP1B1 can have harmful effects such as increased production of genotoxic catechol estrogens, promotion of cancer and potentially other harms.⁸¹ Cancerous tumors overexpress CYP1B1, and research has been done on inhibition of this enzyme as a cancer therapeutic target.

What Causes Increased CYP1B1?

CYP1B1 can be overexpressed by increased leptins/obesity, polyaromatic hydrocarbons, cigarette smoke, inflammatory cytokines, cancerous tissue, UV exposure, biotin supplementation, environmental pollutants, and genetic polymorphisms⁸²⁻⁸⁴

What can inhibit CYP1B1?

The following compounds demonstrate an inhibitory effect on CYP1B1 in pre-clinical studies:

- Stilbene compounds such as:
 - Resveratrol and piceatannol found in grapes, berries, peanuts⁸⁵
- Flavonoids such as:
 - Kaempferol (flavonoid found in foods such as broccoli, apples, strawberries, tea, and beans⁸⁵
 - Quercetin⁸⁵
 - Rutin⁸⁶
 - Apigenin⁸⁶– found in foods such as parsley, celery, basil, chamomile, cilantro, and oregano
 - Naringenin (flavonoid found in citrus, tomatoes, figs)⁸⁵
- Coumarins which are found in cherries, cassia cinnamon, apricots and strawberries
- Anthraquinones which are laxative compounds found in aloe, rhubarb and senna⁸⁷
- Melatonin⁸⁶



What are quinone conjugates and DNA adducts and why are they important?

The catechol estrogens, 2-OHE1 or 2-OHE2 and 4-OHE1 or 4-OHE2 have 2 hydroxyl groups that are generally methylated by the enzyme COMT into 2 and 4 -methoxy forms of estrogens. However, remaining unmethylated hydroxysteroids, if further oxidized can generate compounds called semiquinones which are further oxidized

into quinones. Semiquinones and quinones can become a part of redox cycling processes and generate reactive oxygen species, which can be neutralized, or if not, they can bind to DNA to form stable or depurinating adducts.


The depurinating estrogen adducts from 4-OHE1 or 4-OHE2 are called estrogen (E1 or E2) -3,4-quinones and they are known to become carcinogens. Estrogen-3,4 quinones induce DNA mutations in which errors form in the DNA repair process. Depurinating estrogen adducts can be initiators of many types of cancers, such as breast, ovarian, thyroid, non-Hodgkin lymphoma, and prostate cancer.⁸⁸

The good news is that cells have several protective mechanisms to maintain estrogen metabolism homeostasis and avoidance of carcinogenic damage to DNA. Transformation of catechol estrogens to methyl estrogens through COMT inactivates carcinogenic potential. Also, estrogen-3,4 quinones can be reduced by NQO1 back to catechol estrogens or be conjugated with glutathione to prevent DNA adduct formation. Thus, increasing methylation, glutathione activity and/or NQO1 pathway activity are all ways to decrease carcinogenicity.⁸⁹

Pre-clinical evidence indicates the following compounds can boost detoxification of quinone conjugates via increased activity of NQO1 and/or glutathione conjugation:

- Resveratrol⁸⁸
- N-Acetyl Cysteine (NAC)⁸⁸
- Glutathione⁸⁸
- Sulforaphane⁸⁹

16 α -OH Estrone – 16 α -OH-E1



The Phase I metabolism of estrone via hydroxylation at C-16 yields the product of 16 α -OH-E1. Reduction of 16 α -OH estrone through 17 β -HSD1 enzyme activity produces estriol. Both 16 α -OH-E1 and estriol are considered estrogenic. 16 α -OH-E1 irreversibly and strongly binds to the estrogen receptor, which promotes estrogenic effects. It has been shown to have a mitogenic, growth enhancing, effect in breast cancer cells.⁹⁰ Pre-clinical evidence shows that 16 α -OH-E1 can cause the formation of depurinating DNA adducts.⁹¹ Similar to estrogen, 16-OH E1 appears help to protect against osteoporosis.⁹²

An epidemiologic study of premenopausal women showed that an increased ratio of 16 α -OH pathway metabolites to parent estrogens was associated with increased risk of breast cancer.⁹³ Contrary to this, more recently in postmenopausal women, higher 16 α -OH levels, relative to estrone and estradiol, were found to be associated with a non-statistically significant, decreased risk of breast cancer.⁹⁴ Varying lab technologies may have confounded the issue and further research is warranted.

↑ Causes of Increased 16a-OH Estrone

- Increased Estrone production
- Increased CYP3A4 metabolism
- Reduced activity of 17B-HSD1 which converts 16-OH to estrinol
- Refer to Table 2.0 in E2 section for lifestyle considerations to reduce excess estrogen if clinically indicated
- Enhance estrogen metabolism via 2-OH pathways if clinically indicated

↓ Causes of Decreased 16-OH Estrone

- Decreased levels of estrogen and/or estrone

2-OH (E1 + E2)/16a-OH E1 Ratio

This test measures the ratio of 2-hydroxy estrogens compared to 16 hydroxy estrone. Given that 2-OH estrogens are believed to be less estrogenic than 16 hydroxyestrogens, it follows that a higher 2-OH ratio would be preferred when considering the broad category of estrogen related diseases. Numerous studies have shown that a low ratio of 2-OH/16a-OH is an independent risk factor for breast cancer.⁹⁵ Numerous studies have also shown that diet and lifestyle can lower this ratio.

↓ Causes of Decreased 2-OH/16-OH Ratio

- Breast cancer risk factors
- Estrogen excess related diseases
- Alcohol⁹⁵
- Obesity and increased BMI⁹⁵
- Consider lifestyle interventions to increase 2-OH pathway metabolism and increase the 2-OH/16OHE1 ratio, if clinically indicated, such as:
 - o Cruciferous vegetables^{96,97}
 - Indole-3-Carbinol
 - DIM – di-indolylmethane (studies reference doses of 100-400mg/day)
 - Sulforaphane
 - o Flax⁹⁸ (study dose: 25grams of flax per day)
 - o Isoflavones⁹⁹ (study dose: 107mg soy isoflavones/day)
 - o Resveratrol¹⁰⁰ (study dose: 1gram/day)
 - o Coffee¹⁰¹ (study dose: ≥4 cups/day)
 - o Weight Loss¹⁰²
 - o Rosemary¹⁰³

↑ Causes of Increased 2-OH/16-OH Ratio


- Generally considered favorable unless estrogen deficiency related diseases, such as osteoporosis, are clinically relevant

Phase 2 Estrogen Metabolism

Metabolism of catechol estrogens in Phase 2 metabolism follows the hydroxylation step in Phase 1. Biochemical processes of methylation, sulfation, glucuronidation and glutathione conjugation help to biotransform estrogen metabolites to inactive states and into hydrophilic compounds in preparation for excretion. The actions of enzymes COMT, GST and NQO1 are notably protective in detoxification of catechol estrogens as well as protection from carcinogenic adducts.¹⁰⁴ COMT is present in many human tissues, notably the liver and breast, but also the kidney and red blood cells. Genetic polymorphisms in the COMT gene can result in up to a 4-fold decrease in enzyme activity. This is meaningful as it is hypothesized that there is an increased risk for breast cancer in individuals with low COMT activity, due to an increase in accumulation of harmful estrogen intermediates.¹⁰⁵

Phase 2 Estrogen Metabolism

The 2 Methoxyestrogens - 2-MeO Estradiol, 2-MeO Estrone



Following CYP1A1 conversion to 2-OH metabolites, the next metabolic step for these intermediates is Phase 2 detoxification and methylation via the COMT enzyme pathway. 2-OH metabolites are converted into the 2-methoxyestrogens, 2-MeO E2 and 2-MeO E1. 2-MeO E2 appears to have a low affinity for the estrogen receptor and has shown anti-cancer effects in preclinical studies through promotion of apoptosis and anti-angiogenesis.^{106,107} These studies have generated much interest in this molecule as a potential anti-cancer agent. 2-MeO E2 also has demonstrated protective cardiovascular effects and has been shown to have a negative feedback inhibition of CYP1B1 activity.¹⁰⁸

↑ Causes of Increased 2-MeO Metabolites

- High levels of parent estrogens and/or 2-OH metabolites
- Increased COMT activity
- Generally considered favorable, however can be a signal to reduce estrogen therapy if consistent with additional estrogen elevations (E1, E2, 2-OH metabolites)

↓ Causes of Decreased 2-MeO Metabolites

- Low levels of parent estrogens or 2-OH metabolites
- Reduced COMT activity/methylation
 - Consider addition of methyl donors (See COMT Support section)




COMT Support

Catechol-O-methyltransferase (COMT) is an important phase II enzyme that inactivates the harmful potential of catechol estrogens through the methylation and formation of methoxy estrogens. Support for COMT function can be gained through optimization of nutrient cofactors and methyl donors, such as increasing intake of foods high in methionine, methyl vitamin B12 (methylcobalamin), pyridoxal-5-phosphate (P5P), betaine, choline, folate (5-methyltetrahydrofolate), and magnesium.¹⁰⁹ Also notable is that pre-clinical research suggests that a high sucrose diet may inhibit COMT activity.¹¹⁰

Table 5. Food Sources of Nutrients to Aid Methylation and COMT Enzyme Activity

Nutrients	Food Sources
Methionine	Turkey, pork, chicken, beef, fish, egg, Brazil nuts, soybeans, sesame and pumpkin seeds, teff, oats, spirulina ¹¹¹
Vitamin B12	Meats and meat products (especially liver and kidney), poultry, fish, shellfish, eggs ¹¹¹
Vitamin B6	Poultry, pork, beef, lentils, sunflower seeds, amaranth, garlic, prunes ¹¹¹
Betaine	Quinoa, beets, spinach, whole grains (especially rye, kamut, bulgur, amaranth, barley, and oats) sweet potato, meats, poultry ¹¹¹
Choline	Poultry, soy, egg, salmon, quinoa, wheat germ, milk, cauliflower, peas, tofu, broccoli ¹¹²
Folate	Spinach, legumes, liver, sunflower seeds, quinoa, asparagus, avocados, mustard greens, artichokes ¹¹¹
Magnesium	Nuts, seeds, legumes, whole grains ¹¹¹

The 4 Methoxyestrogens - 4-MeO Estradiol and 4-MeO Estrone



As indicated previously, 4-OH metabolites from CYP1B1 enzyme activity are considered an unfavorable estrogen intermediate due to potential for quinone conjugates which can cause DNA damage and promote carcinogenicity. The solution for inactivation of 4-OH metabolites is rapid methylation of these products in Phase 2 with the formation of the 4-methoxyestrogens, 4 Me-O estradiol and 4-Me-O estrone. The methylation of 4-OH intermediates is governed by COMT and requires substrates such as Vitamins B6, B12, folate and betaine. COMT activity can be affected by genetic polymorphisms and nutrient status and thus can vary greatly among individuals. 4-Methoxyestrogens have been shown to have neutralizing activity on the harmful effects of 4-OH metabolites via prevention of oxidative damage.¹¹³ There have also been other benefits shown by 4-MeO estrogens. In a pre-clinical study, 4-MEO's showed a strong neuroprotective effect against neuronal oxidative damage, stronger than parent estrogens.¹¹⁴

↑ Potential Causes of Increased 4-MeO Metabolites

- Generally considered beneficial
- Increased 4-OH intermediates resulting in efficient methylation
- Increased COMT activity from genetic polymorphisms or increased demand

↓ Potential Causes of Decreased 4-MeO Metabolites

- Decreased 4-OH intermediates
- Decreased or inefficient COMT enzyme activity from genetic polymorphism or decreased demand
 - Consider addition of methyl donors to increase 4-MeO metabolites (See above COMT Support section)
- Use ratio of 4-MeO/4-OH estrogens to aid in assessment of COMT methylation activity

Ratio of 4-MeO E1/4-OH E1 and 4-MeO E2/4-OH E2

↑ Causes of Increased 4-MeO/4-OH Ratios

Methylation of the 4-hydroxyestrogens is beneficial as this renders them inert and prevents estrogen quinone formation. A higher ratio is related to a lower risk of breast cancer.

- Elevations of the ratio of methylated estrogens to catechol estrogens generally indicate increased COMT activity and methylation.
 - Thus, the conversion of catechol estrogen into methylated versions is considered optimal for reduction of carcinogenic potential.
- Causes are similar to causes of increased 4-MeO metabolites, refer to 4-MeO section.

↓ Causes of Decreased 4-MeO/4-OH Ratios

- Low ratios of methoxylated estrogens compared to catechol estrogens indicate reduced or inadequate methylation. This may confer risks of increased DNA damage from oxidative stress of catechol estrogens and quinone conjugates leading to increased breast, and potentially additional, cancer risks.¹¹⁵
- Causes of a reduced ratio are similar to causes of decreased 4-MeO metabolites, refer to MeO section

Additional Estrogen Related Markers

Bisphenol A

Bisphenol A (BPA) is an industrial chemical which is produced in large volumes primarily for polycarbonate plastics and epoxy containing products. Humans have become ubiquitously exposed to BPA through food packaging, paper receipts, sports equipment, medical and dental supplies, and subsequent contamination in dust, food, and water. BPA is a known 'endocrine disruptor'. BPA has a low binding affinity with estrogen receptors, however, is thought to be able to stimulate estrogenic effects at low doses. BPA monitoring and avoidance is advised for conditions where estrogen excess related illness is a concern.¹¹⁶ In addition to estrogenic concerns, studies show BPA exposure promotes the development of obesity and may play a role in the development of diabetes, hypertension, male and female infertility, hormone dependent tumors, PCOS, liver inflammation, COPD, cardiovascular disease as well as increased all-cause mortality.^{117,118}

↑ Potential Causes of Increased BPA levels

Increased levels of BPA can indicate recent and/or ongoing exposure to plastic or epoxy resins. Items that frequently contain BPA include the following:


- Water bottles (#3, 6, or 7), children's toys, food can linings, food packaging, food storage containers, eyeglass lenses, baby bottles and nipples, water supply pipes, medical equipment, sports safety equipment, electronic devices, CD/DVD discs, thermal paper receipts, and dental sealants¹¹⁹
- Lifestyle considerations for elevated BPA
 - o BPA avoidance is advised wherever possible.
 - o Antioxidants
 - Curcumin, NAC, melatonin, selenium, alpha lipoic acid, vitamin A, quercetin, lycopene, and vitamin E have shown mitigation of BPA harms in pre-clinical evidence.^{118,120}
 - o Increased elimination
 - Sauna has shown the ability to increase excretion of BPA through the skin¹²¹
 - Wheat sprout juice (*Triticum aestivum*) (100ml/day) consistently reduced urine BPA levels in a cohort of young women, suggesting potential detoxification of BPA-toxicity via antioxidative mechanisms and interference of absorption and distribution.¹²²

Total Estrogen

The total estrogens are the sum of all estrogens and estrogen metabolites. High or low levels of these values can be used in conjunction with individual estrogen values as part of an overall estrogen dominance or deficiency evaluation.


Androgens

Androgens	DHEA	Etiocholanolone
	Testosterone	5 α -DHT
	Epi-Testosterone	5 α ,3 α -Androstanediol
	Androstenedione	5 α -Androstanediol
	Androsterone	DHEA-S



Total Estrogen

DHEA and DHEA-S are quantitatively the most abundant steroid hormones in circulation in humans. DHEA is considered slightly androgenic while DHEA-S is relatively inactive at the steroid receptor. However, DHEA, and inactive DHEA-S converted to active DHEA, serve as a precursor to androgens and estrogen in the periphery.¹²³ Generally, the level of DHEA and androgen metabolites should trend in similar directions. If a single androgen metabolite is markedly elevated or decreased, then consideration can be given to possible inhibition or upregulation of the enzymatic pathway regulating the metabolite.



DHEA

DHEA begins with the conversion of cholesterol into pregnenolone by the mitochondrial enzyme p450_{scc}. Pregnenolone is then converted into 17-OH pregnenolone by a 17 α -hydroxylase reaction. The 17,20-lyase reaction follows which converts 17-OH pregnenolone to DHEA.

DHEA can be either converted to DHEA-S or go through reactions to produce other androgens, such as androstenedione, etiocholanolone and androsterone, as well as estrone. DHEA also plays a role in the physiology of adipose and studies indicated it can reduce adipose tissue mass, inhibit the proliferation of adipocytes and may protect against obesity.¹²⁴

DHEA-S

DHEA-S occurs by the sulfation of DHEA and this process is catalyzed by the enzyme hydroxysteroid sulfotransferase (HST, *SULT2A1*), commonly known as DHEA sulfotransferase. DHEA-S can also be converted back into DHEA by steroid sulfatase (STS). DHEA-S is the sulphated, most abundant, version of DHEA. Compared to DHEA, DHEA-S it has a longer half-life, does not have diurnal variation, and provides a stable circulating pool from which to measure adrenal androgen activity. Thus, DHEA-S it is often used to measure effects of DHEA supplementation. As mentioned, DHEA and DHEA-S can interconvert through enzymatic actions.

While DHEA-S has little androgenic capacity, it is known to be active as a neurosteroid and a buffer from the effects of oxidation and glucocorticoids.¹²⁵

Typically, females are referred for DHEA/DHEA-S testing for virilization and/or PCOS evaluation. Males are referred for DHEA/DHEA-S testing for congenital adrenal hyperplasia, primary or secondary adrenal insufficiency, adrenal tumors hypertension and alopecia.¹²⁶


↑ Potential Causes of High DHEA/DHEA-S

- DHEA and DHEA-S can be mildly elevated for idiopathic reasons
- Exogenous supplementation
- Androgen secreting adrenal tumor
- Elevated cortisol
- PCOS
- Steroid sulfatase (STS) deficiency (in cases of DHEA-S elevations)¹²⁷
- Precocious puberty
- Congenital adrenal hyperplasia¹²⁸
- If DHEA is high but DHEA-S is low – consider polymorphisms or inflammation as a cause of reduced DHEA sulfotransferase enzyme activity
- Lifestyle considerations for androgen excess can be seen in Table 6.0 in the Testosterone Section.

↓ Potential Causes of Low DHEA/DHEA-S

- Biological aging in males and females
- Chronic stress
- Chronic inflammation¹²⁹
- Primary and secondary adrenal insufficiency
- Hypothyroidism
- Anencephaly in a developing fetus
- Drugs include carbamazepine, dexamethasone, opioids, phenytoin.¹²⁸
- Lifestyle considerations for androgen deficiency, if clinically relevant can be seen in Table 7.0 in the Testosterone Section.

Androstenedione (A4)



Androstenedione is a key intermediate and hormone precursor of estrone and testosterone. Androstenedione is synthesized from DHEA and then either converted into testosterone through the action of 17 β -hydroxysteroid dehydrogenase (17 β -HSD), or to estrone via the aromatase enzyme (AR) complex. It is primarily secreted by the adrenal gland and governed by ACTH. It can be secondarily secreted in the testes or ovaries from DHEA-S. Males rapidly convert androstenedione to testosterone in the testes. Females use androstenedione as a precursor for all estrone production and about half of testosterone production.¹³⁰ Androstenedione is often measured as a diagnostic aid for hyperandrogenism, premature adrenarche, and diagnosis and monitoring of congenital adrenal hyperplasia (CAH). It was once sold as a popular athletic enhancer, however, since 2004, has been classified as a Schedule III-controlled substance and is not available for over the counter or prescription use. In peripheral tissues, testosterone can be converted back to the inactive androstenedione by 17 β -HSD-2, as a means to regulate intracellular androgen levels.¹³¹

↑ Causes of Increased Androstenedione


- Mild elevations are usually idiopathic in men
 - Peripheral conversion to estrogens may produce gynecomastia in men
- May relate to insulin resistance and PCOS in women
- Marked elevations can be caused by adrenal gland disorders
 - Adrenal tumors, congenital adrenal hyperplasia,¹³² Cushing's syndrome
- Androstenedione producing gonadal tumors (testes, ovaries)
- Genetic polymorphisms: CYP19, CYP17A1¹³³
- Exogenous hormone or dietary supplement use
- Elevated 17-beta dehydrogenase-2 which increases peripheral testosterone inactivation to androstenedione
- Lifestyle considerations for androgen excess, if clinically relevant, can be seen in Table 6.0 in the Testosterone Section.

↓ Causes of Decreased Androstenedione

- Aging
- Adrenal insufficiency - low cortisol production
- Osteoporosis

- Medications: such as glucocorticoids and testosterone replacement therapy
- Poor nutrition: A diet that is low in essential vitamins and minerals can also lead to reduced androstenedione
- Rare causes include 17Beta HSD-3 deficiency, STAR (steroidogenic acute regulatory protein) deficiency and 17- α -hydroxylase deficiency
- Lifestyle considerations for androgen deficiency, if clinically relevant can be seen in Table 7.0 in the Testosterone Section

The Metabolites of Androstenedione: Androsterone and Etiocholanolone



Assessment of androstenedione metabolites can provide insight around androsterone and etiocholanolone levels as well as insight into 5 α and 5 β -reductase activities. 5 α -reductase, as a precursor to androgenic metabolites, is often a therapeutic target for androgen excess conditions, while 5 β -reductase yields less androgenic effects.

Androsterone

Androsterone is a metabolite of androstenedione, governed by enzymatic effects of 5 α -reductase. It is one indicator of 5 α -reductase activity. Androsterone is thought to be devoid of androgenic activity. Androsterone has been found to have neurosteroid effects and anticonvulsant properties through modulation of the GABA receptor.¹³⁴ Androstenedione, the parent compound of androsterone and etiocholanolone, along with both androsterone and etiocholanolone, have been found to be increased in androgenic alopecia.¹³⁵

↑ Causes of Increased Androsterone

- Increased levels of DHEA
- Increased levels of androstenedione
- Increased activity of 5 α -reductase

↓ Causes of Decreased Androsterone

- Decreased levels of DHEA
- Decreased levels of androstenedione
- Decreased activity of 5 α -reductase
 - Finasteride is a 5 α -reductase inhibiting drug, used for alopecia and benign prostatic hypertrophy (BPH), which can cause decreased androsterone¹³⁶

Etiocholanolone

Etiocholanolone is a metabolite of androstenedione, governed by the enzymatic actions of 5-beta reductase. It is thought to be one indicator of 5-beta reductase activity. It does not have androgenic action and is classified as a ketosteroid. It has a role in fever production through stimulation of interleukin-1 and other pro-inflammatory cytokines from leukocytes.¹³⁷

Both etiocholanolone and androsterone have been found to have neurosteroid effects and anticonvulsant properties through modulation of the GABA receptor.¹³⁴ Also, both androsterone and etiocholanolone, along with androstenedione, have been found to be increased in androgenic alopecia.¹³⁸

↑ Causes of Increased Etiocholanolone

- Increased levels of DHEA
- Increased levels of androstenedione
- Increased activity of 5-beta reductase

↓ Causes of Decreased Etiocholanolone

- Decreased levels of DHEA
- Decreased levels of androstenedione
- Decreased activity of 5-beta reductase

If either metabolite of androstenedione is high while the other is low, then the enzymatic pathway leading to that metabolite can reflect an upregulation in that preferred metabolic pathway. For example: high androsterone level with low etiocholanolone level may indicate that 5-alpha-reductase pathway is more upregulated or active. The 5-alpha-reductase pathway also governs DHT metabolite production, which is androgenic, so looking at these two hormones, alongside glucocorticoid metabolites, may provide helpful information (see insert on 5 α -reductase enzymes below).

💡 Spotlight on 5 α -Reductase Enzymes

The 5 α -reductase family of NADPH dependent enzymes reduces several steroids to more potent forms. Much is known about the 5 α -reductase pathway due to its role in converting testosterone to the more potent dihydrotestosterone (DHT), primarily in target tissues such as skin, liver, and prostate. Pharmaceutical targeting of this pathway is common for treatment of BPH and lower urinary tract symptoms in men and androgen stimulated skin/hair disorders in both genders.

Less is commonly known about the role of 5 α -reductase reduction activities in progesterone, glucocorticoid, and mineral corticosteroid pathways. Research on finasteride, a prescription 5 α -reductase inhibitor, has indicated several systemic side effects resulting from actions across multiple hormone targets during and after 5 α -reductase inhibition therapy. In addition to expected androgen related sexual health side effects, effects also include cognitive complaints, gynecomastia, depression and suicidality, anxiety, and headaches, among other complaints.^{139,140}

Natural compounds with 5 α -reductase inhibitory activity have this action through diverse plant constituents and mechanisms. Consequently, many natural compounds with 5 α -reductase activity do not appear to have the same severity of side effects when used in moderate dosages.¹⁴¹

How Do You Evaluate 5 α -Reductase Levels?

Evaluation of the ratio of metabolites produced from 5 β -reductase and 5 α -reductase pathways can highlight which may be more metabolically active. Vibrant's Urinary Hormones Panel measures the following metabolites for comparison.

- Androgen Metabolites: Androsterone vs. Etiocholanolone, Testosterone vs. 5 α -DHT
- Glucocorticoid Metabolites: Ratio of α -Tetrahydrocortisol to β -Tetrahydrocortisol¹⁴²
- Progesterone Metabolites: 5- α Pregnenediol vs. β -Pregnenediol

What is Associated with High Levels of 5 α -Reductase?

- Obesity¹⁴³
- Insulin resistance
- Elevated DHEA
- PCOS
- Hirsutism


Natural Compounds with 5 α -Reductase Inhibition Activity

- Flax Lignan¹⁴⁴ (Typical dose 10–30g/day)
- Saw Palmetto^{145,146} (Typical dose 160mg/2x per day)
- Stinging Nettle^{145*}
- Free Fatty Acids^{145*} (i.e., GLA, DHA)
- Green Tea, Quercetin, and other Polyphenols^{147*}
- Reishi Mushroom^{148*}
- Zinc^{149*}
- Pumpkin Seed Oil^{150*}

*Dose which increases 5 α -reductase inhibition activity is not well established

Testosterone (T)

Testosterone is produced primarily by the gonads, i.e., testes and ovaries, in men and women. Secondly, it is produced by the adrenal glands, as well as in peripheral tissues from metabolites originating in gonads and adrenal glands. The main precursor is DHEA, which is converted to androstenedione or androstenediol by the enzyme 3 β -HSD, and then to testosterone via 17 β -HSD or 3 β -HSD2 respectively.



Testosterone is the predominant androgen produced in males and females. In both males and females, testosterone supports reproductive function and libido, maintains muscle mass and bone structure, supports cardiac health and promotes optimal brain function.¹⁵¹ In males, at puberty, testosterone promotes the development of the male sexual organs and secondary sex characteristics such as deep voice, body hair and libido. It also contributes to the anabolic status of tissues such as red blood cells, muscle mass, linear growth, and bone density.¹⁵²

In females, ovarian production of testosterone increases during the follicular phase of the menstrual cycle and reaches the highest levels at ovulation and the luteal phase. In advancing age, ovarian production of testosterone decreases gradually throughout age, rather than suddenly in a menopausal transition.¹⁵³ Nonetheless, by menopause total testosterone levels in women aged 65–74 years is approximately one-third that observed in women who are aged 20 years.¹⁵¹

The benefit of urinary testosterone testing is that it is more closely correlated with serum free testosterone rather than total testosterone levels. Levels of serum total testosterone are often influenced by sex hormone binding globulin levels which are affected by numerous conditions such as obesity and metabolic syndrome, thyroid disorders, steroid use, PCOS, pregnancy, etc.¹⁵⁴

↑ Potential Causes of High Testosterone

- Endogenously high testosterone in males is uncommon, however can be a result of precocious puberty, adrenal hyperplasia or tumor, testicular tumor, and CNS lesion.
- Exogenous supplementation from prescriptions or adulterated over the counter “libido boosting” or “male enhancing” supplements
- In females, high testosterone can be a result of idiopathic hirsutism, polycystic ovary syndrome, abnormal menstrual cycles, congenital adrenal hyperplasia, ovarian tumors, and intersex physical characteristics.
- Cigarette smoking due to a nicotine metabolite which inhibits androgen breakdown
- Pollutant exposure (polychlorobiphenyls, hexachlorobenzene, etc.)
- Eating disorders
- Drugs - anticonvulsants, atrial natriuretic hormone, barbiturates, cimetidine, clomiphene, estrogens, gonadotropin (males), kaliuretic hormone, oral contraceptives, and vessel dilator hormone¹⁵⁵

↓ Potential Causes of Low Testosterone

- There can be many common causes of low testosterone such as: Increased body mass, heavy alcohol use, hypopituitarism, hyperprolactinemia, hypothyroidism, and late-onset hypogonadism (andropause).¹⁵⁶ In women, oophorectomy and menopause are also common contributing causes.¹⁵⁵
- Other causes can be the following: cirrhosis, COPD (moderate to severe), Klinefelter syndrome, Down syndrome, obstructive sleep apnea, end-stage renal disease, adrenal insufficiency, epilepsy, trauma to gonads or head, hemochromatosis, human immunodeficiency virus, and male hypogonadism.
- Drugs which contribute to low testosterone include anabolic steroids, cyproterone, dexamethasone, diethylstilbestrol, digitalis, digoxin (males), estrogen therapy (increases SHBG), ethyl alcohol, glucose, glucosteroids, gonadotropin-releasing hormone analogs, finasteride, halothane, ketoconazole, metoprolol, metyrapone, opioids, phenothiazines, spironolactone, and tetracycline.¹⁵⁵
- It is important to note that individuals with a polymorphism of the UGT2B17 gene have low or undetectable urinary testosterone levels, thus in these individuals’ serum or saliva testing of testosterone will provide a more accurate result in free testosterone levels.

Treatment Considerations for High Testosterone and Low Testosterone

Testosterone/Androgen Excess

Treatment Objectives: Consider which of the following single or combination of therapeutic approaches is the best fit for the patient with symptoms related to androgen excess.

- Reduce circulating free androgens
- Increase insulin sensitivity
- Increase sex hormone binding globulin (SHBG) to decrease circulating free androgens
- Weight loss if obesity is present

Table 6. Holistic Treatment Considerations for Androgen Excess

Lifestyle Considerations	Reduce Free Androgens/Improve Insulin Sensitivity	Interventions to Increase SHBG (Sex Hormone Binding Globulin)
Specific Supplements	<p>NAC in PCOS patients 600mg/3x per day¹⁵⁷</p> <p>Vitamin D 4000 IU/day¹⁵⁸ Additionally, repletion in deficient women (to blood levels above 32ng/mL) ↑SHBG and ↓bioavailable E2 and testosterone¹¹</p> <p>Myo-inositol 4g/day or D-chiro inositol 1g/day¹⁵⁹⁻¹⁶¹</p> <p>Soy isoflavones 50g/day¹⁶¹</p> <p>Berberine 400mg/TID – improves insulin sensitivity¹⁶²</p> <p>Chromium* Improve insulin resistance, may lower DHEA levels in adipose tissue in PCOS patients¹⁶³</p> <p>Zinc, Magnesium and Selenium to counter deficiencies in diet¹⁶⁴</p>	<p>DIM 300mg/day¹⁶⁵</p> <p>Refer to Insulin Modulating interventions increased insulin decreases SHBG¹⁶⁶</p>
Diet	<p>DASH Diet reduce insulin resistance promotes weight loss in PCOS¹⁶⁷</p> <p>Mediterranean Diet reduces insulin¹⁶⁷</p> <p>Calorie Restriction Diets reduce insulin resistance promotes weight loss in PCOS¹⁶⁷</p> <p>Low Glycemic/Ketogenic diet¹⁶⁴</p> <p>Almonds 46 g/d¹⁶⁸</p> <p>Walnuts 36 g/d¹⁶⁸</p> <p>Prebiotics/ Probiotics* -Increased microbiome diversity¹⁶⁹</p> <p>Soy isoflavones in diet 50g/d¹⁶¹</p> <p>Spearmint tea¹⁶⁴</p> <p>Green tea¹⁶⁴</p> <p>Licorice tea¹⁶⁴</p> <p>Flax seed meal 30g/day¹⁷⁰</p>	<p>Coffee, Tea and other caffeine containing beverages^{171*}</p> <p>Walnuts 36g/day¹⁶⁸</p> <p>Soy milk 30g/day¹⁷¹</p> <p>Olive Oil^{173*}</p> <p>High Fiber diet, vegetarian diet¹⁷³</p> <p>Insulin Modulating interventions as increased insulin decreases SHBG¹⁷⁴</p>

Exercise/Sauna	Vigorous aerobic exercise – minimum of 120 minutes per week alone or with resistance training reduce insulin resistance in PCOS, aid body composition ¹⁶⁴	Moderate intensity exercise 225 min/week ¹⁷⁵
Misc. Factors	Optimizing sleep aids with regulation of dysglycemia ¹⁶⁴	Drugs that increase SHBG include Tamoxifen and other selective estrogen-receptor modulators (SERM's), oral estrogen, metformin, and anti-seizure medications

* Dosage information is variable or not well established for this specific benefit; refer to empiric or published standards of safe use

Additional treatment options for androgen excess:

- If hirsutism is present, and signs of elevated 5 α -reductase upregulation are shown, consider modulation of this pathway, see 5 α -reductase insert.
- Pharmacologic therapy such as spironolactone, or spironolactone in combination with licorice¹³²

Further testing considerations for androgen excess:

- **Saliva Hormone** test to determine bioavailable and tissue levels of hormones
- **Environmental Toxins and Heavy Metals** tests to assess xenoestrogen compounds
- **Cardiac Health and Diabetes Panel** to assess cardiometabolic function
- **Thyroid Panel**
- **IGF-1**
- **17-OH Progesterone** – for non-classical congenital adrenal hyperplasia
- **Prolactin** (can be ordered alone)
- **Liver function tests** (included in Comprehensive Metabolic Panel, cannot be ordered alone)
- **SHBG** (included in Serum Hormones Panel, can be ordered alone)

Testosterone/Androgen Deficiency

Treatment Objectives: Consider which of the following single or combination of therapeutic approaches is the best fit for the patient symptoms of testosterone/androgen deficiency.

- Reduce aromatase activity
- Support symptoms of low androgens
- Reduce obesity and insulin resistance to reduce aromatase production
- Reduce exposure to phthalates and heavy metals

Table 7. Holistic Treatment Considerations for Testosterone Deficiency

Lifestyle Interventions	Symptomatic Support/Increase Testosterone (T)	Decrease Aromatase Activity
Supplements	<p>Withania somnifera (Ashwagandha) studied to boost testosterone in men (i.e KSM-66 300mg/BID or TID), Shoden brand 21mg withanolide glycoside/d or 5g root powder/day^{176,177}</p> <p>Trigonella foenum-graceum (Fenugreek) 500-600mg/day seed extract boosts testosterone & libido in males¹⁷⁸</p> <p>Mucuna pruriens 5g/day seed powder improves testosterone in infertile men and boosts semen quality¹⁷⁹</p> <p>Eurycoma longifolia (Long jack) 200-400mg/day increased libido, testosterone, ED in men,¹⁷⁸ conflicting studies on efficacy</p> <p>Korean Red Ginseng 1500mg/BID improves serum testosterone in men with metabolic syndrome,¹⁸⁰ conflicting studies on efficacy¹⁷⁶</p> <p>DHEA >50g/day^{181**} may boost testosterone levels, especially in women</p> <p>Nigella Sativa improves semen parameters in infertile men¹⁷⁸</p> <p>Tribulus terrestris 750mg/day improved sexual enhancement in women,¹⁷⁸ several studies show no increase in testosterone in men,¹⁷⁶ may enhance libido</p> <p>Maca 1500mg root/BID improvement in testosterone/libido/sexual function in women on SSRI,²⁹ several studies in men show no increase in testosterone¹⁷⁸</p> <p>Saw palmetto Increased T, with decreased DHT and Estrogen in men at 800mg-2000mg/day¹⁴</p>	<p>Numerous flavonoids show preclinical evidence of decreased aromatase activity: *</p> <p>Resveratrol¹³</p> <p>Grape Seed extract¹⁴</p> <p>Citrus Peels/citrus flavonoids¹⁸²</p> <p>Apigenin¹³</p> <p>Chrysin (Honey/propolis)¹⁶</p>
Diet	<p>Diets that aid weight loss - In obese men and/or men with metabolic syndrome, weight loss through a variety of diet and exercise or bariatric surgery programs consistently increase testosterone¹⁸³</p>	<p>White button mushrooms^{14*}</p> <p>Foods with apigenin and/or resveratrol: grapes, berries, celery, parsley, onions, oranges, chamomile, thyme, basil, oregano, red wine¹⁸⁴ Kale, collards¹⁴</p>
Exercise/Sauna	<p>Exercise programs that aid weight loss - In obese men and/or men with metabolic syndrome, multiple, varied diet/exercise programs that resulted in weight loss significantly improve testosterone levels.¹⁸³</p>	<p>150 minutes moderate activity per week¹⁸⁵</p>


Stress Management	<p>Optimize sleep, sleep deprivation (<5hours/night) linked to decreased testosterone¹⁸³</p> <p>High occupational stress is linked to lower testosterone levels.¹⁸³ Stress management interventions may be helpful.</p>	
Miscellaneous Support	<p>Aromatase inhibitors, SERM's and hCG have been shown to raise testosterone¹⁸⁶</p> <p>Varicocele repair significantly improves testosterone levels¹⁸³</p> <p>Testosterone replacement therapy (TRT) is an option for hypogonadal men. It has shown to improve sexual function, muscle strength, bone density, and mood and cognition***</p>	Metformin may aid increases in testosterone in men with metabolic syndrome ¹⁸⁷

* Dosage information for this indication is not well established, refer to empiric guidelines of safe and effective use.

** May increase estrogen and other hormone levels, use under physician supervision and hormonal lab monitoring.

*** Fully research the individual patient's risk/benefit analysis regarding hormone related cancer and other side effects of testosterone replacement therapy (TRT) prior to recommending. Side effects of TRT include erythrocytosis, male infertility, testicular atrophy, and gynecomastia.¹⁸⁶ There are also concerns but inconclusive data regarding the role of exogenous testosterone with cardiovascular risks, prostate cancer risks and venous thromboembolism.¹⁸⁶

Epi-Testosterone (Epi-T)



Epi-testosterone is an isomer of testosterone that is endogenously produced by the gonads in response to LH. It is a competitive antagonist of the androgen receptor and as such, demonstrates a weak anti-androgen effect. It is also a strong inhibitor of 5-alpha-reductase activity in studies.¹⁸⁸ The ratio of testosterone to epi-testosterone is used to monitor anabolic drug use in athletic competitions.

↑ Causes of Increased Epitestosterone

- Highest levels are seen in adolescent males
- Levels can be elevated in exogenous epi-testosterone use

↓ Causes of Decreased Epitestosterone

- Exogenous use of testosterone or other androgens
- Oral contraceptive use
- Levels can vary with menstrual cycle fluctuations

T/Epi-T Ratio (Testosterone/Epi-Testosterone Ratio)

This ratio is a well-known biomarker of androgen abuse or doping. A typical T/Epi-T ratio is about 1:1. Use of exogenous synthetic testosterone or testosterone precursors results in elevations in urinary testosterone with a corresponding decrease in urinary epitestosterone. While some factors can alter the ratio aside from doping, a urinary T/Epi-T ratio value which is ≥ 4 is suspicious for exogenous use of androgens. In evaluation of androgen abuse, serial measurement of individual T/Epi-T ratios are preferred, rather than single point in time testing.


↑ Causes of Increased T/Epi-T ratio

- Exogenous androgen use
- Increased alcohol consumption, particularly in women

↓ Causes of Decreased T/Epi-T Ratio

- Ketoconazole use¹⁸⁹
- Genetic polymorphism of UGT2B17 (common in Asian populations)
 - o Lowers urinary levels of testosterone, 5 α -DHT and 5 β -Androstanediol

T/Epi-T Ratio (Testosterone/Epi-Testosterone Ratio)



5 α -dihydrotestosterone (DHT) is an androgenic hormone that is formed from testosterone through the action of the enzyme 5 α -reductase in the gonads, as well as other tissues such as prostate, skin, and liver. DHT can also be indirectly formed through “backdoor pathways” and the action of 5 α -reductase on other androgen intermediates.¹⁹⁰

DHT is a more potent androgen than testosterone due to increased affinity for the androgen receptor. Biologically, DHT plays important roles in the maturity of male genitalia and secondary sex characteristics during male development stages, as well as a primary role in prostate growth and function. In females, DHT is known to contribute to puberty and the development of secondary sex characteristics, however other roles are somewhat unclear. Metabolites of DHT have beneficial neurosteroid and cognitive effects. Unlike testosterone, DHT cannot be aromatized into estradiol.¹⁹¹ DHT levels are currently used to monitor 5 α -reductase deficiency, male pattern baldness, and monitoring of 5 α -reductase therapy.


↑ Potential Causes of Increased 5 α -DHT

- PCOS and hirsutism in females
- Androgenic alopecia
- BPH
- Increased 5 α -reductase activity
 - Refer to Spotlight on 5 α -Reductase Enzymes section
 - Obesity, insulin resistance, DHEA supplementation
- Lifestyle considerations for androgen excess, if clinically relevant, can be seen in Table 6.0 in the Testosterone Section.

↓ Potential Causes of Decreased 5 α -DHT

- Reduced testosterone
- Prostatectomy
- 5 α -reductase inhibition or deficiency
 - Refer to Spotlight on 5 α -Reductase Enzymes section
- Hypopituitarism
- Lifestyle considerations for androgen deficiency, if clinically relevant can be seen in Table 7.0 in the Testosterone Section

DHT Metabolites: 5 α ,3 α -Androstanediol and 5 β -Androstanediol



5 α ,3 α -Androstanediol

5 α -DHT is reduced to 5 α , 3 α -androstanediol via the activity of the enzyme 3 α -HSD. 5 α , 3 α -androstanediol (often referred to as 3 α -diol) can be a useful marker for peripheral androgen activity in conditions such as PCOS, hirsutism and acne. 5 α , 3 α -androstanediol may be a more sensitive marker than DHT for peripheral androgen production because DHT is often bound to sex hormone binding globulin and it is metabolized rapidly.^{192,193} Also considered a neurosteroid, 5 α , 3 α -androstanediol has shown to have anxiolytic effects through actions on the GABA receptor as well as estrogenic effects from binding to estrogen receptor beta (ER- β).¹⁹⁴

↑ Causes of Increased 5 α , 3 α -androstanediol (known as 3 α -diol)

- Increases in 5 α -DHT
 - PCOS and hirsutism in females
 - Androgenic alopecia
 - Benign prostatic hypertrophy
 - Increased 5 α -reductase activity
 - Obesity, insulin resistance, high DHEA
- Increased 3 α -HSD activity

↓ Causes of Decreased 5 α , 3 α -androstanediol (known as 3 α -diol)

- Decreases in 5 α -DHT
 - Reduced testosterone
 - Prostatectomy
 - 5 α -reductase inhibition or deficiency
 - Hypopituitarism


5 β -Androstanediol

Testosterone can also be metabolized by 5 β -reductase into 5 β -DHT, a metabolite with insignificant androgenic action. 5 β -DHT is reduced to 5 β -androstanediol (also known as etiocholanediol or 5 β -androstane-3 α ,17 β -diol) through the action of 3 β -HSD. Clinical implications of high and low levels of 5 β -androstanediol have not been well studied. High and low levels of 5 β -androstanediol are affected by upstream hormone and metabolite levels such as DHEA, testosterone, and 5 β -DHT values, as well as 5 β -reductase and 3 β -HSD enzyme activity.

Progestogens

Progestogens	Progesterone	20 α -Dihydroprogesterone
	Allopregnanolone	β -Pregnanediol
	Allopregnanediol	α -Pregnanediol
	3 α -Dihydroprogesterone	

The Progesterone Metabolites



Progesterone and its metabolites are vital to normal health and reproduction and have roles across many systems in both men and women. While the role of progesterone in normal menstruation and pregnancy is well known, there are other systemic roles that progesterone fulfills. For example, progesterone is a precursor to critical steroid hormones such as aldosterone, cortisol, testosterone, and estradiol. Also, progesterone and its metabolites have important roles in the cardiovascular, renal, and musculoskeletal systems. Further, progesterone is a neuroprotectant, neuromodulator and aids with relaxation and sleep. Importantly, progesterone plays a role in immune support and cancer protection against endometrial, colorectal cancers and potentially others. In addition to these benefits, in men, progesterone also influences spermiogenesis and testosterone biosynthesis.

Progesterone metabolizes into isomers of dihydroprogesterone, pregnanolone and pregnanediol. All of these metabolites then undergo phase I and phase II reactions and result in significantly more progesterone metabolites.¹⁹⁵

Beta(β)-Pregnanediol

Progesterone itself cannot be measured directly in urine, therefore, beta (β) and alpha (α) pregnanediol, the main progesterone metabolites in urine, are used in urine as surrogate progesterone markers. β-pregnanediol is considered the primary metabolite resulting from progesterone metabolism and is inactive. β-pregnanediol is produced through 5β-reductase enzyme activity. Factors which elevate or decrease progesterone will correspondingly increase or decrease β-pregnanediol. In addition, upregulation or down regulation of 5β-reductase enzyme activity can influence overall levels of β-pregnanediol.

↑ Potential Causes of Increased Progesterone/β-Pregnanediol

- Pregnancy
- Exogenous progesterone supplementation (non-topical forms)
- Ovarian cysts
- Ovarian tumors and testicular tumors
- Adrenal hyperplasia
- Stress and caffeine have also been linked to slightly elevated levels.^{196,197}

↑ Potential Causes of Decreased Progesterone/ β -Pregnanediol

- Irregular or anovulatory cycles, including during breastfeeding
- PCOS
- Aging and menopause
- Thyroid disorders
- Obesity
- Over exercise
- Hyperprolactinemia
- Anorexia
- Long term use of NSAIDs¹⁹⁸
- Oral contraceptives¹⁹⁹
- Endometriosis
- Environmental toxicants such as phthalates, pesticides, herbicides, etc.,²⁰⁰
- Low progesterone in pregnancy can be caused by ectopic pregnancy or complications or failure with the fetus or placenta
 - Later in pregnancy, low progesterone can be caused by toxemia or pre-eclampsia of pregnancy
- In men, while sudden dramatic reductions in progesterone are uncommon, waning progesterone can occur as a sequelae of reduced androgens that occur gradually over time after the 4th decade in men²⁰¹

Table 8. Lifestyle Considerations to Support Reduced Progesterone or PG/E2 Ratio

1. Vitex agnus castus 40mg per day²⁰²
2. Vitamin C 750mg/day²⁰³
3. B6 pyridoxine 200mg per day²⁰⁴
4. White peony root daily tea or 3-5 g of root/day²⁰⁵
5. Evening Primrose Oil²⁰⁶
6. Increase foods high in zinc, B6, Vitamin C, Magnesium²⁰³
7. Consider Progesterone replacement therapy with oral micronized, transdermal, or compounded prescription.
8. Assess risk/benefit analysis for the specific individual prior to giving hormone replacement therapy

Alpha-Pregnanediol (5 α -Pregnanediol)

5 α -pregnanediol has been studied to stimulate GABA receptor action and promote relaxation and sleep.²⁰⁷ 5 α -pregnanediol can be a substrate of allopregnanolone and can convert into allopregnanolone through the action of 20-hydroxysteroid dehydrogenase.²⁰⁸

↕ Potential Causes of Increased and Decreased 5 α -Pregnanediol

- Elevations and decreases in 5 α -pregnanediol correspond to those listed for general progesterone and β -pregnanediol, (see above) however may independently vary due to increases or decreases in 5 α -reductase activity.
- Medications which block 5 α -reductase activity, such as those used for enlarged prostate and male pattern baldness, can correspondingly cause decreases in 5 α -pregnanediol.

Allopregnanolone

Allopregnanolone (3 α -THP) is formed when progesterone is reduced by 5 α -reductase to 5 α -dihydroprogesterone and this is further reduced by 3 α -hydroxysteroid dehydrogenase to create allopregnanolone. Allopregnanolone is a regulator of female reproductive function and lactation, but also has been greatly studied for its neurosteroid effects. Allopregnanolone is a modulator of the GABA^A receptor.

↑ Potential Causes of Increased Allopregnanolone

- Typically increases or decreases in allopregnanolone relate to increased or decreased progesterone levels.
- Refer to the Beta (β)-Pregnanediol section for a complete listing of factors affecting progesterone levels.
- Common scenarios causing increased allopregnanolone include the luteal phase of the menstrual cycle, oral progesterone supplementation and pregnancy.
- It has been found that premenopausal women have higher allopregnanolone than postmenopausal women.²⁰⁹

↓ Potential Causes of Decreased Allopregnanolone

- Common scenarios causing decreased allopregnanolone include situations of low progesterone such as perimenopause and menopause, luteal phase defect or anovulation, and pharmaceutical interventions such as finasteride or oral contraceptives.
- Variations in enzyme activity of 5 α -reductase, such as those seen with pharmaceuticals, as well as 3 α -hydroxysteroid dehydrogenase levels, can also affect allopregnanolone levels.
- Decreased production of allopregnanolone or fluctuations in levels have been associated with increased depression and anxiety, premenstrual dysphoric disorder, increased seizure activity in epileptics, Parkinson's disease and Alzheimer's disease.^{210,211,212}


💡 A Note About 4-Pregnanes and 5 α -Pregnanes

Progesterone metabolites that are formed from 5 α reductase activity have been termed 5 α -pregnanes. As indicated, these metabolites have been well studied for calming, anti-seizure, and sleep benefits due to stimulation of GABA^A receptors. Metabolites that have been formed through other enzymatic reactions have been termed 4-pregnanes, such as 3 α -dihydroprogesterone and 20 α -dihydroprogesterone.

Contrary to the calming benefits of 5 α -pregnanes, some pre-clinical studies have demonstrated a potential role 5 α -pregnanes may have in proliferation and cancer promotion effects in breast cancer cells as well as colon, ovarian and glioblastoma tumor tissue.^{213,214} In opposition to this, the 4-pregnanes, such as 3 α -dihydroprogesterone and 20 α -dihydroprogesterone have demonstrated cancer inhibiting effects. Early studies show that in the microenvironment of a tumor cell, the ratio of 5 α -pregnanes to 4-pregnanes is

much higher than a normal cell, which is supportive of tumorigenesis.^{215,216} While more research is needed, it may be pragmatic to evaluate overall levels of 5 α -pregnanes in comparison to 4-pregnenes regarding questions surrounding potential oral progesterone supplementation. This would be most relevant in clinical situations where recurrence of breast cancer and/or other cancer types is a concern.

3 α -Dihydroprogesterone and 20 α -Dihydroprogesterone



3 α -Dihydroprogesterone

3 α -dihydroprogesterone is created from progesterone through the enzymatic actions of 3 α -hydroxysteroid oxidoreductase (3 α -HSD) in tissues. It has been identified as a metabolite which has anti-proliferative effects in tumor cells.²¹⁷ Variations in activity of 3 α -HSD and levels of precursor hormones can affect levels of 3 α -dihydroprogesterone.

20 α -Dihydroprogesterone

20 α -dihydroprogesterone is created from progesterone through the enzymatic action of 20 α -hydroxysteroid oxidoreductase (20 α -HSD) enzymes. 20 α -dihydroprogesterone has a lower affinity for the progesterone receptor and is significantly less potent than progesterone.²¹⁸ 20 α -has been studied to have drug metabolizing, anti-aromatase, and anti-proliferative effects.²¹⁹

Variations in 20 α -HSD enzymes and precursor hormones can contribute to elevations and reductions in 20 α -dihydroprogesterone.

β -Pregnanediol/E2 (β -pregnanediol/estradiol ratio)

The β -pregnanediol/estradiol (β -pregnanediol/E2) ratio is used similarly to the progesterone/E2 ratio. This was commonly used empirically as a marker of "Estrogen Dominance," developed, and popularized by John Lee, MD.²²⁰ Conventionally, the Pg/E2 ratio is also used in IVF research for pregnancy rates and has been studied for assorted characteristics of menstrual cycles and fibrocystic breast disease.^{221,222}

↑ Causes of Increased β -Pregnanediol/E2 – High Progesterone Relative to Estradiol


- Exogenous oral progesterone use (most common)
- Increased endogenous progesterone production (i.e., common in pregnancy)
- Decreased estradiol production

↓ Causes of Decreased β -Pregnanediol/E2 – Low Progesterone Relative to Estradiol

- Natural aging in men and women – common in 4th and 5th decades of life and later
- Sub-optimal clearance of estradiol and/or estradiol metabolites
 - Detoxification pathways
 - Environmental xenoestrogen
 - Exposure

Glucocorticoids and Glucocorticoid Metabolites

Glucocorticoids	Cortisol	α -Tetrahydrocortisol (α -THF)
	Cortisone	β -Tetrahydrocortisone (β -THE)
	β -Tetrahydrocortisol (β -THF)	Corticosterone



Cortisol is widely known as “the stress hormone,” and nearly all tissues of the body have glucocorticoid receptors to respond to the actions of cortisol. Cortisol plays a significant role in maintaining glucose and protein homeostasis, mediation of the stress and immune response, and suppression of inflammation.²²³

The hypothalamus-pituitary-adrenal axis (HPA axis) regulates production and secretion of cortisol. It does this through release of CRH, corticotropin releasing hormone, from the hypothalamus which signals ACTH, adrenocorticotropic hormone, release from the pituitary gland to the adrenal cortex, which then releases cortisol. After cortisol is released in response to these signals, cortisol sends a negative feedback loop to suppress further production of ACTH and CRF. See Figure A.

The HPA axis function follows a diurnal pattern of release, therefore cortisol levels are highest in the morning after waking and lowest at night around bedtime.²²³ Aside from diurnal secretion of hormones, and measured pulsatile releases of hormones, HPA axis function and release of cortisol are also triggered by stressors, both acute and chronic.





Figure A.

Total Cortisol



Total cortisol is a measure of the sum of free cortisol and conjugated cortisol excreted in the urine. Cortisol undergoes steroidogenesis similarly to other steroid hormones, originating from cholesterol and progressing to progesterone. From 17-OH progesterone, cortisol goes through two hydroxylation steps to arrive at 11-deoxycortisol which is then further hydroxylated by 11 β -hydroxylase to arrive at cortisol. In the tissues, the glucocorticoids (cortisol, cortisone, and corticosterone) are regulated by 11 β -hydroxysteroid dehydrogenases, type 1 and type 2.

↑ Table 9. Potential Causes of Increased Cortisol

Cushing's Disease/Syndrome	Infectious Disease	Psychological Stress
Adrenal Neoplastic Disease	Crohn's Disease	Obesity
Hyperpituitarism or Hyperthyroidism	Stress: Heat/Cold/Trauma/Pain	Night Shift Work/Circadian Disruption
Diabetes Mellitus (uncontrolled)	Burns	Alcoholism and Alcohol Withdrawal
Shock	Surgery/Post Op Recovery	Anorexia/Hypoglycemia
Cigarette Smoking	Chronic Renal Disease	HIV
Pregnancy	Virilism	Smoking
Hypertension	Exercise or Chronic Overexercise	Glucocorticoid Resistance
Eclampsia	Severe Hepatic Disease	
Opioid Withdrawal	Depression/Neuropsychiatric disease	


Drugs which can increase cortisol levels include glucocorticoids, caffeine, nicotine, corticotropin, estrogens, oral contraceptives, yohimbine, and vasopressin.²²⁴ Initial or short-term use of marijuana can increase cortisol, while opposite effects may result from chronic or heavy use.²²⁵

↓ Table 10. Potential Causes of Decreased Cortisol

Primary adrenal insufficiency- due to Addison's disease	Primary adrenal insufficiency- due to adrenal destruction from other causes (malignancy, surgery, etc.)
Primary adrenal insufficiency due to genetic conditions (i.e., CAH - congenital adrenal hyperplasia)	Secondary adrenal insufficiency (i.e., withdrawal from glucocorticoid therapy)
Liver disease	Long term use or recent withdrawal from corticosteroid therapy
Pituitary adenoma	Craniopharyngioma/hypophysectomy
	Hypopituitarism
Postpartum pituitary necrosis	Rheumatoid arthritis
Early traumatic experiences, chronic stress, PTSD	Hypothyroidism
Recurrent infectious processes	Waterhouse-Friderichsen syndrome


Drugs which can decrease cortisol levels include opioids,²²⁶ ketoconazole, rifampin, phenytoin, dexamethasone, dexamethasone acetate, and dexamethasone sodium phosphate. There is some evidence that marijuana affects cortisol response, namely that heavy/chronic use can blunt cortisol response.²²⁷

Total Cortisone



Total cortisone is a measure of the sum of free cortisone and conjugated cortisone excreted in the urine. Cortisol is converted to cortisone, an inactive metabolite by 11-β-hydroxysteroid dehydrogenase type 2 (11β-HSD2). 11β-HSD2 is highly expressed in aldosterone-selective target tissues such as the distal nephron, colon, skin, and the salivary glands.²²⁸

The function of 11β-HSD2 to convert cortisol to cortisone is critical in tissues to protect from overexposure to active cortisol and corticosterone.²²⁹ Otherwise, these active hormones would occupy mineralocorticoid receptors and produce a cascade effect of mineralocorticoid excess symptoms such as sodium retention, hypertension, and hypokalemia.²³⁰ In contrast, in metabolically active tissues such as the liver (20-40% of daily production) adipose tissue, and skeletal muscle, 11β-Hydroxysteroid dehydrogenase, type 1 (11β-HSD1), reversibly catalyzes the 11β-reduction of cortisone to cortisol, and regenerates cortisol within these tissues.²³⁰



Which is More Accurate: Cortisol or Cortisone?

There are specific clinical situations in which there would be expected differing results between cortisol and cortisone. In these clinical situations, cortisone may reflect a greater accuracy than cortisol.

- Systemic glucocorticoid deficiency or excess
- Use of oral glucocorticoid therapy
- Directly following acute stress

While cortisone itself is an inactive metabolite, it has been shown, like cortisol, to directly compare to serum free cortisol levels.^{231,232} In general, cortisone values tend to parallel cortisol values. Thus, similar situational or condition triggers cause elevations or depressions in cortisone values as cortisol values. Diurnal cortisone values have not yet been studied for normative value comparisons.²³¹



11 β -HSD Enzyme Activity

11 β -HSD1 and 11 β -HSD2 enzyme activity affect cortisol and cortisone levels. These enzymes affect local target tissues in an intracrine manner, as in 11 β -HSD2 upregulating in the distal nephron and 11 β -HSD1 upregulating in adipose tissue and the brain. In these scenarios, the effect is intracrine, i.e., local, and does not affect systemic cortisol production.²³³

However, these enzymes also work in an endocrine manner, and may affect cortisone and cortisol levels systemically as well. For example, 11 β -HSD1 in the splanchnic bed generates 30–40% of the total daily production of cortisol in humans, while 11 β -HSD2 in the kidney deactivates a similar percentage.²³³ Thus, up or downregulation of these enzymes can play a role in cortisone and cortisol dynamic systemic levels.

Careful evaluation of cortisol/cortisone ratio (below) and levels of terminal cortisol and cortisone metabolites, 5 β THF, 5 α -THF, and 5 β -THE, can yield insights about systemic 11 β -HSD enzyme activity. For example, the trend of increased cortisol terminal metabolites (β -THF + α -THF) compared to tetrahydrocortisone from cortisone, can indicate increased 11 β -HSD1 activity. Similarly, an examination of the cortisol/cortisone ratio can provide insight on 11 β -HSD2 activity.

Cortisol/Cortisone Ratio

This value measures the ratio of total cortisol to total cortisone.



Potential Causes of Increased Cortisol/Cortisone Ratio


- Upregulated 11 β -HSD1 activity which is linked to the following:
 - Inflammatory diseases such as joint disorders, neurodegenerative disorders, diabetes,²³⁴ metabolic syndrome, cognitive decline, and obesity²³⁵
 - Drugs, foods, and herbs such as grapefruit juice, licorice, glucocorticoids, progesterone²³⁶
 - Synthetic endocrine disruptors such as phthalates, organotins (tin compounds widely used as pesticide, biocides and in PVC materials), alkylphenols²³³



Potential Causes of Decreased Cortisol/Cortisone Ratio

- Upregulated 11 β -HSD2 activity which is linked to the following:
 - Hormones such as estrogen, progesterone, 7-keto-DHEA, thyroid²³³
 - Pregnancy upregulates cortisone and 11 β -HSD2²³⁸⁷
 - PCOS upregulates cortisone and 11 β -HSD2²³⁸
- Supplements, drugs, and foods that upregulate 11 β -HSD2 activity
 - Green tea, coffee, holy basil, curcumin, vitamin A, pravastatin, bile acids.^{239,240,241}

Cortisol and Cortisone Metabolites



Evaluation of cortisol metabolites can provide insight into trends of enzyme activity and show patterns in chronic diseases such as PCOS, obesity, apparent mineralocorticoid excess, hyper or hypothyroidism, anorexia, and others.^{242, 243} In the metabolism process, cortisol is reduced by either 5 β -reductase or 5 α -reductase, often in the liver, to form dihydrocortisol intermediates. These are further reduced by 3 α -hydroxysteroid dehydrogenases to form 5 β -tetrahydrocortisol (5 β -THF) and 5 α -tetrahydrocortisol (5 α -THF), the terminal metabolites of cortisol. Cortisone is reduced by 5 β -reductase to produce tetrahydrocortisone (5 β -THE), the final product of cortisone metabolism.

Cortisol Metabolites: β -Tetrahydrocortisol (β -THF or THF) and α -Tetrahydrocortisol (α -THF or Allo-THF)

These metabolites can be used alongside cortisol measures to aid assessment of glandular cortisol production as well as shine a spotlight on overall cortisol metabolism in the body. β -THF is the metabolite from active cortisol which has been reduced by 5 β -reductase. α -THF is the metabolite from active cortisol which has been reduced by 5 α -reductase.

These markers also offer insight into respective 5 α -reductase or 5 β -reductase enzyme activity, along with other metabolites produced from the same enzymatic pathway. See the Spotlight on 5 α -Reductase Enzymes section for further information regarding specific enzymatic influences.

↑ Potential Causes of Increased β -THF and α -THF

- If cortisol metabolites, β -THF and α -THF, are elevated it can reflect excess cortisol production or increased cortisol metabolism. Conditions such as stress, insulin resistance, obesity and hyperthyroidism can cause increased cortisol production and/or metabolism.
- Refer to Table 9., Potential Causes of Increased Cortisol, for further possible causes.

↓ Potential Causes of Decreased β -THF and α -THF

- If cortisol metabolites are decreased it can potentially illustrate reduced glandular output, but also depending on other markers, reflect decreased metabolism and clearance as seen in adrenal insufficiency, suboptimal liver function and hypothyroid.
- Refer to Table 10., Potential Causes of Decreased Cortisol, for further causes of decreased cortisol metabolites.

Cortisone Metabolite: β -Tetrahydrocortisone (β -THE)

β -THE is the metabolite from cortisone which has been reduced by 5β -reductase. This metabolite, in conjunction with β -THF and α -THF can aid in evaluation of overall cortisol production and metabolism in the body. If this marker is elevated or decreased, it can reflect situations of cortisol excess or deficiency such as those listed in Tables 9. and 10.

This marker can also offer insight into 5β – reductase enzyme activity, along with other metabolites produced from the same pathway. See the Spotlight on 5α -Reductase Enzymes section for further information regarding specific enzymatic influences.

Metabolized Cortisol (THF+THE)

Metabolized cortisol reflects the sum of α -THF, β -THF, and β -THE. Metabolized cortisol is the best measure of overall cortisol production and clearance through the liver. It does not include the fraction of cortisol that is free and bioavailable for use, generally about 5-10% of total production.²⁴⁴

The Mineralocorticoids: Deoxycorticosterone (11-Deoxycorticosterone) and Corticosterone

Deoxycorticosterone (11-Deoxycorticosterone)

Deoxycorticosterone (DOC) is an intermediate steroid hormone synthesized from progesterone through the action of 21α -hydroxylase and is a precursor for the synthesis of cortisol and aldosterone. It has negligible mineralocorticoid activity. Typically, this measure is helpful for diagnosis of disorders of steroid synthesis such as congenital adrenal hyperplasia.²⁴⁵

↑ Potential Causes of Increased Deoxycorticosterone

- Congenital adrenal hyperplasia due to 11 beta-hydroxylase 2 (CYP11B2) deficiency
- Pregnancy and newborns
- DOC producing adrenal tumor
- Cushing's disease
- Primary cortisol resistance syndrome.

↓ Potential Causes of Decreased Deoxycorticosterone

- Congenital adrenal hyperplasia due to 21 -hydroxylase deficiency, 3 -beta-steroid dehydrogenase deficiency, or steroidogenic acute regulatory protein (StAR) deficiency.
- Primary or secondary adrenal insufficiency

Corticosterone

Corticosterone is an intermediate metabolite produced from deoxycorticosterone through 11β -hydroxylase activity in the pathway of aldosterone formation. It circulates at 10–20-fold lower concentrations than cortisol.²⁴⁶ In other species, such as rodents and birds, corticosterone acts as the main glucocorticoid, however, in humans, the role of corticosterone, other than as an intermediate, is somewhat undefined. There is ongoing research into the role of corticosterone in insulin sensitivity, HPA axis function and mineralocorticoid effects.²⁴⁷

Testing levels of corticosterone plays a role, in conjunction with deoxycorticosterone, for diagnosis of disorders of steroid synthesis, especially relating to congenital adrenal hyperplasia and 11 -beta-hydroxylase activity.

Oxidative Stress

Other	Melatonin 8-OHdG
-------	---------------------

8-OHdG

8-hydroxy-2'-deoxyguanosine (8-OHdG, or 8-oxodG) is a marker of oxidative damage to DNA from reactive oxygen species (ROS). The guanine molecule, one of 4 nucleobases in DNA, has a high susceptibility to ROS oxidative damage. When oxidative damage to guanine occurs, 8-OHdG is formed. When this damage is repaired by the cell, the 8-OHdG is released into the urine.

8-OHdG formation is considered to have a close relationship with tumors, cell aging, exposure to xenobiotics, cardiovascular disease, and other degenerative diseases. A meta-analysis of 8-OHdG in solid tumor cancers found that in most solid tumors, except breast cancer, elevated 8-OHdG was associated with poorer survival outcomes.²⁴⁸ Some researched associations to elevated 8-OHdG are as follows:

↑ Causes of Increased 8-OHdG

- Increased cortisol²⁴⁹
- Smoking²⁵⁰
- Pollutants such as benzene, styrene, toluene, PAH, asbestos, BPA, heavy metals²⁵¹
- Cancer (numerous solid tumor types and leukemia)²⁵²
- Occupational radiation exposure²⁵³
- Atherosclerosis²⁵⁴
- Hypertension²⁵⁴
- Diabetes²⁵⁵
- COPD
- Recent, acute exercise²⁵⁶ (moderate exercise over time shows reductions in 8-OHdG²⁵⁷)

The Good News about Diet and 8-OHdG

The good news is that a variety of food-based interventions with antioxidant effects have been shown to reduce 8-OHdG levels.

- Diets rich in colorful fruits, vegetables and polyphenols have been shown to reduce 8-OHdG. The inclusion of foods such as onions, almonds, green tea, lycopene, curcumin, and red wine have independently shown the ability to reduce 8-OHdG.¹¹⁸
- Simply increasing the number of fruits and vegetables per day in one study, i.e., up to 12 servings, has shown decreased 8-OHdG markers in study participants.²⁵⁸
- The Mediterranean diet is a cohesive diet that brings these elements together and has shown reductions in 8-OHdG in metabolic syndrome patients.¹¹⁸

Diurnal Values of Free Cortisol and Cortisone

Urinary Free Cortisol

Urinary Free Diurnal Cortisol	Free Cortisol (1st morning) Free Cortisol (2nd morning) Free Cortisol (Evening) Free Cortisol (Night) Free Cortisol (Pooled)
-------------------------------	--

Free cortisol is a measure of unbound cortisol. Serum testing typically measures the carrier or protein-bound version of hormone that is travelling through the blood. While free and albumin-bound hormones can enter tissues, only unbound hormones can exert actions on target cells. Typically, about 5-10% of cortisol is thought to be free and 90-95% bound to binding proteins.^{259,260}

Urine free cortisol has been shown to strongly correlate with mean serum-free cortisol in cortisol excess situations.²⁶¹ Measurement of free cortisol is helpful for determination of HPA axis function. It is also relevant for patients who are on drugs which affect binding proteins of cortisol, such as contraceptives, low dose prednisone and others.²⁶²

↑ Potential Causes of Elevated Timed Cortisol Samples

Much of the research on diurnal patterns of cortisol has been conducted on salivary cortisol, rather than urine, however, summaries of those findings are listed here.

↑ **Free Cortisol (1st morning)** – Elevated waking cortisol values can specifically reflect depression, early stages of burnout, pain, glycemic dysregulation, or job-related stress.²⁶³ Severe obstructive sleep apnea has been shown to produce irregular AM cortisol values in studies (both high and low).^{264, 265}

↑ **Free Cortisol (2nd morning)** – As a solo elevation, this generally reflects sporadic or situational stress triggers. One recent study showed that elevations specifically of noon or night cortisol accurately predicted job stress and/or burnout better than frequently studied cortisol awakening response (CAR) values (not measured with this test) which show conflicting results regarding burnout.²⁶⁶

↑ **Free Cortisol (Evening)** – In addition to situational stress triggers, high evening cortisol can be related to glycemic dysregulation due to late-day, prior to dinner collection time. In the literature it has also been associated with autism spectrum disorder, adolescence, and home related stressors such as divorce and financial strain.^{267, 268}

↑ **Free Cortisol (Night)** – This value reflects the expected nadir of cortisol output and therefore, baseline cortisol levels. Increased levels are related to insomnia, situational stressors, inflammation, and disease conditions such as the ones in the above table. Markedly increased bedtime cortisol values should prompt consideration of Cushing's syndrome or disease.²⁶⁹

↑ **Free Cortisol (Pooled)** - Pooled cortisol reflects overall diurnal cortisol output, collected from all samples given. Thus, it can be a measure of general HPA axis functionality. Some evidence shows that elevations in pooled cortisol reflect chronic stress.²⁶⁸

↓ Potential Causes of Decreased Timed Cortisol Samples

- ↓ **Free Cortisol (1st morning)** – Severe exhaustion can result in a decreased AM cortisol, in addition to severe obstructive sleep apnea and conditions of adrenal insufficiency such as Addison’s Disease.^{270, 265}
- ↓ **Free Cortisol (2nd morning)** – Generally, decreased values represent incomplete recovery from acute stress exposures, exhaustion from chronic stress and/or conditions related to low cortisol output.
- ↓ **Free Cortisol (Evening)** – Any of the conditions listed in the Causes of Decreased Cortisol Table listed earlier in this section can produce a decreased pooled cortisol value. Chronic stress has been linked to a flattened or attenuated diurnal output of cortisol, however the opposite has also been shown. This discrepancy has been theorized to be related to changes in HPA axis output related to long-term exposure to chronic stress.²⁷¹

Urinary Free Cortisol

Urinary Free Diurnal Cortisol	Free Cortisol (1st morning) Free Cortisol (2nd morning) Free Cortisol (Evening) Free Cortisol (Night) Free Cortisol (Pooled)
-------------------------------	--

Diurnal cortisone values have not been well studied for normative value comparisons to diurnal cortisol. Nonetheless, while urinary cortisone itself is an inactive metabolite, it has been shown, like urinary cortisol, to directly compare to serum free cortisol levels. Thus, elevations or depressions in urinary cortisone values tend to parallel cortisol values. However, differences between urinary free cortisone and cortisol can result from activity in 11-βHSD enzyme function. Evaluation of the ratio of cortisol/cortisone values and cortisol and cortisone metabolites can shed light on these enzymatic actions. Refer to the 11β-HSD Activity section under Total Cortisone for further information.

Interpretation of Cortisol Circadian Rhythm – HPA Axis Dysregulation

Our collective understanding of the stress response and the HPA axis largely comes from the pioneering work of Hans Selye. His theory of General Adaptation Syndrome described a limited, yet functional, model of stress response and adaptation to persistent stress. It is based on a stressor producing 3 stages of response.

Stage 1 – The first stage that one goes through in response to a stressor is termed “alarm” and indicates a release in catecholamines along with a transient increase in corticosteroids.

Stage 2 – The second stage of stress response is termed “resistance” and describes a heightened stress response with persistently elevated cortisol and other physiological adaptations necessary to resist the stressor.

Stage 3 – The third stage is termed “exhaustion” and results in persistent hypocortisolism and depletion of other biological resources needed to maintain the stress response.²⁷²

As research evolves it can be said that numerous factors influence the stress response and subsequent HPA axis result, and the issue is more complex than originally modeled. Selye and others have noted,²⁷³ timing is a critical element. Hormonal activity is elevated at stressor onset but reduces as time passes. However, how much time passes, and how much the cortisol response is affected, is widely variable among persons. This heterogeneity may relate to HPA axis linked genetic variation among other factors.²⁷¹ Second, ongoing stressors that threaten physical integrity, involve trauma, and are uncontrollable can elicit a high, flat diurnal profile of cortisol secretion.²⁷³ However, it is demonstrated that this is not the case in every situation. For example, it is lower in people with posttraumatic stress, and it may be lower after many years of persistent stress. Third, regarding hypocortisolemia, situations that provoke repeated or sustained cortisol elevations may lead to a breakdown in the negative feedback system of cortisol secretion, ultimately resulting in low flattened slopes.²⁷³

HPA Axis Function

A normal cortisol circadian rhythm is that cortisol is elevated in the morning, peaks within 30-45 minutes following waking, known as the cortisol awakening response or CAR, then begins a steady downward trend until it reaches its nadir at bedtime.²⁷¹

Questions to Ask in an HPA Axis Evaluation

1. Does the Diurnal rhythm peak in the morning and then trend downward from AM to PM?
2. Is there a single elevation only, that can be explained by a situational trigger or chronic condition?
3. Is there a combination of either high or low elevations with loss of diurnal rhythm?




Figure 1. Normal Variation of HPA Axis Curve

The graph above has values within an expected HPA axis configuration. Although a typical HPA axis curve shows all values trending downward from the morning cortisol value, in this test, the second morning cortisol, taken 2 hours after waking, is higher than the first morning collection. In urine testing this can be a normal variant. As urine captures hormone levels accumulated in the urine since the previous urine sample, it is possible that the increased hormone produced from the cortisol awakening response (CAR) has been captured in this 2nd morning sample. Clinical correlation would be helpful to rule in or out situational triggers as well. The evening and night values proceed in a typical downward trend with the expected nadir at the night sample.




Figure 2. Flattened Diurnal Slope

A flattened diurnal slope is found when there is a pattern of blunted response to stress throughout the day, typically with hypocortisolemia. This flattened diurnal slope pattern represents a more advanced state of HPA axis dysregulation. This pattern has been studied in relationship to chronic stress, early childhood adversity, burnout, etc. and has associations with many clinical conditions as seen on the right.²⁷⁴



Clinical Tip

A Flattened HPA Axis Curve is an Important Finding!

The most studied variation in diurnal cortisol slope pattern is a flattened or blunted slope. The following conditions have been associated with a flatter diurnal cortisol slope²⁷⁴

- Immune dysregulation
- Inflammatory dysregulation
- Chronic fatigue syndrome
- Chronic stress
- Breast cancer mortality
- Obesity/BMI/adiposity
- Depression and negative affect
- Cardiovascular disease
- Recent withdrawal from corticosteroid therapy

Lifestyle Treatment Considerations for Glucocorticoids

Hypercortisolism – Treatment Objectives

Consider which of the following single or combination of therapeutic approaches is the best fit for the patient with hypercortisolism. Transient elevated cortisol values may be corrected by remediation of situational triggers. Further modifiable HPA axis triggers for persistent hypercortisolism include dysglycemic and inflammatory conditions as listed below.

- Identify stressful situational triggers and offer lifestyle remediation long commutes, caffeine, jet lag, stressful life event, hypoglycemia, insomnia
- Identify and treat glycemic dysregulation diabetes, hypoglycemia, stress/comfort eating, eating disorder
- Identify and treat chronic inflammatory conditions which may promote hyper signaling of HPA axis such as: autoimmune disease, arthritis or chronic musculoskeletal pain, chronic sinusitis, CVD, GI inflammation

Table 11. Holistic Treatment Considerations for High Cortisol Effects

Lifestyle Considerations	To Aid Cortisol Reduction/Mitigate Effects of High Cortisol
Supplements	<p>Ashwagandha 300m/ BID KSM-66 brand (Chandrasekhar), also + studies for Sensoril and Shoden brands²⁷⁵</p> <p>Phosphatidyl serine/phosphatidic acid 400mg²⁷⁶</p> <p>L-theanine 200mg²⁷⁷ anti-stress and cortisol lowering effect</p> <p>Curcumin 500mg BID lowered AM cortisol²⁷⁸</p> <p>Curserin brand (Phosphatidyl serine/curcumin blend -800mg/BID) lowers AM cortisol²⁷⁹</p> <p>Panax Ginseng 1.5g/ Korean ginseng BID reduced cortisol in metabolic syndrome¹⁸⁰</p> <p>Rhodiola 300mg/BID²⁷⁵</p> <p>Relora Magnolia/phellodendron 250mg BID²⁸⁰</p> <p>DHEA* 10-450mg/d in meta-analysis of multiple studies reduces cortisol, considered safe under doctor monitoring (avg. dose 50 or 100mg), moderate's effects of high cortisol in target tissues²⁸¹</p>
Diet	<p>Low sucrose diet/carbohydrate restriction May modulates 11B-HSD enzyme activity to reduce truncal obesity²⁸²</p> <p>Prebiotics, high fiber diet may reduce cortisol levels, anxiety, and related cognitive function²⁸³</p> <p>Cherry juice 125 ml/BID²⁸⁴enhances mood, reduces cortisol</p> <p>In Cushing's syndrome and prolonged hypercortisolism consider the following dietary measures: Low sodium diet, (excess cortisol can lead to hypertension), higher calcium diet (to offset osteoporosis), higher protein diet (offset muscle wasting)</p>

Exercise	<p>Yoga heterogenous mix of styles, is associated with decreased evening and/or waking cortisol, improved regulation of HPA axis system²⁸⁵</p> <p>Habitual exercise (low or moderate intensity) regulates HPA axis function, buffers anxiety response and rumination time^{286,287} however extreme exercise is associated with prolonged cortisol elevations and decreased immunity</p> <p>Weight bearing exercise to offset bone loss from chronic hypercortisolism</p>
Stress Management	<p>Meditation Numerous styles of meditation are associated with decreases in salivary cortisol and anxiety, rumination^{288,289}</p> <p>Forest Bathing Decreases cortisol after forest walking, anticipation of stress relief is also helpful²⁹⁰</p> <p>Biofeedback 5 minutes once a day decreases stress and cortisol levels²⁹¹</p> <p>QiGong Reduces cortisol, raises serotonin and brain derived neurotropic factor (BDNF)²⁹²</p> <p>Sleep Maximize sleep – at least 7 hours each night for prevention of high cortisol</p>

*Not advised to use DHEA without physician supervision, monitoring cortisol values, and monitoring for excess androgenic side effects/symptoms

Further Testing Considerations for Elevated Cortisol

- **Cardiac Health Panel** and **Diabetes Panel**
- **Methylation Panel**
- **Tickborne** or Infections tests to explore occult infections
- **Gut Zoomer** for inflammation if cause of persistent inflammation is unknown or clinical symptoms correlate
- **hs-CRP** for inflammation
- **Thyroid Panel**
- **DEXA scan** for osteoporosis
- **Blood pressure** monitoring
- Extreme elevations of cortisol require a more **extensive endocrine workup** including workup for Cushing’s syndrome, adrenal tumor/hyperplasia and/or pituitary tumor etc. Endocrinologist referral is advised.

Hypocortisolism – Treatment Objectives

- Differentiate low to moderate hypocortisolism, which is more common, from adrenal insufficiency which is more uncommon.
- Hypocortisolism, sub-optimal adrenal output, and/or HPA axis dysregulation are all terms which indicate maladaptation to chronic stress and are associated with conditions such as burnout, PTSD, impacts of early life trauma, chronic fatigue syndrome, etc.
- Adrenal insufficiency (Primary, secondary, or tertiary) is extreme hypocortisolism; it can lead to life threatening events and necessitates glucocorticoid and possibly mineralocorticoid replacement therapy. Similarly, a blunted or flattened circadian release of cortisol should prompt a workup for adrenal insufficiency.
- Maximize sleep – at least 7 hours each night for restoration and repair
- Address chronic life/job stressors/major depression and provide relevant support
- Re-establish circadian rhythm through morning light, early, daytime exercise and sleeping in total darkness
- Consider nutrient and herbal support for adrenal nourishment and stress buffering
- Regular macronutrient balanced meals to stabilize blood glucose levels

Table 12. Holistic Treatment Considerations for Low Cortisol Effects

Lifestyle Considerations	Improve Low Cortisol/Re-establish Diurnal Rhythm
Supplements	<p>Phosphatidyl serine/omega 3 fatty acid blend 300mg normalization of blunted cortisol response²⁹³</p> <p>St John's Wort 600mg/d increased cortisol at 600mg/d in 2 studies, null cortisol effect in others at varying doses²⁷⁵</p> <p>Licorice Inhibits 11β-HSD2 – increases cortisol: cortisone ratio in target tissues. dosage varies based on glycyrrhizin content of product, genetics, gut flora, blood pressure. Licorice may promote pseudoaldosteronism effects if dose is too high – monitor BP²⁹⁴</p> <p>Adaptogens such as Rhodiola rosea, Eleutherococcus senticosus and Schisandra chinensis show benefits for fatigue, endurance and cognitive performance in several studies, however, no evidence found for increased cortisol with use²⁹⁵</p> <p>Hydrocortisone (Cortef) for adrenal insufficiency*</p>
Diet	Grapefruit Juice - inhibits 11β-HSD2 – can increase cortisol: cortisone ratio
Exercise	Conflicting studies - Low intensity exercise only, to patient tolerance
Stress Management	Mindfulness meditation 12-minute meditation scan reduces blunted cortisol in cancer patients ²⁹⁶

*Hydrocortisone/Cortef is a prescription glucocorticoid replacement used for adrenal insufficiency. Lower doses (2.5-10mg daily) are sometimes used by integrative physicians for low or moderate hypocortisolism or HPA axis dysregulation with blunted cortisol levels. Risk reward considerations, frequent cortisol/DHEAS monitoring, evaluation for side effects and established tapering programs are advised with any glucocorticoid replacement plan.

Further Testing Considerations for Low Cortisol

- **Cardiac Health Panel**
- **Neurotransmitter test**
- **Total Tox Burden** (Heavy Metals, Environmental Toxins, and Mycotoxins tests)
- **Thyroid Panel**
- Extreme hypocortisolism requires a more extensive endocrine workup including ACTH stimulation test, etc., for adrenal insufficiency. Endocrinologist referral is advised.

Urinary Free Melatonin - 6-Hydroxymelatonin Sulfate (MT6s)

Urinary Diurnal Melatonin Mt6s	Melatonin (1st morning) Melatonin (2nd morning) Melatonin (Evening) Melatonin (Night) Melatonin (pooled)
--------------------------------	--

Melatonin is a hormone synthesized from tryptophan and secreted by the pineal gland. Its principal role is in maintaining circadian rhythm and sleep induction; however, it also has antioxidant, anti-inflammatory, immunomodulatory, and anticancer effects.²⁹⁷ Disruption in melatonin secretion is positively correlated with certain clinical conditions such as major depressive disorder, metabolic syndrome, Alzheimer's disease, Parkinson's disease, fibromyalgia, hypertension, and cancer.²⁹⁸⁻³⁰¹

Melatonin is metabolized to 6-hydroxymelatonin, primarily in the liver, and then conjugated with sulfuric acid to form 6-sulfatoxymelatonin sulfate (MT6s). This is the primary melatonin metabolite found in the urine.³⁰² Urinary MT6s strongly correlates with plasma melatonin rhythm and is an established marker for evaluation of plasma melatonin.³⁰³

Melatonin is primarily secreted in response to environmental light cues. Melatonin production begins to increase about 2 hours prior to an established bedtime and increases throughout the sleep cycle. During peak darkness, the highest amount of melatonin is produced, typically around 3-4am in a typical sleep cycle. As light enters the room in the early daytime, melatonin levels drop off. Melatonin production is at its lowest at midday.




Figure 3. Diurnal Melatonin

This graph is an average result of melatonin urine results from 4x a day testing. The peak morning result is the collection of excreted overnight urine, which typically has the highest amount of melatonin. Epidemiological studies have shown a significant relationship between nocturnal plasma melatonin and morning urinary MT6 levels.³⁰²

Melatonin levels typically drop off at midday and results of that urine collection containing metabolized melatonin are seen in the evening sample. Melatonin metabolism studies have shown that there is a time lag seen between plasma melatonin values and urine values that is expected to be about 2 hours due to melatonin metabolism in the liver and subsequent renal clearance.³⁰⁴ This timing is affected by individual renal variability. Towards the end of day, melatonin levels rise again, and this initial rise is seen in the night urine sample.

↑ Potential Causes of Increased Melatonin

- Use of supplemental melatonin
- Use of SSRI antidepressants and MAO inhibitors or supplements such as 5-HTP or tryptophan
- Caffeine)³⁰⁵
- Oral contraceptives (mixed studies)³⁰⁵

↓ Potential Causes of Decreased Melatonin

- Age – In the 40's melatonin production begins to decline, and declines increase with age related calcification of the pineal gland
- Exposure to blue light from computer screens in the evening and electromagnetic frequencies
- Exercise at night can delay or reduce melatonin secretion
- Comorbid health conditions such as: renal disease, hepatic disease, ophthalmic disease, spinal cord disease
- Medications: β 1-adrenergic blockers (i.e., propranolol), α 2 adrenergic agonists (i.e., clonidine), NSAIDS, benzodiazepines^{302,305}
- Consider melatonin supplementation if clinically relevant

Urinary Creatinine




Figure 4. Diurnal Creatinine

Urine creatinine, a metabolite of creatine from muscle and protein metabolism, is often used as a biomarker of renal function and/or can reflect a number of systemic conditions. Creatinine elimination is stable, largely independent of diet and hydration fluctuations. Thus, creatinine values in the urine, in addition to use as a disease biomarker, can also be used to adjust for hydration status in urine lab testing.³⁰⁶ All the hormones tested on Vibrant Wellness Urine Hormones test are divided by urinary creatinine concentration values to improve accuracy of results which may be otherwise skewed by hydration variability. Individual creatinine concentrations are charted here for additional reference.

Disclaimer

The Urine Hormone test is a general wellness test. This interpretive guide is intended to offer information for improvement of functions associated with a general state of health while making references to diseases or conditions. The content in this guide is not meant to diagnose, treat, or cure any disease or condition. The clients who receive Vibrant Wellness Urine Hormone test results are advised to consult their physician, and/or health care provider team for diagnosis and further follow up care, including but not limited to additional testing, prescription medication, and any treatment interventions including diet, exercise, or lifestyle management.

References

1. Ziegler RG, Fuhrman BJ, Moore SC, Matthews CE. Epidemiologic studies of estrogen metabolism and breast cancer. *Steroids*. 2015;99(Pt A):67-75. doi:10.1016/j.steroids.2015.02.015
2. Files JA, Ko MG, Pruthi S. Bioidentical hormone therapy. *Mayo Clin Proc*. 2011;86(7):673-680. doi:10.4065/mcp.2010.0714
3. Wise PM, Suzuki S, Brown CM. Estradiol: A hormone with diverse and contradictory neuroprotective actions. *Dialogues Clin Neurosci*. 2009;11(3):297-303. doi:10.31887/dcns.2009.11.3/pmwise
4. Kumar A, Banerjee A, Singh D, et al. Estradiol: A Steroid with Multiple Facets. *Horm Metab Res*. 2018;50(5):359-374. doi:10.1055/s-0044-100920
5. Patel S, Homaei A, Raju AB, Meher BR. Estrogen: The necessary evil for human health, and ways to tame it. *Biomed Pharmacother*. 2018;102(March):403-411. doi:10.1016/j.biopha.2018.03.078
6. Parida S, Sharma D. The Microbiome-Estrogen Connection and Breast Cancer Risk. *Cells*. 2019;8(12). doi:10.3390/cells8121642
7. X. Chai, S. Zeng WX. Estrogen-Metabolizing Enzymes in Systemic and Local Liver Injuries: A Case Study of Disease-Drug Interaction. In: *Drug Metabolism in Diseases*. ; 2017:241-255.
8. Davis SR, Martinez-Garcia A, Robinson PJ, et al. Estrone is a strong predictor of circulating estradiol in women age 70 years and older. *J Clin Endocrinol Metab*. 2020;105(9):E3348-E3354. doi:10.1210/clinem/dgaa429
9. Williams DE. Indoles Derived From Glucobrassicin: Cancer Chemoprevention by Indole-3-Carbinol and 3,3'-Diindolylmethane. *Front Nutr*. 2021;8(October):1-14. doi:10.3389/fnut.2021.734334
10. Yang L, Zahid M, Liao Y, et al. Reduced formation of depurinating estrogen-DNA adducts by sulforaphane or KEAP1 disruption in human mammary epithelial MCF-10A cells. *Carcinogenesis*. 2013;34(11):2587- 2592. doi:10.1093/carcin/bgt246
11. Mason C, De Dieu Tapsoba J, Duggan C, Imayama I, Wang CY, Korde LA, Stanczyk F MA. Effects of vitamin D supplementation during weight loss on sex hormones in postmenopausal women. *Meno-pause*. 2016;Jun;23((6)):645-652.
12. Tao J, Li S, Gan RY, Zhao CN, Meng X, Li H Bin. Targeting gut microbiota with dietary components on cancer: Effects and potential mechanisms of action. *Crit Rev Food Sci Nutr*. 2020;60(6):1025-1037. doi:10.1080/10408398.2018.1555789
13. Balunas M, Su B, Brueggemeier R, Kinghorn A. Natural Products as Aromatase Inhibitors. *Anticancer Agents Med Chem*. 2013;8(6):646-682. doi:10.2174/187152008785133092
14. Zamir, A.; Ben-Zeev T., Hoffman J. Manipulation of Dietary Intake on Circulating Changes in Circulating- Testosterone Concentrations. *Nutrients*. 2021;13(3375). doi:10.3390/ nu13103375
15. El-Kersh DM, Ezzat SM, Salama MM, et al. Anti-estrogenic and anti-aromatase activities of citrus peels major compounds in breast cancer. *Sci Rep*. 2021;11(1):1-14. doi:10.1038/s41598-021-86599-z
16. Balam FH, Ahmadi ZS, Ghorbani A. Inhibitory effect of chrysin on estrogen biosynthesis by suppression of enzyme aromatase (CYP19): A systematic review. *Heliyon*. 2020;6(3):e03557. doi:10.1016/j.heliyon.2020. e03557
17. Yi B, Kasai H, Lee HS, Kang Y, Park JY YM. Inhibition by wheat sprout (*Triticum aestivum*) juice of bisphenol A-induced oxidative stress in young women. *Mutat Res*. 2011;Sep 18;724(1-2):64-68
18. Rock CL, Flatt SW, Thomson CA, et al. Effects of a high-fiber, low-fat diet intervention on serum concentrations of reproductive steroid hormones in women with a history of breast cancer. *J Clin Oncol*. 2004;22(12):2379-2387. doi:10.1200/JCO.2004.09.025

19. Calado A, Neves PM, Santos T, Ravasco P. The Effect of Flaxseed in Breast Cancer: A Literature Review. *Front Nutr.* 2018;5(February):1-7. doi:10.3389/fnut.2018.00004
20. Salehi B, Venditti A, Sharifi-Rad M, et al. The therapeutic potential of Apigenin. *Int J Mol Sci.* 2019;20(6). doi:10.3390/ijms20061305
21. Schmitz KH, Williams NI, Kontos D, et al. Dose-response effects of aerobic exercise on estrogen among women at high risk for breast cancer: a randomized controlled trial. *Breast Cancer Res Treat.* 2015;154(2):309-318. doi:10.1007/s10549-015-3604-z
22. Ticinesi A, Lauretani F, Tana C, Nouvenne A, Ridolo E, Meschi T. Exercise and immune system as modulators of intestinal microbiome: Implications for the gut-muscle axis hypothesis. *Exerc Immunol Rev.* 2019;25(96):84-95
23. Mnif W, Hassine AIH, Bouaziz A, Bartegi A, Thomas O, Roig B. Effect of endocrine disruptor pesticides: A review. *Int J Environ Res Public Health.* 2011;8(6):2265-2303. doi:10.3390/ijerph8062265
24. Colleluori G, Chen R, Turin CG, et al. Aromatase Inhibitors Plus Weight Loss Improves the Hormonal Profile of Obese Hypogonadal Men Without Causing Major Side Effects. *Front Endocrinol (Lausanne).* 2020;11(May):1-12. doi:10.3389/fendo.2020.00277
25. Morris PG, Hudis CA, Giri D, et al. Breast Tissue of Obese Women With Breast Cancer. *Cancer Prev Res (Phila).* 2011;4(7):1021-1029. doi:10.1158/1940-6207.CAPR-11-0110.
26. Saran S, Gupta BS, Philip R, Singh KS, Bende SA, Agroiya P, Agrawal P. Effect of hypothyroidism on female reproductive hormones. *Indian J Endocrinol Metab.* 2016 Jan-Feb;20(1):108-13. doi: 10.4103/2230-8210.172245. PMID: 26904478; PMCID: PMC4743370.
27. Osmer R, Friede M, Liske E, Schnitker J, Freudenstein J HZH. Efficacy and safety of isopropanolic black cohosh extract for climacteric symptoms. *Obs Gynecol.* 2005;May;105(5 Pt 1):1074-1083
28. Chen LR, Chen KH. Utilization of isoflavones in soybeans for women with menopausal syndrome: An overview. *Int J Mol Sci.* 2021;22(6):1-23. doi:10.3390/ijms22063212
29. Dording CM, Schettler PJ, Dalton ED, Parkin SR, Walker RS, Fehling KB, Fava M MD. A double-blind placebo-controlled trial of maca root as treatment for antidepressant-induced sexual dysfunction in women. *Evid Based Complement Altern Med.* 2015;949036.
30. Suwanvesh N, Manonai J, Sophonsritsuk A, Cherdshewasart W. Comparison of Pueraria mirifica gel and conjugated equine estrogen cream effects on vaginal health in postmenopausal women. *Menopause.* 2017;24(2):210-215. doi:10.1097/GME.0000000000000742
31. Kim JE, Jeong H, Hur S, Lee J, Kwon O. Efficacy and safety of kudzu flower–mandarin peel on hot flashes and bone markers in women during the menopausal transition: A randomized controlled trial. *Nutrients.* 2020;12(11):1-13. doi:10.3390/nu12113237
32. Myers SP, Vigar V. Effects of a standardised extract of *Trifolium pratense* (Promensil) at a dosage of 80 mg in the treatment of menopausal hot flashes: A systematic review and meta-analysis. *Phytomedicine.* 2017;24:141-147. doi:10.1016/j.phymed.2016.12.003
33. So SH, Lee JW, Kim YS, Hyun SH, Han CK. Red ginseng monograph. *J Ginseng Res.* 2018;42(4):549-561. doi:10.1016/j.jgr.2018.05.002
34. Gopal S, Ajgaonkar A, Kanchi P, Kaundinya A, Thakare V, Chauhan S LD. Effect of an ashwagandha (*Withania Somnifera*) root extract on climacteric symptoms in women during perimenopause: A randomized, double-blind, placebo-controlled study. *J Obs Gynaecol Res.* 2021;Dec;47(12):4414-4425
35. De Franciscis P, Colacurci N, Riemma G, et al. A nutraceutical approach to menopausal complaints. *Med.* 2019;55(9):16-18. doi:10.3390/medicina55090544
36. Steels E, Steele ML, Harold M CS. Efficacy of a Proprietary *Trigonella foenum-graecum* L. De-Husked Seed Extract in Reducing Menopausal Symptoms in Otherwise Healthy Women: A Double-Blind, Randomized, Placebo-Controlled Study. *Phytother Res.* 2017;Sep;31(9):1316-1322

37. Soleymani M, Siassi F, Qorbani M, et al. Dietary patterns and their association with menopausal symptoms: A cross-sectional study. *Menopause*. 2019;26(4):365-372. doi:10.1097/GME.0000000000001245
38. Pugliese G. Mediterranean diet as tool to manage obesity in menopause: A narrative review. *Nutrition*. 2020;Nov-Dec;79(110991)
39. Stojanovska L, Apostolopoulos V, Polman R BE. To exercise, or, not to exercise, during menopause and beyond. *Maturitas*. 2014;Apr;77(4):318-323
40. Goldstein KM, Coeytaux RR, Mcduffie JR, et al. Improving vasomotor symptoms; psychological symptoms; and health-related quality of life in peri- or post-menopausal women through yoga: An umbrella systematic review and meta_analysis. *Complement Ther Med*. Published online 2020:156-164. doi:10.1016/j.ctim.2017.08.011.
41. Guthrie KA, Larson JC, Ensrud KE, et al. Effects of pharmacologic and nonpharmacologic interventions on insomnia symptoms and self-reported sleep quality in women with hot flashes: A pooled analysis of individual participant data from four MsFLASH trials. *Sleep*. 2018;41(1). doi:10.1093/sleep/zsx190
42. Johnson A, Roberts L, Elkins G. Complementary and Alternative Medicine for Menopause. *J Evidence-Based Integr Med*. 2019;24:1-14. doi:10.1177/2515690X19829380
43. Carmody JF, Crawford S, Salmoirago-Blotcher E, Leung K, Churchill L ON. Mindfulness training for coping with hot flashes: results of a randomized trial. *Menopause*. 2011;Jun;18(6):611-620
44. Chiu HY, Pan CH, Shyu YK, Han BC TP. Effects of acupuncture on menopause-related symptoms and quality of life in women in natural menopause: a meta-analysis of randomized controlled trials. *Menopause*. 2015;Feb;22(2):234-244
45. Jeyakumar M, Carlson KE, Gunther JR KJ. Exploration of dimensions of estrogen potency: parsing ligand binding and coactivator binding affinities. *J Biol Chem*. 2011;286(15):12971-12982.
46. MA. JKA. "Estrogens." In: Meyler's Side Effects of Drugs. Elsevier Publishing; 2015:122-151
47. Reed MJ, Beranek PA, Ghilchik MW JV. Conversion of estrone to estradiol and estradiol to estrone in postmenopausal women. *Obs Gynecol*. 1985;Sept(66):361-5.
48. Karim R, Mack WJ, Hodis HN, Roy S SF. Influence of age and obesity on serum estradiol, estrone, and sex hormone binding globulin concentrations following oral estrogen administration in postmenopausal women. *J Clin Endocrinol Metab*. 2009;Nov;94(11):4136-43. doi:10.1210/jc.2009-0643
49. Jasonni VM, Bulletti C, Bolelli GF, Franceschetti F, Bonavia M, Ciotti P FC. Estrone sulfate, estrone and estradiol concentrations in normal and cirrhotic postmenopausal women. *Steroids*. 1983;May;41(5):569-73. doi:10.1016/0039-128x(83)90022-3
50. Southren L. E al. The Conversion of Androgens to Estrogens in Hyperthyroidism. *J Clin Endocrinol Metab*. 1974;38(2):207-214
51. Trabert B, Brinton LA E al. Circulating Estrogens and Postmenopausal Ovarian Cancer Risk in the Women's Health Initiative Observational Study. *Cancer Epidemiol Biomarkers Prev*. 2016;Apr;25(4):648-656. doi:10.1158/1055-9965
52. Cavlan D, Bharwani N GA. Androgen- and estrogen-secreting adrenal cancers. *Semin Oncol*.2010;-Dec;37(6):638-648
53. Ali ES, Mangold C, Peiris AN. Estriol: Emerging clinical benefits. *Menopause*. 2017;24(9):1081-1085. doi:10.1097/GME.0000000000000855
54. Falah N, Torday J, K. Quinney S, M. Haas D. Estriol review: Clinical applications and potential biomedical importance. *Clin Res Trials*. 2015;1(2):29-33. doi:10.15761/crt.1000109
55. Mucci LA, Lagiou P, Tamimi RM, Hsieh C, Adami H. Pregnancy estriol , estradiol , progesterone and prolactin in relation to birth weight and other birth size variables (United States). Published online 2003:311-318.

56. Minsart A, Onderbergen A Van, Jacques F, Kurt C, Gillerot Y. Indication of prenatal diagnosis in pregnancies complicated by undetectable second-trimester maternal serum estriol levels. Published online 2008:27-30.
57. Lemon HM, Wotiz HH, Parsons L, Mozden PJ. Reduced Estriol Excretion in Patients with Breast Cancer Prior to Endocrine Therapy. *JAMA J Am Med Assoc.* 1966;196(13):1128-1136. doi:10.1001/jama.1966.03100260066020
58. Yonemoto RH. Breast Cancer in Japan and the United States. *Arch Surg.* 1980;115(9):1056. doi:10.1001/archsurg.1980.01380090036009
59. Stoddard FR, Brooks AD, Eskin BA, Johannes GJ. Iodine alters gene expression in the MCF7 breast cancer cell line: Evidence for an anti-estrogen effect of iodine. *Int J Med Sci.* 2008;5(4):189-196. doi:10.7150/ijms.5.189
60. Manjer J, Sandsveden M, Borgquist S. Serum iodine and breast cancer risk: A prospective nested case-control study stratified for selenium levels. *Cancer Epidemiol Biomarkers Prev.* 2020;29(7):1335-1340. doi:10.1158/1055-9965.EPI-20-0122
61. Farebrother J, Zimmermann MB AM. Excess iodine intake: sources, assessment, and effects on thyroid function. *Ann N Y Acad Sci.* 2019;Jun(1446(1)):44-65
62. Sepkovic DW, Bradlow HL. Estrogen hydroxylation--the good and the bad. *Ann N Y Acad Sci.* 2009 Feb;1155:57-67. doi: 10.1111/j.1749-6632.2008.03675.x. PMID: 19250192.
63. Moore SC, Matthews CE, Ou Shu X, Yu K, Gail MH, Xu X, Ji BT, Chow WH, Cai Q, Li H, Yang G, Ruggieri D, Boyd-Morin J, Rothman N, Hoover RN, Gao YT, Zheng W, Ziegler RG. Endogenous Estrogens, Estrogen Metabolites, and Breast Cancer Risk in Postmenopausal Chinese Women. *J Natl Cancer Inst.* 2016 May 18;108(10):djw103. doi: 10.1093/jnci/djw103. PMID: 27193440; PMCID: PMC5858156.
64. Taioli E, Im A, Xu X, Veenstra TD, Ahrendt G, Garte S. Comparison of estrogens and estrogen metabolites in human breast tissue and urine. *Reprod Biol Endocrinol.* 2010 Aug 2;8:93. doi: 10.1186/1477-7827-8-93. PMID: 20678202; PMCID: PMC2922211.
65. Ziegler RG, Fuhrman BJ, Moore SC, Matthews CE. Epidemiologic studies of estrogen metabolism and breast cancer. *Steroids.* 2015 Jul;99(Pt A):67-75. doi: 10.1016/j.steroids.2015.02.015. Epub 2015 Feb 26. PMID: 25725255; PMCID: PMC5722219.
66. Yager JD. Mechanisms of estrogen carcinogenesis: The role of E2/E1-quinone metabolites suggests new approaches to preventive intervention--A review. *Steroids.* 2015 Jul;99(Pt A):56-60. doi: 10.1016/j.steroids.2014.08.006. Epub 2014 Aug 24. PMID: 25159108; PMCID: PMC4339663.
67. Cavalieri E, Rogan E. The 3,4-Quinones of Estrone and Estradiol Are the Initiators of Cancer whereas Resveratrol and N-acetylcysteine Are the Preventers. *International Journal of Molecular Sciences.* 2021; 22(15):8238. doi.org/10.3390/ijms22158238
68. Romilly E, Hodges, Deanna M, Minich, "Modulation of Metabolic Detoxification Pathways Using Foods and Food-Derived Components: A Scientific Review with Clinical Application", *Journal of Nutrition and Metabolism*, vol. 2015, Article ID 760689, 23 pages, 2015. doi.org/10.1155/2015/760689
69. Dalessandri KM, Firestone GL, Fitch MD, Bradlow HL, Bjeldanes LF. Pilot study: effect of 3,3'-diindolylmethane supplements on urinary hormone metabolites in postmenopausal women with a history of early-stage breast cancer. *Nutr Cancer.* 2004;50(2):161-7. doi: 10.1207/s15327914nc5002_5. PMID: 15623462.
70. Brooks JD, Ward WE, Lewis JE, Hilditch J, Nickell L, Wong E, Thompson LU. Supplementation with flaxseed alters estrogen metabolism in postmenopausal women to a greater extent than does supplementation with an equal amount of soy. *Am J Clin Nutr.* 2004 Feb;79(2):318-25. doi: 10.1093/ajcn/79.2.318. PMID: 14749240.
71. Chow HH, Garland LL, Heckman-Stoddard BM, Hsu CH, Butler VD, Cordova CA, Chew WM, Cornelison TL. A pilot clinical study of resveratrol in postmenopausal women with high body mass index: effects on systemic sex steroid hormones. *J Transl Med.* 2014 Aug 14;12:223. doi: 10.1186/s12967-014-0223-0. PMID: 25115686; PMCID: PMC4243716.

72. Sisti JS, Hankinson SE, Caporaso NE, Gu F, Tamimi RM, Rosner B, Xu X, Ziegler R, Eliassen AH. Caffeine, coffee, and tea intake and urinary estrogens and estrogen metabolites in premenopausal women. *Cancer Epidemiol Biomarkers Prev.* 2015 Aug;24(8):1174-83. doi: 10.1158/1055-9965.EPI-15-0246. Epub 2015 Jun 10. PMID: 26063478; PMCID: PMC4526325.
73. Lucille R. Marchand MD, BSN, James A. Stewart MD, FACP. "Breast Cancer". *Integrative Medicine (Fourth Edition)*, 2018
74. Zhu BT, Loder DP, Cai MX, Ho CT, Huang MT, Conney AH. Dietary administration of an extract from rosemary leaves enhances the liver microsomal metabolism of endogenous estrogens and decreases their uterotrophic action in CD-1 mice. *Carcinogenesis.* 1998 Oct;19(10):1821-7. doi: 10.1093/carcin/19.10.1821. PMID: 9806165.
75. Zahid M, Kohli E, Saeed M, Rogan E, Cavalieri E. The greater reactivity of estradiol-3,4-quinone vs estradiol-2,3-quinone with DNA in the formation of depurinating adducts: implications for tumor-initiating activity. *Chem Res Toxicol.* 2006 Jan;19(1):164-72. doi: 10.1021/tx050229y. PMID: 16411670.
76. Cavalieri, E.L., Rogan, E.G. and Zahid, M. (2017), Critical depurinating DNA adducts: Estrogen adducts in the etiology and prevention of cancer and dopamine adducts in the etiology and prevention of Parkinson's disease. *Int. J. Cancer*, 141: 1078-1090. <https://doi.org/10.1002/ijc.30728>
77. Miao S, Yang F, Wang Y, Shao C, Zava DT, Ding Q, Shi YE. 4-Hydroxy estrogen metabolite, causing genomic instability by attenuating the function of spindle-assembly checkpoint, can serve as a biomarker for breast cancer. *Am J Transl Res.* 2019 Aug 15;11(8):4992-5007. PMID: 31497216; PMCID: PMC6731443.
78. Romilly E. Hodges, Deanna M. Minich, "Modulation of Metabolic Detoxification Pathways Using Foods and Food-Derived Components: A Scientific Review with Clinical Application", *Journal of Nutrition and Metabolism*, vol. 2015, Article ID 760689, 23 pages, 2015. <https://doi.org/10.1155/2015/760689>
79. Li Yang, Muhammad Zahid, et al. Reduced formation of depurinating estrogen–DNA adducts by sulforaphane or KEAP1 disruption in human mammary epithelial MCF-10A cells, *Carcinogenesis*, Volume 34, Issue 11, November 2013, Pages 2587–2592, <https://doi.org/10.1093/carcin/bgt246>
80. Falero-Perez J, Song YS, Sorenson CM, Sheibani N. CYP1B1: A key regulator of redox homeostasis. *Trends Cell Mol Biol.* 2018;13:27-45. PMID: 30894785; PMCID: PMC6422340.
81. Liu X, Huang T, Li L, Tang Y, Tian Y, Wang S, Fan C. CYP1B1 deficiency ameliorates obesity and glucose intolerance induced by high fat diet in adult C57BL/6J mice. *Am J Transl Res.* 2015 Apr 15;7(4):761-71. PMID: 26064443; PMCID: PMC4455350.
82. Jacob A, Hartz AM, Potin S, Coumoul X, Yousif S, Scherrmann JM, Bauer B, Declèves X. Aryl hydrocarbon receptor-dependent upregulation of Cyp1b1 by TCDD and diesel exhaust particles in rat brain microvessels. *Fluids Barriers CNS.* 2011 Aug 25;8:23. doi: 10.1186/2045-8118-8-23. PMID: 21867498; PMCID: PMC3180256.
83. Alsubait A, Aldossary W, Rashid M, Algamdi A, Alrfaei BM. CYP1B1 gene: Implications in glaucoma and cancer. *J Cancer.* 2020 May 19;11(16):4652-4661. doi: 10.7150/jca.42669. PMID: 32626511; PMCID: PMC7330686.
84. Rodriguez-Melendez R, Griffin JB, Zemleni J. Biotin supplementation increases expression of the cytochrome P450 1B1 gene in Jurkat cells, increasing the occurrence of single-stranded DNA breaks. *J Nutr.* 2004 Sep;134(9):2222-8. doi: 10.1093/jn/134.9.2222. PMID: 15333708.
85. Shah BR, Xu W, Mraz J. Cytochrome P450 1B1: role in health and disease and effect of nutrition on its expression. *RSC Adv.* 2019 Jul 4;9(36):21050-21062. doi: 10.1039/c9ra03674a. PMID: 35515562; PMCID: PMC9065998.
86. Li F, Zhu W, Gonzalez FJ. Potential role of CYP1B1 in the development and treatment of metabolic diseases. *Pharmacol Ther.* 2017 Oct;178:18-30. doi: 10.1016/j.pharmthera.2017.03.007. Epub 2017 Mar 16. PMID: 28322972; PMCID: PMC5600638.

87. Meng X, An X, Zhou L, Fu B, Jia L. The isomers, aloe-emodin and emodin, possess differential inhibitory activities against CYP1B1 enzyme. *Steroids*. 2022 Sep;185:109055. doi: 10.1016/j.steroids.2022.109055. Epub 2022 Jun 2. PMID: 35661798.
88. Cavalieri EL, Rogan EG. Depurinating estrogen-DNA adducts, generators of cancer initiation: their minimization leads to cancer prevention. *Clin Transl Med*. 2016 Mar;5(1):12. doi: 10.1186/s40169-016-0088-3. Epub 2016 Mar 15. PMID: 26979321; PMCID: PMC4792821.
89. Li Yang, Muhammad Zahid, et al. Reduced formation of depurinating estrogen–DNA adducts by sulforaphane or KEAP1 disruption in human mammary epithelial MCF-10A cells, *Carcinogenesis*, Volume 34, Issue 11, November 2013, Pages 2587–2592, <https://doi.org/10.1093/carcin/bgt246>
90. Chow HH, Garland LL, Heckman-Stoddard BM, Hsu CH, Butler VD, Cordova CA, Chew WM, Cornelison TL. A pilot clinical study of resveratrol in postmenopausal women with high body mass index: effects on systemic sex steroid hormones. *J Transl Med*. 2014 Aug 14;12:223. doi: 10.1186/s12967-014-0223-0. PMID: 25115686; PMCID: PMC4243716.
91. Sisti JS, Hankinson SE, Caporaso NE, Gu F, Tamimi RM, Rosner B, Xu X, Ziegler R, Eliassen AH. Caffeine, coffee, and tea intake and urinary estrogens and estrogen metabolites in premenopausal women. *Cancer Epidemiol Biomarkers Prev*. 2015 Aug;24(8):1174-83. doi: 10.1158/1055-9965.EPI-15-0246. Epub 2015 Jun 10. PMID: 26063478; PMCID: PMC4526325.
92. Nicola Napoli, Sirisha Donepudi, Sharmin Sheikh, Giovam Batista Rini, Reina Armamento-Villareal, Increased 2-Hydroxylation of Estrogen in Women with a Family History of Osteoporosis, *The Journal of Clinical Endocrinology & Metabolism*, Volume 90, Issue 4, 1 April 2005, Pages 2035–2041, <https://doi.org/10.1210/jc.2004-1425>
93. Eliassen AH, Spiegelman D, Xu X, Keefer LK, Veenstra TD, Barbieri RL, Willett WC, Hankinson SE, Ziegler RG. Urinary estrogens and estrogen metabolites and subsequent risk of breast cancer among premenopausal women. *Cancer Res*. 2012 Feb 1;72(3):696-706. doi: 10.1158/0008-5472.CAN-11-2507. Epub 2011 Dec 5. PMID: 22144471; PMCID: PMC3271178.
94. Fuhrman BJ, Schairer C, Gail MH, et al. Estrogen metabolism and risk of breast cancer in postmenopausal women. *J Natl Cancer Inst*. 2012;104(4):326-339. doi:10.1093/jnci/djr531
95. Im A, Vogel VG, Ahrendt G, Lloyd S, Ragin C, Garte S, Taioli E. Urinary estrogen metabolites in women at high risk for breast cancer. *Carcinogenesis*. 2009 Sep;30(9):1532-5. doi: 10.1093/carcin/bgp139. Epub 2009 Jun 5. PMID: 19502596; PMCID: PMC2736301.
96. Dalessandri KM, Firestone GL, Fitch MD, Bradlow HL, Bjeldanes LF. Pilot study: effect of 3,3'-diindolylmethane supplements on urinary hormone metabolites in postmenopausal women with a history of early-stage breast cancer. *Nutr Cancer*. 2004;50(2):161-7. doi: 10.1207/s15327914nc5002_5. PMID: 15623462.
97. Williams DE. Indoles Derived From Glucobrassicin: Cancer Chemoprevention by Indole-3-Carbinol and 3,3'-Diindolylmethane. *Front Nutr*. 2021;8:734334. Published 2021 Oct 1. doi:10.3389/fnut.2021.734334
98. Brooks JD, Ward WE, Lewis JE, Hilditch J, Nickell L, Wong E, Thompson LU. Supplementation with flaxseed alters estrogen metabolism in postmenopausal women to a greater extent than does supplementation with an equal amount of soy. *Am J Clin Nutr*. 2004 Feb;79(2):318-25. doi: 10.1093/ajcn/79.2.318. PMID: 14749240.
99. Lu LJ, Cree M, Josyula S, Nagamani M, Grady JJ, Anderson KE. Increased urinary excretion of 2-hydroxyestrone but not 16 α -hydroxyestrone in premenopausal women during a soya diet containing isoflavones. *Cancer Res*. 2000;60:1299–305
100. Chow HH, Garland LL, Heckman-Stoddard BM, Hsu CH, Butler VD, Cordova CA, Chew WM, Cornelison TL. A pilot clinical study of resveratrol in postmenopausal women with high body mass index: effects on systemic sex steroid hormones. *J Transl Med*. 2014 Aug 14;12:223. doi: 10.1186/s12967-014-0223-0. PMID: 25115686; PMCID: PMC4243716.

101. Sisti JS, Hankinson SE, Caporaso NE, Gu F, Tamimi RM, Rosner B, Xu X, Ziegler R, Eliassen AH. Caffeine, coffee, and tea intake and urinary estrogens and estrogen metabolites in premenopausal women. *Cancer Epidemiol Biomarkers Prev.* 2015 Aug;24(8):1174-83. doi: 10.1158/1055-9965.EPI-15-0246. Epub 2015 Jun 10. PMID: 26063478; PMCID: PMC4526325.
102. Lucille R. Marchand MD, BSN, James A. Stewart MD, FACP. "Breast Cancer". *Integrative Medicine (Fourth Edition)*, 2018
103. Zhu BT, Loder DP, Cai MX, Ho CT, Huang MT, Conney AH. Dietary administration of an extract from rosemary leaves enhances the liver microsomal metabolism of endogenous estrogens and decreases their uterotrophic action in CD-1 mice. *Carcinogenesis.* 1998 Oct;19(10):1821-7. doi: 10.1093/carcin/19.10.1821. PMID: 9806165.
104. Miao S, Yang F, Wang Y, Shao C, Zava DT, Ding Q, Shi YE. 4-Hydroxy estrogen metabolite, causing genomic instability by attenuating the function of spindle-assembly checkpoint, can serve as a biomarker for breast cancer. *Am J Transl Res.* 2019 Aug 15;11(8):4992-5007. PMID: 31497216; PMCID: PMC6731443
105. Qin X, Peng Q, Qin A, Chen Z, Lin L, Deng Y, Xie L, Xu J, Li H, Li T, Li S, Zhao J. Association of COMT Val-158Met polymorphism and breast cancer risk: an updated meta-analysis. *Diagn Pathol.* 2012 Oct 8;7:136. doi: 10.1186/1746-1596-7-136. PMID: 23039364; PMCID: PMC3543196.
106. Lee JS, Kim YK, Yang H, Kang HY, Ahn C and Jeung EB: Two faces of the estrogen metabolite 2-methoxyestradiol in vitro and in vivo. *Mol Med Rep* 12: 5375-5382, 2015
107. Samavat H, Kurzer MS. Estrogen metabolism and breast cancer. *Cancer Lett.* 2015 Jan 28;356(2 Pt A):231-43. doi: 10.1016/j.canlet.2014.04.018. Epub 2014 Apr 28. PMID: 24784887; PMCID: PMC4505810.
108. Choi, H.J., Lee, A.J., Kang, K.S. et al. 4-Hydroxyestrone, an Endogenous Estrogen Metabolite, Can Strongly Protect Neuronal Cells Against Oxidative Damage. *Sci Rep* 10, 7283 (2020). <https://doi.org/10.1038/s41598-020-62984-y>
109. Romilly E. Hodges, Deanna M. Minich, "Modulation of Metabolic Detoxification Pathways Using Foods and Food-Derived Components: A Scientific Review with Clinical Application", *Journal of Nutrition and Metabolism*, vol. 2015, Article ID 760689, 23 pages, 2015. <https://doi.org/10.1155/2015/760689>
110. Busserolles J, Zimowska W, Rock E, Rayssiguier Y, Mazur A. Rats fed a high sucrose diet have altered heart antioxidant enzyme activity and gene expression. *Life Sci.* 2002;71(11):1303-1312. doi:10.1016/s0024-3205(02)01846-5
111. Hodges RE, Minich DM. Modulation of Metabolic Detoxification Pathways Using Foods and Food-Derived Components: A Scientific Review with Clinical Application. *J Nutr Metab.* 2015;2015:760689. doi:10.1155/2015/760689
112. Zeisel SH, da Costa KA. Choline: an essential nutrient for public health. *Nutr Rev.* 2009;67(11):615-623. doi:10.1111/j.1753-4887.2009.00246.x
113. Samavat H, Kurzer MS. Estrogen metabolism and breast cancer. *Cancer Lett.* 2015 Jan 28;356(2 Pt A):231-43. doi: 10.1016/j.canlet.2014.04.018. Epub 2014 Apr 28. PMID: 24784887; PMCID: PMC4505810.
114. Choi, H.J., Lee, A.J., Kang, K.S. et al. 4-Hydroxyestrone, an Endogenous Estrogen Metabolite, Can Strongly Protect Neuronal Cells Against Oxidative Damage. *Sci Rep* 10, 7283 (2020). <https://doi.org/10.1038/s41598-020-62984-y>
115. Brinton LA, Trabert B, Anderson GL, Falk RT, Felix AS, Fuhrman BJ, Gass ML, Kuller LH, Pfeiffer RM, Rohan TE, Strickler HD, Xu X, Wentzensen N. Serum Estrogens and Estrogen Metabolites and Endometrial Cancer Risk among Postmenopausal Women. *Cancer Epidemiol Biomarkers Prev.* 2016 Jul;25(7):1081-9. doi: 10.1158/1055-9965.EPI-16-0225. Epub 2016 Apr 12. PMID: 27197275; PMCID: PMC4930692.
116. Wan MLY, Co VA, El-Nezami H. Endocrine disrupting chemicals and breast cancer: a systematic review of epidemiological studies. *Crit Rev Food Sci Nutr.* 2022;62(24):6549-6576. doi: 10.1080/10408398.2021.1903382. Epub 2021 Apr 5. PMID: 33819127.

117. Bao W, Liu B, Rong S, Dai SY, Trasande L, Lehmler HJ. Association Between Bisphenol A Exposure and Risk of All-Cause and Cause-Specific Mortality in US Adults. *JAMA Netw Open*. 2020 Aug 3;3(8):e2011620. doi: 10.1001/jamanetworkopen.2020.11620. PMID: 32804211; PMCID: PMC7431989.
118. Crinnion W, Pizzorno J. *Clinical Environmental Medicine: Identification and Natural Treatment of Diseases Caused by Common Pollutants*. 1st Edition. St. Louis. Elsevier. 2018.
119. National Institute of Environmental Health Sciences. 2022. Bisphenol A (BPA). [online] Available at: <<https://www.niehs.nih.gov/health/topics/agents/syabpa/index.cfm>> (Accessed 1 October 2022)
120. Amjad S, Rahman MS, Pang MG. Role of Antioxidants in Alleviating Bisphenol A Toxicity. *Biomolecules*. 2020;10(8):1105. Published 2020 Jul 25. doi:10.3390/biom10081105
121. Genuis SJ, Beesoon S, Birkholz D, Lobo RA. Human excretion of bisphenol A: blood, urine, and sweat (BUS) study. *J Environ Public Health*. 2012;2012:185731. doi:10.1155/2012/185731
122. Yi B, Kasai H, Lee HS, Kang Y, Park JY, Yang M. Inhibition by wheat sprout (*Triticum aestivum*) juice of bisphenol A-induced oxidative stress in young women. *Mutat Res*. 2011;724(1-2):64-68. doi:10.1016/j.mrgentox.2011.06.007
123. Dušková M, Kolátorová L, Stárka L. Androgens in women - critical evaluation of the methods for their determination in diagnostics of endocrine disorders. *Physiol Res*. 2018;67:S379-S390. doi:10.33549/physiol-res.933964
124. Karbowska J, Kochan Z. Effects of DHEA on metabolic and endocrine functions of adipose tissue. *Horm Mol Biol Clin Investig*. 2013;14(2):65-74. doi:10.1515/hmbci-2013-0009
125. Matulevicius V, Urbanavicius V, Lukosevicius S, Banisaukaite I, Donielaite G, Galkine A. Importance of dehydroepiandrosterone sulfate assessment with special attention for adrenal tumours and arterial hypertension. *Acta Endocrinol (Copenh)*. 2021;17(1):68-76. doi:10.4183/aeb.2021.68
126. Chernecky C. DHEAS. In: *Laboratory Tests and Diagnostic Procedures*. 6th ed. Saunders; 2013:432-459.
127. Idkowiak J, Taylor AE, Subtil S, et al. Steroid Sulfatase Deficiency and Androgen Activation Before and After Puberty. *J Clin Endocrinol Metab*. 2016;101(6):2545-2553. doi:10.1210/jc.2015-4101
128. Chernecky C. DHEAS. In: *Laboratory Tests and Diagnostic Procedures*. 6th ed. Saunders; 2013:432-459.
129. Straub RH, Lehle K, Herfarth H, et al. Dehydroepiandrosterone in relation to other adrenal hormones during an acute inflammatory stressful disease state compared with chronic inflammatory disease: Role of interleukin-6 and tumour necrosis factor. *Eur J Endocrinol*. 2002;146(3):365-274. doi:10.1530/eje.0.1460365
130. Hindmarsh P, Geertsma K. Other Hormones and Their Roles. *Congenital Adrenal Hyperplasia*. Academic Press; 2017; Pages 309-325. Accessed 2/2023. doi.org/10.1016/B978-0-12-811483-4.00025-8
131. Sherbet, D.P., Auchus, R.J. (2007). Peripheral Testosterone Metabolism. In: Payne, A.H., Hardy, M.P. (eds) *The Leydig Cell in Health and Disease*. Contemporary Endocrinology. Humana Press. https://doi.org/10.1007/978-1-59745-453-7_12
132. Witchel SF. Congenital Adrenal Hyperplasia. *J Pediatr Adolesc Gynecol*. 2017 Oct;30(5):520-534. doi: 10.1016/j.jpag.2017.04.001. Epub 2017 Apr 24. PMID: 28450075; PMCID: PMC5624825.
133. Zhang XL, Zhang CW, Xu P, Liang FJ, Che YN, Xia YJ, Cao YX, Wu XK, Wang WJ, Yi L, Gao Q, Wang Y. SNP rs2470152 in CYP19 is correlated to aromatase activity in Chinese polycystic ovary syndrome patients. *Mol Med Rep*. 2012 Jan;5(1):245-9. doi: 10.3892/mmr.2011.616. Epub 2011 Oct 3. PMID: 21972004.
134. Kaminski RM, Marini H, Kim WJ, Rogawski MA. Anticonvulsant activity of androsterone and etiocholalone. *Epilepsia*. 2005 Jun;46(6):819-27. doi: 10.1111/j.1528-1167.2005.00705.x. PMID: 15946323; PMCID: PMC1181535.
135. Poór V, Juricskay S, Telegdy E. Urinary steroids in men with male-pattern alopecia. *J Biochem Biophys Methods*. 2002 Oct-Nov;53(1-3):123-30. doi: 10.1016/s0165-022x(02)00100-8. PMID: 12406594.

136. Thevis M, Geyer H, Mareck U, Flenker U, Schänzer W. Doping-control analysis of the 5 α -reductase inhibitor finasteride: determination of its influence on urinary steroid profiles and detection of its major urinary metabolite. *Ther Drug Monit.* 2007 Apr;29(2):236-47. doi: 10.1097/FTD.0b013e31803bb85d. PMID: 17417080.
137. Hampl R, Stárka L, Janský L. Steroids and thermogenesis. *Physiol Res.* 2006;55(2):123-131. doi: 10.33549/physiolres.930758. Epub 2005 May 24. PMID: 15910167.
138. Poór V, Juricskay S, Telegdy E. Urinary steroids in men with male-pattern alopecia. *J Biochem Biophys Methods.* 2002 Oct-Nov;53(1-3):123-30. doi: 10.1016/s0165-022x(02)00100-8. PMID: 12406594.
139. Caruso D, Abbiati F, Giatti S, Romano S, Fusco L, Cavaletti G, Melcangi RC. Patients treated for male pattern hair with finasteride show, after discontinuation of the drug, altered levels of neuroactive steroids in cerebrospinal fluid and plasma. *J Steroid Biochem Mol Biol.* 2015 Feb;146:74-9. doi: 10.1016/j.jsbmb.2014.03.012. Epub 2014 Apr 6. PMID: 24717976.
140. Melcangi RC, Santi D, Spezzano R, Grimoldi M, Tabacchi T, Fusco ML, Diviccaro S, Giatti S, Carrà G, Caruso D, Simoni M, Cavaletti G. Neuroactive steroid levels and psychiatric and andrological features in post-finasteride patients. *J Steroid Biochem Mol Biol.* 2017 Jul;171:229-235. doi: 10.1016/j.jsbmb.2017.04.003. Epub 2017 Apr 10. PMID: 28408350.
141. Evron E, Juhasz M, Babadjouni A, Mesinkovska NA. Natural Hair Supplement: Friend or Foe? Saw Palmetto, a Systematic Review in Alopecia. *Skin Appendage Disord.* 2020;6(6):329-337. doi:10.1159/000509905
142. Remer, Thomas; Maser-Gluth, Christiane; Wudy, Stefan (2008). Glucocorticoid Measurements in Health and Disease - Metabolic Implications and the Potential of 24-h Urine Analyses. *Mini-Reviews in Medicinal Chemistry*, 8(2), 153–170. doi:10.2174/138955708783498096
143. Akalestou E, Genser L, Rutter GA. Glucocorticoid Metabolism in Obesity and Following Weight Loss. *Front Endocrinol (Lausanne).* 2020;11:59. Published 2020 Feb 20. doi:10.3389/fendo.2020.00059
144. Evans BA, Griffiths K, Morton MS. Inhibition of 5 α -reductase in genital skin fibroblasts and prostate tissue by dietary lignans and isoflavonoids. *Journal of Endocrinology.* 1995;147:295–302.
145. Azizi A, Mumin NH, Shafqat N. Phytochemicals With Anti 5-alpha-reductase Activity: A Prospective For Prostate Cancer Treatment. *F1000Res.* 2021;10:221. Published 2021 Jul 6. doi:10.12688/f1000research.51066.3
146. Pais P, Villar A, Rull S. Determination of the potency of a novel saw palmetto supercritical CO₂ extract (SPSE) for 5 α -reductase isoform II inhibition using a cell-free in vitro test system. *Res Rep Urol.* 2016;8:41-49. Published 2016 Apr 21. doi:10.2147/RRU.S96576
147. Cicero AFG, Allkanjari O, Busetto GM, et al. Nutraceutical treatment and prevention of benign prostatic hyperplasia and prostate cancer. *Arch Ital Urol Androl.* 2019;91(3):10.4081/aiua.2019.3.139. Published 2019 Oct 2. doi:10.4081/aiua.2019.3.139
148. Grant P, Ramasamy S. An update on plant derived anti-androgens. *Int J Endocrinol Metab.* 2012;10(2):497-502. doi:10.5812/ijem.3644
149. Hosking AM, Juhasz M, Atanaskova Mesinkovska N. Complementary and Alternative Treatments for Alopecia: A Comprehensive Review. *Skin Appendage Disord.* 2019;5(2):72-89. doi:10.1159/000492035
150. Dhariwala, Maria Yusuf, and Padmini Ravikumar. "An overview of herbal alternatives in androgenetic alopecia." *Journal of cosmetic dermatology* 18.4 (2019): 966-975.
151. Bianchi VE, Bresciani E, Meanti R, Rizzi L, Omeljaniuk RJ, Torsello A. The role of androgens in women's health and wellbeing. *Pharmacol Res.* 2021;171(June):105758. doi:10.1016/j.phrs.2021.105758
152. Jameson JL. Testosterone. In: *Endocrinology: Adult and Pediatric.* 17th ed. Saunders; 2016
153. Burger N, Johnson J. Androgen production in women. *Androg Reprod Aging.* 2005;77(4):1-4. doi:10.1201/b14632-2
154. Clifton S, Macdowall W, Copas AJ, et al. Salivary testosterone levels and health status in men and women

- in the British general population: Findings from the third national survey of sexual attitudes and lifestyles (natsal-3). *J Clin Endocrinol Metab.* 2016;101(11):3939-3951. doi:10.1210/jc.2016-1669
155. Chernecky CC. Testosterone. In: *Laboratory Tests and Diagnostic Procedures*,. 6th ed. Saunders;2013
 156. Crawford M, Kennedy L. Testosterone replacement therapy: Role of pituitary and thyroid in diagnosis and treatment. *Transl Androl Urol.* 2016;5(6):850-858. doi:10.21037/tau.2016.09.01
 157. Nemati M, Nemati S, Taheri AM HB. Comparison of metformin and N-acetyl cysteine, as an adjuvant to clomiphene citrate, in clomiphene-resistant women with polycystic ovary syndrome. *J Gynecol Obs Hum Reprod.* 2017;Sep;46(7):579-585
 158. Menichini D, Facchinetti F. Effects of vitamin D supplementation in women with polycystic ovary syndrome: a review. *Gynecol Endocrinol.* 2020;36(1):1-5. doi:10.1080/09513590.2019.1625881
 159. Unfer V, Facchinetti F, Orrù B, Giordani B, Nestler J. Myo-inositol effects in women with PCOS: A meta-analysis of randomized controlled trials. *Endocr Connect.* 2017;6(8):647-658. doi:10.1530/EC-17-0243
 160. Pizzo A, Laganà AS, Barbaro L. Comparison between effects of myo-inositol and d-chiro-inositol on ovarian function and metabolic factors in women with PCOS. *Gynecol Endocrinol.* 2014;30(3):205-208. doi:10.3109/09513590.2013.860120
 161. Jamilian M AZ. The Effects of Soy Isoflavones on Metabolic Status of Patients With Polycystic Ovary Syndrome. *J Clin Endocrinol Metab.* 2016;Sep;101((9)):3386-3394
 162. Li L, Li C, Pan P, et al. A single arm pilot study of effects of berberine on the menstrual pattern, ovulation rate, hormonal and metabolic profiles in anovulatory Chinese women with polycystic ovary syndrome. *PLoS One.* 2015;10(12):1-11. doi:10.1371/journal.pone.0144072
 163. Piotrowska A, Pilch W, Czerwińska-Ledwig O, et al. The Possibilities of Using Chromium Salts as an Agent Supporting Treatment of Polycystic Ovary Syndrome. *Biol Trace Elem Res.* 2019;192(2):91-97. doi:10.1007/s12011-019-1654-5
 164. Szczuko M, Kikut J, Szczuko U, et al. Nutrition strategy and life style in polycystic ovary syndrome—narrative review. *Nutrients.* 2021;13(7):1-18. doi:10.3390/nu13072452
 165. Thomson CA, Chow HHS, Wertheim BC, et al. Breast Cancer Biomarker Modulation in Patients Taking Tamoxifen. 2018;165(1):97-107. doi:10.1007/s10549-017-4292-7.
 166. Haffner SM, Katz MS, Stern MP DJ. The relationship of sex hormones to hyperinsulinemia and hyperglycemia. *Metabolism.* 1988;Jul;37(7):683-688
 167. Shang Y, Zhou H, Hu M FH. Effect of Diet on Insulin Resistance in Polycystic Ovary Syndrome. . *J Clin Endocrinol Metab.* 2020;Oct 1;105(10):dgaa425
 168. Kalgaonkar S, Almario RU, Gurusinghe D, et al. Differential effects of walnuts vs almonds on improving metabolic and endocrine parameters in PCOS. *Eur J Clin Nutr.* 2011;65(3):386-393. doi:10.1038/ejcn.2010.266
 169. Yurtdaş G AY. A New Approach to Polycystic Ovary Syndrome: The Gut Microbiota. *J Am Coll Nutr.* 2020;May-Jun;39(4):371-382
 170. Haidari F, Banaei-Jahromi N, Zakerkish M, Ahmadi K. The effects of flaxseed supplementation on metabolic status in women with polycystic ovary syndrome: A randomized open-labeled controlled clinical trial. *Nutr J.* 2020;19(1):1-11. doi:10.1186/s12937-020-0524-5
 171. Nagata C, Kabuto M SH. Association of coffee, green tea, and caffeine intakes with serum concentrations of estradiol and sex hormone-binding globulin in premenopausal Japanese women. *Nutr Cancer.* 1998;30(1):21-24
 172. Pino AM, Valladares LE, Palma MA, Mancilla AM, Yáñez M AC. Dietary isoflavones affect sex hormone-binding globulin levels in postmenopausal women. *J Clin Endocrinol Metab.* 2000;Aug;85(8):2797-2800
 173. Sáez-López C, Soriguer F, Hernandez C, Rojo-Martinez G, Rubio-Martín E, Simó R SD. Oleic acid increases

- hepatic sex hormone binding globulin production in men. *Mol Nutr Food Res*. 2014;Apr;58(4):760-767
174. Haffner SM, Katz MS, Stern MP DJ. The relationship of sex hormones to hyperinsulinemia and hyperglycemia. *Metabolism*. 1988;Jul;37(7):683-688
175. Friedenreich CM, Woolcott CG, McTiernan A, et al. Alberta physical activity and breast cancer prevention trial: Sex hormone changes in a year-long exercise intervention among postmenopausal women. *J Clin Oncol*. 2010;28(9):1458-1466. doi:10.1200/JCO.2009.24.9557
176. Smith SJ, Lopresti AL, Teo SYM, Fairchild TJ. Examining the Effects of Herbs on Testosterone Concentrations in Men: A Systematic Review. *Adv Nutr*. 2021;12(3):744-765. doi:10.1093/advances/nmaa134
177. Wankhede S, Langade D, Joshi K, Sinha SR, Bhattacharyya S. Examining the effect of *Withania somnifera* supplementation on muscle strength and recovery: A randomized controlled trial. *J Int Soc Sports Nutr*. 2015;12(1):1-11. doi:10.1186/s12970-015-0104-9
178. Santos HO, Howell S, Teixeira FJ. Beyond tribulus (*Tribulus terrestris* L.): The effects of phytotherapies on testosterone, sperm and prostate parameters. *J Ethnopharmacol*. 2019;235:392-405. doi:10.1016/j.jep.2019.02.033
179. Shukla KK, Mahdi AA, Ahmad MK, Shankhwar SN, Rajender S, Jaiswar SP. *Mucuna pruriens* improves male fertility by its action on the hypothalamus-pituitary-gonadal axis. *Fertil Steril*. 2009;92(6):1934-1940. doi:10.1016/j.fertnstert.2008.09.045
180. Jung DH, Lee YJ, Kim CB, Kim JY, Shin SH PJ. Effects of ginseng on peripheral blood mitochondrial DNA copy number and hormones in men with metabolic syndrome: A randomized clinical and pilot study. *Complement Ther Med*. 2016;Feb(24):40-46
181. Li Y, Ren J, Li N, Liu J, Tan SC, Low TY MZ. A dose-response and meta-analysis of dehydroepiandrosterone (DHEA) supplementation on testosterone levels: perinatal prediction of randomized clinical trials. *Exp Gerontol*. 2020;Nov;141(111110)
182. El-Kersh DM, Ezzat SM, Salama MM, et al. Anti-estrogenic and anti-aromatase activities of citrus peels major compounds in breast cancer. *Sci Rep*. 2021;11(1):1-14. doi:10.1038/s41598-021-86599-z
183. Lo EM, Rodriguez KM, Pastuszak AW, Khara M. Alternatives to Testosterone Therapy: A Review. *Sex Med Rev*. 2018;6(1):106-113. doi:10.1016/j.sxmr.2017.09.004
184. Salehi B, Venditti A, Sharifi-Rad M, et al. The therapeutic potential of Apigenin. *Int J Mol Sci*. 2019;20(6). doi:10.3390/ijms20061305
185. Colleluori G, Chen R, Turin CG, et al. Aromatase Inhibitors Plus Weight Loss Improves the Hormonal Profile of Obese Hypogonadal Men Without Causing Major Side Effects. *Front Endocrinol (Lausanne)*. 2020;11(May):1-12. doi:10.3389/fendo.2020.00277
186. Raheem OA, Chen T, Akula KP, Greenberg J, Le TV, Chernobylsky D, Sikka SC WT. Efficacy of Non-Testosterone-Based Treatment in Hypogonadal Men: A Review. *Sex Med Rev*. 2021;Jul;9(3):381-392
187. Winter AG, Zhao F, Lee RK. Androgen deficiency and metabolic syndrome in men. *Transl Androl Urol*. 2014;3(1):50-58. doi:10.3978/j.issn.2223-4683.2014.01.04
188. Gasser B, Kurz J, Mohaupt M. Testosterone/Epitestosterone Ratios-Further Hints to Explain Hyperandrogenemia in Children with Autism. *Diseases*. 2021 Feb 1;9(1):13. doi: 10.3390/diseases9010013. PMID: 33535392; PMCID: PMC7931062.
189. Kicman AT, Oftebro H, Walker C, Norman N, Cowan DA. Potential use of ketoconazole in a dynamic endocrine test to differentiate between biological outliers and testosterone use by athletes. *Clin Chem*. 1993 Sep;39(9):1798-803. PMID: 8375050.
190. Stanczyk FZ. Diagnosis of hyperandrogenism: biochemical criteria. *Best Pract Res Clin Endocrinol Metab*. 2006 Jun;20(2):177-91. doi: 10.1016/j.beem.2006.03.007. PMID: 16772150.
191. Hammes SR, Levin ER. Impact of estrogens in males and androgens in females. *J Clin Invest*. 2019 May

- 1;129(5):1818-1826. doi: 10.1172/JCI125755. Epub 2019 May 1. PMID: 31042159; PMCID: PMC6486327.
192. Stanczyk FZ. Diagnosis of hyperandrogenism: biochemical criteria. *Best Pract Res Clin Endocrinol Metab.* 2006 Jun;20(2):177-91. doi: 10.1016/j.beem.2006.03.007. PMID: 16772150.
193. Dhayat NA, Marti N, Kollmann Z, Troendle A, Bally L, Escher G, Grössl M, Ackermann D, Ponte B, Pruijm M, Müller M, Vogt B, Birkhäuser MH, Bochud M, Flück CE; members of the SKIPOGH Study Group. Urinary steroid profiling in women hints at a diagnostic signature of the polycystic ovary syndrome: A pilot study considering neglected steroid metabolites. *PLoS One.* 2018 Oct 11;13(10):e0203903. doi: 10.1371/journal.pone.0203903. PMID: 30308019; PMCID: PMC6181287.
194. Frye CA, Edinger KL, Lephart ED, Walf AA. 3alpha-androstanediol, but not testosterone, attenuates age-related decrements in cognitive, anxiety, and depressive behavior of male rats. *Front Aging Neurosci.* 2010 Apr 8;2:15. doi: 10.3389/fnagi.2010.00015. PMID: 20552051; PMCID: PMC2874398.
195. Trabert B, Sherman ME, Kannan N, Stanczyk FZ. Progesterone and Breast Cancer. *Endocr Rev.* 2020 Apr 1;41(2):320-44. doi: 10.1210/edrev/bnz001. PMID: 31512725; PMCID: PMC7156851.
196. Herrera AY, Nielsen SE, Mather M. Stress-induced increases in progesterone and cortisol in naturally cycling women. *Neurobiol Stress.* 2016;3:96-104. doi:10.1016/j.yjnstr.2016.02.006
197. Kotsopoulos J, Eliassen AH, Missmer SA, Hankinson SE, Tworoger SS. Relationship between caffeine intake and plasma sex hormone concentrations in premenopausal and postmenopausal women. *Cancer.* 2009;115(12):2765-2774. doi:10.1002/cncr.24328
198. B.Q S, A.M DA-Z, S.S DS. Effects of Some Non Steroidal Anti-inflammatory Drugs on Ovulation in Women with Mild Musculoskeletal Pain (A Clinical Study). *IOSR J Pharm Biol Sci.* 2014;9(4):43-49. doi:10.9790/3008-09444349
199. Fleischman DS, Navarrete CD, Fessler DMT. Oral contraceptives suppress ovarian hormone production. *Psychol Sci.* 2010;21(5):750-752. doi:10.1177/0956797610368062
200. Piazza MJ, Urbanetz AA. Environmental toxins and the impact of other endocrine disrupting chemicals in women's reproductive health. *J Bras Reprod Assist.* 2019;23(2):154-164. doi:10.5935/1518-0557.20190016
201. Martelli M, Zingaretti L, Salvio G, Bracci M, Santarelli L. Influence of work on andropause and menopause: A systematic review. *Int J Environ Res Public Health.* 2021;18(19):1-20. doi:10.3390/ijerph181910074
202. Van Die MD, Burger H, Teede H, Bone K. Vitex agnus-castus extracts for female reproductive disorders: A systematic review of clinical trials. *Planta Med.* 2013;79(7):562-575. doi:10.1055/s-0032-1327831
203. Henmi H, Endo T, Kitajima Y, Manase K, Hata H, Kudo R. Effects of ascorbic acid supplementation on serum progesterone levels in patients with a luteal phase defect [2]. *Fertil Steril.* 2003;80(2):459-461. doi:10.1016/S0015-0282(03)00657-5
204. Abraham GE. Nutritional factors in the etiology of the premenstrual tension syndromes. *J Reprod Med.* 1983;Jul;28(7):446-464
205. Kashani L, Akhondzadeh S. Herbal medicine in the treatment of polycystic ovary syndrome. *J Med Plants.* 2016;15(59):1-5
206. Atteia HH, Alzahrani S, El-Sherbeeney NA, et al. Evening Primrose Oil Ameliorates Hyperleptinemia and Reproductive Hormone Disturbances in Obese Female Rats: Impact on Estrus Cyclicity. *Front Endocrinol* 2020;10(January). doi:10.3389/fendo.2019.00942
207. Chantilis S, Dombroski R, Shackleton CH, Casey ML, MacDonald PC. Metabolism of 5 alpha-dihydroprogesterone in women and men: 3 beta- and 3 alpha-, 6 alpha-dihydroxy-5 alpha-pregnan-20-ones are major urinary metabolites. *J Clin Endocrinol Metab.* 1996 Oct;81(10):3644-9. doi: 10.1210/jcem.81.10.8855816. PMID: 8855816.
208. Zamora-Sánchez CJ, Camacho-Arroyo I. Allopregnanolone: Metabolism, Mechanisms of Action, and Its Role in Cancer. *Int J Mol Sci.* 2022 Dec 29;24(1):560. doi: 10.3390/ijms24010560. PMID: 36614002; PMCID: PMC9820109.
209. Kimball A, Dichtel LE, Nyer MB, et al. The allopregnanolone to progesterone ratio across the menstrual cy-

- cle and in menopause. *Psychoneuroendocrinology*. 2020;112:104512. doi:10.1016/j.psyneuen.2019.104512
210. Ben Dor R, Marx CE, Shampine LJ, Rubinow DR, Schmidt PJ. DHEA metabolism to the neurosteroid androsterone: a possible mechanism of DHEA's antidepressant action. *Psychopharmacology (Berl)*. 2015 Sep;232(18):3375-83. doi: 10.1007/s00213-015-3991-1. Epub 2015 Jun 25. Erratum in: *Psychopharmacology (Berl)*. 2015 Oct;232(19):3683. PMID: 26105109; PMCID: PMC6309885.
211. Zamora-Sánchez CJ, Camacho-Arroyo I. Allopregnanolone: Metabolism, Mechanisms of Action, and Its Role in Cancer. *Int J Mol Sci*. 2022 Dec 29;24(1):560. doi: 10.3390/ijms24010560. PMID: 36614002; PMCID: PMC9820109.
212. Hantsoo L, Epperson CN. Allopregnanolone in premenstrual dysphoric disorder (PMDD): Evidence for dysregulated sensitivity to GABA-A receptor modulating neuroactive steroids across the menstrual cycle. *Neurobiol Stress*. 2020;12:100213. Published 2020 Feb 4. doi:10.1016/j.yinstr.2020.100213
213. Trabert B, Bauer DC, Buist DSM, Cauley JA, Falk RT, Geczik AM, Gierach GL, Hada M, Hue TF, Lacey JV Jr, LaCroix AZ, Tice JA, Xu X, Dallal CM, Brinton LA. Association of Circulating Progesterone With Breast Cancer Risk Among Postmenopausal Women. *JAMA Netw Open*. 2020 Apr 1;3(4):e203645. doi: 10.1001/jamanetworkopen.2020.3645. PMID: 32329771; PMCID: PMC7182797.
214. Wiebe JP. Progesterone metabolites in breast cancer. *Endocr Relat Cancer*. 2006 Sep;13(3):717-38. doi: 10.1677/erc.1.01010. PMID: 16954427.
215. Wiebe JP, Zhang G, Welch I, Cadieux-Pitre HA. Progesterone metabolites regulate induction, growth, and suppression of estrogen- and progesterone receptor-negative human breast cell tumors. *Breast Cancer Res*. 2013 May 11;15(3):R38. doi: 10.1186/bcr3422. PMID: 25927181; PMCID: PMC3706910.
216. Zamora-Sánchez CJ, Camacho-Arroyo I. Allopregnanolone: Metabolism, Mechanisms of Action, and Its Role in Cancer. *Int J Mol Sci*. 2022 Dec 29;24(1):560. doi: 10.3390/ijms24010560. PMID: 36614002; PMCID: PMC9820109.
217. Trabert B, Bauer DC, Buist DSM, et al. Association of Circulating Progesterone with Breast Cancer Risk Among Postmenopausal Women. *JAMA Netw Open*. 2020 Apr 1;3(4):e203645. doi: 10.1001/jamanetworkopen.2020.3645. PMID: 32329771; PMCID: PMC7182797.
218. Beranič, N., Gobec, S., & Rižner, T. L. (2011). Progestins as inhibitors of the human 20-ketosteroid reductases, AKR1C1 and AKR1C3. *Chemico-Biological Interactions*, 191(1–3), 227–233. <https://doi.org/10.1016/J.CBI.2010.12.012>
219. Traven, K., Sinreih, M., Stojan, J., Seršen, S., Kljun, J., Bezenšek, J., ... Rižner, T. L. (2015). Ruthenium complexes as inhibitors of the aldo-keto reductases AKR1C1–1C3. *Chemico-Biological Interactions*, 234, 349–359. <https://doi.org/10.1016/J.CBI.2014.11.005>
220. Lee, J., Hanley, J., Hopkins V. *What Your Doctor May Not Tell You about Menopause: The Breakthrough Book on Natural Hormone Balance*. Grand Central Publishing; 2004.
221. Grant LK, Gooley JJ, St Hilaire MA, et al. Menstrual phase-dependent differences in neurobehavioral performance: The role of temperature and the progesterone/estradiol ratio. *Sleep*. 2020;43(2):1-11. doi:10.1093/sleep/zsz227
222. Shalom-Paz E, Aslih N, Samara N, Michaeli M, Ellenbogen A. Late follicular progesterone to estradiol ratio is not influenced by protocols or gonadotropins used. *Reprod Biol Endocrinol*. 2015;13(1):1-7. doi:10.1186/s12958-015-0116-y
223. Thau L, Gandhi J SS. *Physiology, Cortisol*. Stat Pearls Publishing; 2021. https://www.ncbi.nlm.nih.gov/books/NBK538239/#_NBK538239_pubdet_
224. Chernecky B. Cortisol. In: *Laboratory Tests and Diagnostic Procedures*. 6th ed. Saunders; 2013
225. Cservenka A, Lahanas S, Dotson-Bossert J. Marijuana use and hypothalamic pituitary-adrenal axis functioning in humans. *Front Psychiatry*. 2018;9(OCT). doi:10.3389/fpsy.2018.00472
226. Aloisi AM, Aurilio C, Bachiocco V, Biasi G FP. Endocrine consequences of opioid therapy. *Psychoneuroendocrinology*. Dec; 34(Suppl 1:S):162-168.

227. Cservenka A, Lahanas S, Dotson-Bossert J. Marijuana use and hypothalamic pituitary-adrenal axis functioning in humans. *Front Psychiatry*. 2018;9(OCT). doi:10.3389/fpsy.2018.00472
228. Blair J, Adaway J, Keevil B, Ross R. Salivary cortisol and cortisone in the clinical setting. *Curr Opin Endocrinol Diabetes Obes*. 2017;24(3):161-168. doi:10.1097/MED.0000000000000328
229. Zhou C, Ye F, Wu H, Ye H, Chen Q. Recent advances in the study of 11 β -Hydroxysteroid dehydrogenase type 2 (11 β -HSD2) inhibitors. *Environ Toxicol Pharmacol*. 2017;52:47-53. doi:10.1016/j.etap.2017.02.021
230. Seckl JR. 11 β -hydroxysteroid dehydrogenases: Changing glucocorticoid action. *Curr Opin Pharmacol*. 2004;4(6):597-602. doi:10.1016/j.coph.2004.09.001
231. Bakusic J, De Nys S, Creta M, Godderis L, Duca RC. Study of temporal variability of salivary cortisol and cortisone by LC-MS/MS using a new atmospheric pressure ionization source. *Sci Rep*. 2019;9(1):1-12. doi:10.1038/s41598-019-55571-3
232. Bae YJ, Reinelt J, Netto J, et al. Salivary cortisone, as a biomarker for psychosocial stress, is associated with state anxiety and heart rate. *Psychoneuroendocrinology*. 2019;101(October 2018):35-41. doi:10.1016/j.psyneuen.2018.10.015
233. Chapman K, Holmes M, Seckl J. 11 β -Hydroxysteroid Dehydrogenases Intracellular Gate-Keepers of Tissue Glucocorticoid Action. *Physiol Rev*. 2013;93(3):1139-1206. doi:10.1152/physrev.00020.2012
234. Gant CM, Minovic I, Binnenmars H, et al. Lower Renal Function Is Associated With Derangement of 11- β Hydroxysteroid Dehydrogenase in Type 2 Diabetes. *J Endocr Soc*. 2018;2(7):609-620. Published 2018 May 22. doi:10.1210/js.2018-00088
235. Dube S, Norby BJ, Pattan V, Carter RE, Basu A BR. 11 β -hydroxysteroid dehydrogenase types 1 and 2 activity in subcutaneous adipose tissue in humans: implications in obesity and diabetes. *J Clin Endocrinol Metab*. 2015;Jan;100(1):E70-6.
236. Methlie P, Husebye EE, Hustad S, Lien EA LK. Grapefruit juice and licorice increase cortisol availability in patients with Addison's disease. *Eur J Endocrinol*. 2011;Nov;165(5):761-769
237. Shearer FJG, Wyrwoll CS, Holmes MC. The Role of 11 β -Hydroxy Steroid Dehydrogenase Type 2 in Glucocorticoid Programming of Affective and Cognitive Behaviours. *Neuroendocrinology*. 2019;109(3):257-265. doi:10.1159/000499660
238. Vantighem MC, Marcelli-Tourvieille S, Defrance F WJ. 11beta-hydroxysteroid dehydrogenases. Recent advances. *Ann Endocrinol*. 2007;Oct;68(5):349-356
239. Hintzpeter J, Stapelfeld C, Loerz C, Martin HJ, Maser E. Green tea and one of its constituents, epigallocatechin-3-gallate, are potent inhibitors of human 11 β -hydroxysteroid dehydrogenase type 1. *PLoS One*. 2014;9(1):2-10. doi:10.1371/journal.pone.0084468
240. Hu GX, Lin H, Lian QQ, et al. Curcumin as a Potent and Selective Inhibitor of 11 β -Hydroxysteroid Dehydrogenase 1: Improving Lipid Profiles in High-Fat-Diet-Treated Rats. *PLoS One*. 2013;8(3):1-7. doi:10.1371/journal.pone.0049976
241. Marissal-Arvy N, Hamiani R, Richard E, Moisan MP P V. Vitamin A regulates hypothalamic-pituitary-adrenal axis status in LOU/C rats. *J Endocrinol*. 2013;Sep 6;219(1):21-27.
242. Cuzzola A, Mazzini F, Petri A. A comprehensive study for the validation of a LC-MS/MS method for the determination of free and total forms of urinary cortisol and its metabolites. *J Pharm Biomed Anal*. 2014 Jun;94:203-9. doi: 10.1016/j.jpba.2014.02.003. Epub 2014 Feb 16. PMID: 24631839.
243. Taniyama M, Honma K, Ban Y. Urinary cortisol metabolites in the assessment of peripheral thyroid hormone action: application for diagnosis of resistance to thyroid hormone. *Thyroid*. 1993 Fall;3(3):229-33. doi: 10.1089/thy.1993.3.229. PMID: 8257864.
244. Angus Cameron, David Henley, Robin Carrell, Aiwu Zhou, Anthony Clarke, Stafford Lightman, Temperature-Responsive Release of Cortisol from Its Binding Globulin: A Protein Thermocouple, *The Journal of*

245. Asla Q, Sardà H, Lerma E, et al. 11-Deoxycorticosterone Producing Adrenal Hyperplasia as a Very Unusual Cause of Endocrine Hypertension: Case Report and Systematic Review of the Literature. *Front Endocrinol (Lausanne)*. 2022;13:846865. Published 2022 Mar 31. doi:10.3389/fendo.2022.846865
246. Raubenheimer, Peter J.; Young, Elizabeth A.; Andrew, Ruth; Seckl, Jonathan R. (2006). "The role of corticosterone in human hypothalamic– pituitary–adrenal axis feedback." *Clinical Endocrinology* 65(1): 22-26
247. Dlamini SN, Lombard Z, Micklesfield LK, Crowther N, Norris SA, Snyman T, Crawford AA, Walker BR, Goedecke JH. Glucocorticoids associate with cardiometabolic risk factors in black South Africans. *Endocr Connect*. 2021 Aug 2;10(8):873-884. doi: 10.1530/EC-21-0195. PMID: 34261039; PMCID: PMC8346194.
248. Qing, X., Shi, D., Lv, X. et al. Prognostic significance of 8-hydroxy-2'-deoxyguanosine in solid tumors: a meta-analysis. *BMC Cancer* 19, 997 (2019). <https://doi.org/10.1186/s12885-019-6189-9>
249. Joergensen A, Broedbaek K, Weimann A, Semba RD, Ferrucci L, Joergensen MB, Poulsen HE. Association between urinary excretion of cortisol and markers of oxidatively damaged DNA and RNA in humans. *PLoS One*. 2011;6(6):e20795. doi: 10.1371/journal.pone.0020795. Epub 2011 Jun 7. PMID: 21687734; PMCID: PMC3110199.
250. Valavanidis A, Vlachogianni T, Fiotakis C. 8-hydroxy-2' -deoxyguanosine (8-OHdG): A critical biomarker of oxidative stress and carcinogenesis. *J Environ Sci Health C Environ Carcinog Ecotoxicol Rev*. 2009 Apr;27(2):120-39. doi: 10.1080/10590500902885684. PMID: 19412858.
251. Li J, Lu S, Liu G, Zhou Y, Lv Y, She J, Fan R. Co-exposure to polycyclic aromatic hydrocarbons, benzene and toluene and their dose-effects on oxidative stress damage in kindergarten-aged children in Guangzhou, China. *Sci Total Environ*. 2015 Aug 15;524-525:74-80. doi: 10.1016/j.scitotenv.2015.04.020. Epub 2015 Apr 15. PMID: 25889546.
252. Jelic MD, Mandic AD, Maricic SM, Srdjenovic BU. Oxidative stress and its role in cancer. *J Cancer Res Ther*. 2021 Jan-Mar;17(1):22-28. doi: 10.4103/jcrt.JCRT_862_16. PMID: 33723127.
253. Gao Y, Wang P, Wang Z, et al. Serum 8-Hydroxy-2'-Deoxyguanosine Level as a Potential Biomarker of Oxidative DNA Damage Induced by Ionizing Radiation in Human Peripheral Blood. *Dose-Response*. 2019;17(1). doi:10.1177/1559325818820649
254. Di Minno A, Turnu L, Porro B, Squellerio I, Cavalca V, Tremoli E, Di Minno MN. 8-Hydroxy-2-Deoxyguanosine Levels and Cardiovascular Disease: A Systematic Review and Meta-Analysis of the Literature. *Antioxid Redox Signal*. 2016 Apr 1;24(10):548-55. doi: 10.1089/ars.2015.6508. PMID: 26650622; PMCID: PMC4827317.
255. Wang XB, Cui NH, Liu X, Liu X. Mitochondrial 8-hydroxy-2'-deoxyguanosine and coronary artery disease in patients with type 2 diabetes mellitus. *Cardiovasc Diabetol*. 2020 Feb 19;19(1):22. doi: 10.1186/s12933-020-00998-6. PMID: 32075646; PMCID: PMC7029479.
256. Zare M, Shateri Z, Nouri M, Sarbakhsh P, Eftekhari MH, Pourghassem Gargari B. Association between urinary levels of 8-hydroxy-2-deoxyguanosine and F2a-isoprostane in male football players and healthy non-athlete controls with dietary inflammatory and antioxidant indices. *Front Nutr*. 2023 Jan 24;9:1101532. doi: 10.3389/fnut.2022.1101532. PMID: 36761985; PMCID: PMC9902597.
257. Narukawa, Teruaki, et al. "Effects of Body Mass Index (BMI), Dietary Intake and Serum Antioxidant Vitamin Concentration on Urinary 8-hydroxydeoxyguanosine and F2-isoprostane Excretions." *ANTI-AGING MEDICINE* 8.1 (2011): 1-6.
258. Thompson HJ, Heimendinger J, Haegele A, et al. Effect of increased vegetable and fruit consumption on markers of oxidative cellular damage. *Carcinogenesis*. 1999;20(12):2261-2266. doi:10.1093/carcin/20.12.2261
259. Talbott SM, Talbott JA, Pugh M. Effect of *Magnolia officinalis* and *Phellodendron amurense* on cortisol and psychological mood state in moderately stressed subjects. *J Int Soc Sports Nutr*. 2013;10(1):1. doi:10.1186/1550-2783-10-37

260. Newman M, Curran DA, Mayfield BP. Dried urine and salivary profiling for complete assessment of cortisol and cortisol metabolites. *J Clin Transl Endocrinol.* 2020;22:100243. Published 2020 Nov 27. doi:10.1016/j.jcte.2020.100243
261. El-Farhan N, Rees DA, Evans C. Measuring cortisol in serum, urine and saliva – are our assays good enough? *Annals of Clinical Biochemistry.* 2017;54(3):308-322. doi:10.1177/0004563216687335
262. Imrich, Richard; Rovenský, Jozef (2010). Hypothalamic-Pituitary-Adrenal Axis in Rheumatoid Arthritis. 36(4), 721–727. doi:10.1016/j.rdc.2010.09.003
263. Klein M, Weksler N, Gidron Y, Heldman E, Gurski E, Smith OR GG. Do waking salivary cortisol levels correlate with anesthesiologist’s job involvement? *J Clin Monit Comput.* 2012;Dec;26(6):407-413
264. Mohammadi H, Rezaei M, Sharafkhaneh A, Khazaie H GM. Serum testosterone_cortisol ratio in people with obstructive sleep apnea _ Enhanced Reader.pdf. *J Clin Lab Anal.* 2020;Jan;34(1):e23011
265. Bozic J, Galic T, Supe-Domic D, Ivkovic N, Ticinovic Kurir T, Valic Z, Lesko J DZ. Morning cortisol levels and glucose metabolism parameters in moderate and severe obstructive sleep apnea patients. *Endocrine.* 2016;Sep;53.(3):730-739
266. Pilger A, Haslacher H, Meyer BM, et al. Midday and nadir salivary cortisol appear superior to cortisol awakening response in burnout assessment and monitoring. *Sci Rep.* 2018;8(1):1-12. doi:10.1038/s41598-018-27386-1
267. Muscatello R et. a. Comparing the effects of age, pubertal development, and symptom profile on cortisol rhythm in children and adolescents with Autism Spectrum Disorder. *Autism Res.* 2018;January;(1):110-120. doi:10.1002/aur.1879
268. Saxbe DE. A field (researcher’s) guide to cortisol: tracking HPA axis functioning in everyday life. *Health Psychol Rev.* 2008;2(2):163-190. doi:10.1080/17437190802530812
269. Papanicolaou DA, Mullen N, Kyrou I, Nieman LK. Nighttime salivary cortisol: A useful test for the diagnosis of Cushing’s syndrome. *J Clin Endocrinol Metab.* 2002;87(10):4515-4521. doi:10.1210/jc.2002-020534
270. Kim YJ, Kim JH, Hong AR, et al. Stimulated Salivary Cortisol as a Noninvasive Diagnostic Tool for Adrenal Insufficiency. *Endocrinol Metab.* 2020;35(3):628-635. doi:10.3803/enm.2020.707
271. Starr LR, Dienes K, Li YI, Shaw ZA. Chronic stress exposure, diurnal cortisol slope, and implications for mood and fatigue: Moderation by multilocus HPA-Axis genetic variation. *Psychoneuroendocrinology.* 2019;100:156-163. doi:10.1016/j.psyneuen.2018.10.003
272. Selye H. The general adaptation syndrome and the diseases of adaptation. *J Allergy.* 1946;17:231
273. Miller GE, Chen E, Zhou ES. If it goes up, must it come down? Chronic stress and the hypothalamic-pituitary-adrenocortical axis in humans. *Psychol Bull.* 2007;133(1):25-45. doi:10.1037/0033-2909.133.1.25
274. Tomas C, Newton J, Watson S. A Review of Hypothalamic-Pituitary-Adrenal Axis Function in Chronic Fatigue Syndrome. *ISRN Neurosci.* 2013;2013:1-8. doi:10.1155/2013/784520
275. Lopresti AL, Drummond PD, Smith SJ. A Randomized, Double-Blind, Placebo-Controlled, Crossover Study Examining the Hormonal and Vitality Effects of Ashwagandha (*Withania somnifera*) in Aging, Overweight Males. *Am J Mens Health.* 2019;13(2). doi:10.1177/1557988319835985
276. Hellhammer J, Hero T, Franz N, Contreras C, Schubert M. Omega-3 fatty acids administered in phosphatidylserine improved certain aspects of high chronic stress in men. *Nutr Res.* 2012;32(4):241-250. doi:10.1016/j.nutres.2012.03.003
277. White DJ, de Klerk S, Woods W, Gondalia S, Noonan C, Scholey AB. Anti-stress, behavioural and magnetoencephalography effects of an l-theanine-based nutrient drink: A randomised, double-blind, placebo-controlled, crossover trial. *Nutrients.* 2016;8(1). doi:10.3390/nu8010053
278. Yu JJ, Pei LB, Zhang Y, Wen ZY YJ. Chronic supplementation of curcumin enhances the efficacy of antide-

- pressants in Major Depressive Disorder: a randomized, double-blind, placebo-controlled pilot study. *J Clin Psychopharmacol.* 2015; Aug;35(4):406-410
279. Cicero AFG, Sahebkar A, Fogacci F, Bove M, Giovannini M BC. Effects of phytosomal curcumin on anthropometric parameters, insulin resistance, cortisolemia and non-alcoholic fatty liver disease indices: a double-blind, placebo-controlled clinical trial. *Eur J Nutr.* 2020;Mar;59(2):477-483
280. Talbott SM, Talbott JA, Pugh M. Effect of *Magnolia officinalis* and *Phellodendron amurense* (Relora) on cortisol and psychological mood state in moderately stressed subjects. *J Int Soc Sports Nutr.* 2013;10(1):1. doi:10.1186/1550-2783-10-37
281. Chen H, Jin Z, Sun C, Santos HO, Kord Varkaneh H. Effects of dehydroepiandrosterone (DHEA) supplementation on cortisol, leptin, adiponectin, and liver enzyme levels: A systematic review and meta-analysis of randomised clinical trials. *Int J Clin Pract.* 2021;75(11):1-11. doi:10.1111/ijcp.14698
282. London E CT. Diet and the role of 11beta-hydroxysteroid dehydrogenase-1 on obesity. *J Nutr Biochem.* Jul;20(7):485-493
283. Schmidt K, Cowen PJ, Harmer CJ, Tzortzis G, Errington S, Burnet PWJ. Prebiotic intake reduces the waking cortisol response and alters emotional bias in healthy volunteers. *Psychopharmacology (Berl).* 2015;232(10):1793-1801. doi:10.1007/s00213-014-3810-0
284. Garrido M, Espino J, Gonzalez-Gomez D, Lozano M, Barriga C, Paredes SD et al. The consumption of a Jerte Valley cherry product in humans enhances mood, and increases 5-hydroxyindoleacetic acid but reduces cortisol levels in urine. *Exp Gerontol.* 2012;Aug;47((8)):573-580
285. Pascoe M. A Systematic Review of Randomised Control Trials Examining the Effects of Mindfulness on Stress and Anxious Symptomatology. 2016;(March):1-23
286. Puterman, E., O; Donovan A. Physical activity moderates stressor-induced rumination on cortisol reactivity. *Psychosom Med.* 2011;73(7):604-611. doi:10.1097/PSY.0b013e318229e1e0
287. Wunsch K, Wurst R, von Dawans B, Strahler J, Kasten N FR. Habitual and acute exercise effects on salivary biomarkers in response to psychosocial stress. *Psychoneuroendocrinology.* 2019;Aug(106):216-225
288. Pascoe M. A Systematic Review of Randomised Control Trials Examining the Effects of Mindfulness on Stress and Anxious Symptomatology. 2016;(March):1-23
289. Sanada K, Montero-Marin J, Díez MA, et al. Effects of mindfulness-based interventions on salivary cortisol in healthy adults: A meta-analytical review. *Front Physiol.* 2016;7(OCT). doi:10.3389/fphys.2016.00471
290. Antonelli, M., Barbieri, G. & Donelli D. Effects of forest bathing (shinrin-yoku) on levels of cortisol as a stress biomarker: a systematic review and meta-analysis. *Int J Biometeorol.* Published online 2019:1117–1134. doi:/10.1007/s00484-019-01717-x
291. Kotozaki Y, Takeuchi H, Sekiguchi A, et al. Biofeedback-based training for stress management in daily hassles: An intervention study. *Brain Behav.* 2014;4(4):566-579. doi:10.1002/brb3.241
292. Lu EY, Lee P, Cai S, So WWY, Ng BFL, Jensen MP, Cheung WM TH. Qigong for the treatment of depressive symptoms: Preliminary evidence of neurobiological mechanisms. *Int J Geriatr Psychiatry.* 2020;Nov;35(11):1393-1401
293. Hellhammer J, Vogt D, Franz N, Freitas U, Rutenberg D. A soy-based phosphatidylserine/ phosphatidic acid complex (PAS) normalizes the stress reactivity of hypothalamus-pituitary-adrenal-axis in chronically stressed male subjects: A randomized, placebo-controlled study. *Lipids Health Dis.* 2014;13(1):1-11. doi:10.1186/1476-511X-13-121
294. Nazari S, Rameshrad M HH. Toxicological Effects of *Glycyrrhiza glabra* (Licorice): A Review. . *Phytother Res.* 2017;Nov;31(11):1635-1650
295. Panossian A, Wikman G. Effects of adaptogens on the central nervous system and the molecular mechanisms associated with their stress - Protective activity. *Pharmaceuticals.* 2010;3(1):188-224. doi:10.3390/ph3010188

296. Black D, Peng C, Sleight A. Mindfulness practice reduces cortisol blunting during chemotherapy: A randomized controlled study of colorectal cancer patients. *Cancer*. 2017;August 15(16):3088-3096. doi:10.1002/cncr.30698.
297. Ohki S, Kunimatsu M, Ogawa S, Takano H, Furihata T, Shibasaki H, Yokokawa A. Development and Validation of an LC-MS/MS-Based Method for Quantifying Urinary Endogenous 6-Hydroxymelatonin. *Chem Pharm Bull (Tokyo)*. 2022;70(5):375-382. doi: 10.1248/cpb.c21-00982. PMID: 35491194.
298. Zaki NF, Sabri YM, Farouk O, et al. Depressive Symptoms, Sleep Profiles and Serum Melatonin Levels in a Sample of Breast Cancer Patients. *Nat Sci Sleep*. 2020;12:135-149. Published 2020 Feb 13. doi:10.2147/NSS.S206768
299. Emet M, Ozcan H, Ozel L, Yayla M, Halici Z, Hacimuftuoglu A. A Review of Melatonin, Its Receptors and Drugs. *Eurasian J Med*. 2016;48(2):135-141. doi:10.5152/eurasianjmed.2015.0267
300. Jehan S, Jean-Louis G, Zizi F, et al. Sleep, Melatonin, and the Menopausal Transition: What Are the Links?. *Sleep Sci*. 2017;10(1):11-18. doi:10.5935/1984-0063.20170003
301. Partsernyak AS, Polyakova VO, Trufanov AG, et al. Melatonin: Manager of psychosomatic and metabolic disorders in polymorbid cardiovascular pathology. *Front Neurosci*. 2022;16:989497. Published 2022 Sep 28. doi:10.3389/fnins.2022.989497
302. Rzepka-Migut B, Paprocka J. Melatonin-Measurement Methods and the Factors Modifying the Results. A Systematic Review of the Literature. *Int J Environ Res Public Health*. 2020 Mar 15;17(6):1916. doi: 10.3390/ijerph17061916. PMID: 32183489; PMCID: PMC7142625.
303. Arendt J, Bojkowski C, Franey C, Wright J, Marks V. Immunoassay of 6-hydroxymelatonin sulfate in human plasma and urine: abolition of the urinary 24-hour rhythm with atenolol. *J Clin Endocrinol Metab*. 1985;60(6):1166-1173. doi:10.1210/jcem-60-6-1166
304. Nowak, Rachel, et al. "The correlation between serum and salivary melatonin concentrations and urinary 6-hydroxymelatonin sulphate excretion rates: two non-invasive techniques for monitoring human circadian rhythmicity." *Clinical endocrinology* 27.4 (1987): 445-452
305. Kennaway DJ. Measuring melatonin by immunoassay. *J Pineal Res*. 2020;69(1):e12657. doi:10.1111/jpi.12657
306. O'Brien KM, Upson K, Buckley JP. Lipid and Creatinine Adjustment to Evaluate Health Effects of Environmental Exposures. *Curr Environ Health Rep*. 2017;4(1):44-50. doi:10.1007/s40572-017-0122-7